

# Electrochemistry accessories



- Analytical cell
- Corrosion cell
- Glassware
- Electrodes
- Rotating electrodes
- Spectroelectrochemistry
- Quartz crystal analyzer
- Battery holders
- Cells for battery testing
- Fuel cell
- Connection accessories
- Local probes accessories

Bio-Logic SAS (founded in 1983) is located near Grenoble in the French Alps.


We design and manufacture high performance, research laboratory instruments and software for electrochemistry, battery testing, fuel cell and material testing.

Bio-Logic's product range is based on a modular and flexible design. We offer electrochemical workstations (EC-Lab) and battery test equipments (BT-Lab) from single channel up to 128-channel systems, integrating EIS capability and high current boosters.

In 2012, Bio-Logic SAS extended its choice of instruments to the scanning electrochemical workstations. The Scan-Lab division is dedicated to high resolution electrochemical mapping tools. These products are developed, designed and manufactured by Uniscan Instruments Ltd., a company based in Macclesfield, UK. For more than 20 years, Uniscan Instruments Ltd. have been at the cutting edge of scanning probe electrochemistry technology.



|                                 |   |    |
|---------------------------------|---|----|
| <b>Analytical cell</b>          | Small volume cells                        | 4  |
|                                 | Large volume cells                        | 6  |
|                                 | Flow cells                                | 7  |
|                                 | InterDigitated Array (IDA) electrodes     | 8  |
|                                 | Ring disk type electrodes                 | 8  |
|                                 | Screen Printed Electrodes (SPE)           | 9  |
| <b>Corrosion cell</b>           | Standard corrosion cells                  | 10 |
|                                 | Avesta cell                               | 11 |
|                                 | Flat cells                                | 12 |
|                                 | Galvanic cells                            | 13 |
|                                 | Plate material evaluating cell            | 13 |
|                                 | Coating cell                              | 13 |
| <b>Glassware</b>                | Small & big volume cell vials             | 14 |
|                                 | Bridge & purge tubes for corrosion cells  | 16 |
|                                 | Small size bridge tubes                   | 16 |
|                                 | Bridge & purge tubes for analytical cells | 17 |
|                                 | Bridge tubes with ceramic junction        | 18 |
| <b>Electrodes</b>               | Working electrodes                        | 18 |
|                                 | Small size reference electrodes           | 20 |
|                                 | King size reference electrodes            | 22 |
|                                 | Counter electrodes                        | 23 |
|                                 | Metallic electrodes                       | 23 |
| <b>Rotating electrodes</b>      | Rotating Ring Disk Electrode (RRDE)       | 24 |
|                                 | Working electrode tips                    | 25 |
|                                 | Rotating Disk Electrode (RDE)             | 26 |
| <b>Spectro-electrochemistry</b> | Spectrometer                              | 28 |
|                                 | Static experiment                         | 30 |
|                                 | Flow experiment                           | 31 |
|                                 | Thin layer spectroelectrochemical cell    | 32 |
| <b>Quartz crystal analyzer</b>  | EQCM                                      | 34 |
|                                 | Quartz resonators                         | 35 |
|                                 | Biosensor QCM                             | 36 |
| <b>Battery accessories</b>      | Test cells                                | 38 |
|                                 | Dilatometers                              | 39 |
|                                 | Rack                                      | 40 |
|                                 | Battery (BH-1) & Coin Cell Holder (CCH)   | 40 |
|                                 | Current Collector (CC5)                   | 41 |
|                                 | Sense Adapter Module (SAM-50)             | 41 |
| <b>Fuel cell</b>                | Membrane Electrode Assemblies (MEA)       | 42 |
|                                 | Educational kit                           | 43 |
| <b>Connection accessories</b>   | Specific cables                           | 44 |
|                                 | Connectors                                | 44 |
|                                 | Multi-electrode investigation cables      | 45 |
|                                 | Faraday cages                             | 45 |
|                                 | Transport cases                           | 46 |
|                                 | External device connection                | 46 |
|                                 | Test Boxes & Dummy cells                  | 46 |
| <b>Local probes accessories</b> | Probes                                    | 47 |
|                                 | Cells                                     | 47 |
|                                 | Video Microscope System                   | 47 |



## Small volume cells

Each voltammetry cell is designed for a specific application (specific working electrode, volume of solution, oxygen free condition...).


For example, for voltammetry investigations using standard working electrode with an outer diameter (OD) of 6 mm the fixed configuration of SVC-3 kit is recommended.

For application requiring other working electrode shape, SVC-2 is more suitable.



If only a small amount of the electroactive compound is available, SVC-2 in the microvolume mode is recommended.

- **SVC-2**, modular
- **SVC-3**, for a volume of 5 to 20 mL, only for working electrode with OD of 6 mm.
- **VC-4**, for a volume of 1 to 3 mL, only for working electrode with OD of 6 mm.
- **Bulk electrolysis cell**, for a volume of 100 mL.

### SVC-2 modularity: microvolume mode



### oxygen free mode




Note that a complete analytical special kit **SK-2** (A-012763) is also proposed, but the reference electrode has to be purchased separately.

This kit includes:

- SVC-3 kit (A-012669),
- PK-3 polishing kit (A-011975) (see page 18),
- one glassy carbon electrode, OD 6.0 mm, ID 3.0 mm (A-002012) (see page 18),
- one platinum electrode, OD 6.0 mm, ID 1.6 mm (A-002013) (see page 18).

| Small volume cells        | Catalog n°        | Catalog n°      | Catalog n°        | Catalog n°                  |
|---------------------------|-------------------|-----------------|-------------------|-----------------------------|
| <b>Products</b>           | <b>SVC-2</b>      | <b>A-012668</b> | <b>SVC-3</b>      | <b>A-012669</b>             |
| <b>Content</b>            |                   |                 |                   |                             |
| Sample vial               | 20 mL (7 pieces)  | <b>A-001056</b> | 20 mL (7 pieces)  | <b>A-001056</b>             |
| Counter electrode (CE)    | 57 mm             | <b>A-002222</b> | 50 mm             | <b>A-002222</b>             |
| PTFE cap                  |                   | <b>A-012670</b> |                   | <b>A-012671</b>             |
| Purge tube (ETFE), 100 mm |                   | -               |                   | -                           |
| Additional items          | adapter 10 to 6   | -               |                   | cell holder <b>A-011227</b> |
| <b>Options</b>            |                   |                 |                   |                             |
| Sample holder             | 9.0 mm (2 pieces) | <b>A-012177</b> | 6.0 mm (2 pieces) | <b>A-012176</b>             |
| Cell holder               | for 20 mL         | <b>A-001209</b> | for 20 mL         | <b>A-001209</b>             |
| Purge tube (ETFE)         | 1 m               | <b>A-010537</b> | 1 m               | <b>A-010537</b>             |
| Working electrodes        |                   | see page 18     |                   | see page 18                 |
| Reference electrodes      |                   | see page 20     |                   | see page 20                 |

| Products                      | Catalog n°      | Catalog n° | Catalog n° | Catalog n° |
|-------------------------------|-----------------|------------|------------|------------|
| <b>Bulk electrolysis cell</b> | <b>A-001197</b> |            |            |            |
| 100 mL (1 piece)              | <b>A-012632</b> |            |            |            |
| 230 mm                        | <b>A-002234</b> |            |            |            |
|                               | <b>A-012551</b> |            |            |            |
|                               | -               |            |            |            |
| porous carbon electrode       | <b>A-010530</b> |            |            |            |
| lid for CE                    | <b>A-001198</b> |            |            |            |
| chamber for CE                | <b>A-001196</b> |            |            |            |
| O-ring                        | <b>A-001236</b> |            |            |            |
| port plug                     | <b>A-009131</b> |            |            |            |
| stirrer bar                   | <b>A-000178</b> |            |            |            |



## Cell geometry

The geometry of the cell should be optimized to reduce the ohmic drop. Working and reference electrodes have to be close to each other. The counter should not be the limiting factor for the electron transfer, so the surface area of the counter electrode should be larger than the surface area of the working electrode.

## Large volume cells



This analytical cell with volumes from 50 to 250 mL is well adapted for macro electrochemistry applications.

Two types of packages are offered:  
- the standard analytical cell kit,  
- the complete analytical cell kit (which allows temperature control and gas purging).

|  | Standard analytical cell kit (80 mL)<br><b>EL-ELECTRO-1</b> | Complete analytical cell kit (80 mL)<br><b>EL-ELECTRO-2</b> | Complete analytical cell kit (150 mL)<br><b>EL-ELECTRO-3</b> | Catalog n°      |
|--|---|---|--|-----------------|
| Glass cell 80 mL                             | ■   |   |  | <b>EL-A-001</b> |
| Double jacketed glass cell 80 mL             |   | ■   |  | <b>EL-A-002</b> |
| Double jacketed glass cell 150 mL            |   |   | ■  | <b>EL-A-020</b> |
| PTFE cap 5 holes                             | ■   | ■   | ■  | <b>EL-A-003</b> |
| PTFE ring, silicon encapsulated, OD 102 mm   | ■   | ■   | ■  | <b>EL-A-004</b> |
| Cell collar with clamp                       | ■   | ■   | ■  | <b>EL-A-005</b> |
| Double purge tube                            |   | ■   | ■  | <b>EL-A-006</b> |
| Bridge tube for reference electrode, OD 6 mm | ■   | ■   | ■  | <b>EL-A-008</b> |
| Platinum counter electrode                   | ■   | ■   | ■  | <b>EL-A-009</b> |
| Purge tube                                   | ■   | ■   | ■  | <b>EL-A-016</b> |
| Reference electrode                          | ■   | ■   | ■  | <b>A-013430</b> |
| Double nut 25 mm and 12 mm diameter          |   | ■   | ■  | <b>EL-A-011</b> |
| Telescopic cell stand                        |   | ■   | ■  | <b>EL-A-012</b> |

| Options   |       |                  |
|---|-------|------------------|
| Electrode bridge extension for electroanalytical cell                       |       | <b>EL-A-022</b>  |
| Bridge tube for reference electrode of OD 8 mm                              |       | <b>EL-A-017</b>  |
| PT100 probe, indicate connector type  |       | <b>EL-C-014</b>  |
| Magnetic stirrer & heater, without PT100 probe                              | 220 V | <b>EL-C-015A</b> |
|   | 110 V | <b>EL-C-015B</b> |
| Aluminum base holder for magnetic stirrer                                   |       | <b>EL-C-018</b>  |
| Set of 10 porous 4 mm glass frit (CoralPor™) with PTFE heat shrink (200 mm) |       | <b>092-VYC4</b>  |

250 mL cells are also available: glass cell EL-A-018 or double jacketed glass cell EL-A-021




## Flow cells

|                   | Cross flow cell  | Catalog n°      | Radial flow cell  | Catalog n°      |
|-------------------|------------------|-----------------|-------------------|-----------------|
| Product           | Cross flow cell* | <b>A-012798</b> | Radial flow cell* | <b>A-012799</b> |
| Schematic diagram |                  |                 |                   |                 |


\* Gasket, working, reference electrodes and connectors are sold separately

| Options                   |   |  |   |                                       |                 |
|---------------------------|---|--|---|---------------------------------------|-----------------|
| Cross flow                | electrode   | glassy carbon                              | Ø 3 mm x 2 (dual type), size 25 x 25 mm | <b>A-001000</b>                       |                 |
|                           |   | gold                                       | Ø 3 mm x 2 (dual type), size 25 x 25 mm | <b>A-001002</b>                       |                 |
|                           |   | platinum                                   | Ø 3 mm x 2 (dual type), size 25 x 25 mm | <b>A-001012</b>                       |                 |
|                           |   | silver                                     | Ø 3 mm x 2 (dual type), size 25 x 25 mm | <b>A-001008</b>                       |                 |
|                           |   | carbon paste                               | Ø 3 mm x 2 (dual type), size 25 x 25 mm | <b>A-001004</b>                       |                 |
|                           |   | nickel                                     | Ø 3 mm x 2 (dual type), size 25 x 25 mm | <b>A-001009</b>                       |                 |
|                           |   | glassy carbon/gold                         | Ø 3 mm x 2 (dual type), size 25 x 25 mm | <b>A-001006</b>                       |                 |
|                           |   | glassy carbon/platinum                     | Ø 3 mm x 2 (dual type), size 25 x 25 mm | <b>A-012583</b>                       |                 |
|                           |   | gasket                                     | TG-2M PTFE (4 pieces)                   | film thickness 12 µm                  | <b>A-001046</b> |
|                           |   |  | TG-5M PTFE (4 pieces)                   | film thickness 25 µm                  | <b>A-001047</b> |
| TG-6M PTFE (4 pieces)     | film thickness 50 µm  |  | <b>A-001048</b>                         |                                       |                 |
| TG-8M PTFE (4 pieces)     | film thickness 100 µm   |  | <b>A-001049</b>                         |                                       |                 |
| TG-11M silicon (4 pieces) | film thickness 500 µm   |  | <b>A-001092</b>                         |                                       |                 |
| TG-12M silicon (4 pieces) | film thickness 1,000 µm   |  | <b>A-001093</b>                         |                                       |                 |
| cross flow cell block     | material PEEK   |  | <b>A-001032</b>                         |                                       |                 |
| Radial flow               | electrode   |  | glassy carbon                           | Ø 6 mm (single type), size 25 x 25 mm | <b>A-001016</b> |
|                           |   | glassy carbon                              | Ø 3 mm (single type), size 25 x 25 mm   | <b>A-012124</b>                       |                 |
|                           |   | platinum                                   | Ø 3 mm (single type), size 25 x 25 mm   | <b>A-009908</b>                       |                 |
|                           |   | gold                                       | Ø 3 mm (single type), size 25 x 25 mm   | <b>A-011155</b>                       |                 |
|                           |   | PFCE carbon                                | Ø 3 mm (single type), size 25 x 25 mm   | <b>A-000999</b>                       |                 |
|                           |   | carbon paste                               | Ø 3 mm (single type), size 25 x 25 mm   | <b>A-010251</b>                       |                 |
|                           |   | gasket                                     | TG-2MR PTFE (4 pieces)                  | film thickness 12 µm                  | <b>A-001146</b> |
|                           |   |  | TG-5MR PTFE (4 pieces)                  | film thickness 25 µm                  | <b>A-001147</b> |
|                           |   |  | TG-6MR PTFE (4 pieces)                  | film thickness 50 µm                  | <b>A-001148</b> |
|                           |   |  | TG-8MR PTFE (4 pieces)                  | film thickness 100 µm                 | <b>A-012802</b> |
| radial flow cell block    |   | material PEEK                              | <b>A-001193</b>                         |                                       |                 |
|                           |   | material PEEK                              | <b>A-001031</b>                         |                                       |                 |
| Cross/radial              | RE-3V Ag/AgCl reference electrode screw type  | size Ø 10 x 48 mm                          | <b>A-012169</b>                         |                                       |                 |
|                           | RE-7V non aqueous reference electrode [Ag/Ag <sup>+</sup> ] screw type, poly acetal resin |  | <b>A-012173</b>                         |                                       |                 |
|                           | RE-3VP Ag/AgCl screw type electrode, PEEK resin   |  | <b>A-012170</b>                         |                                       |                 |
|                           | RE-7VP non aqueous reference electrode [Ag/Ag <sup>+</sup> ] screw type, PEEK resin       |  | <b>A-012174</b>                         |                                       |                 |
|                           | 0.04" single leal connector (2 pieces)  |  | <b>A-012912</b>                         |                                       |                 |
|                           | 0.04" MM connector (10 pieces)  |  | <b>A-013273</b>                         |                                       |                 |
|                           | 1/16" PEEK tube   | ID 0.25 mm, length 3.0 m                   | <b>A-001531</b>                         |                                       |                 |
|                           | dynaseal PEEK fingertight (10 pieces)   | screw for pipe connecting integrated 1/16" | <b>A-004130</b>                         |                                       |                 |





## InterDigitated Array (IDA) electrodes



The passivation membrane is a Novolac resin + naphthoquinore-diazido compounds.

| IDA electrode                         | Width | Interval | Length | N° of feet | Film thickness | Catalog n°      |
|---------------------------------------|-------|----------|--------|------------|----------------|-----------------|
| Gold                                  | 2 µm  | 2 µm     | 2 mm   | 65 pairs   | 90 nm          | <b>A-012257</b> |
| Platinum                              | 2 µm  | 2 µm     | 2 mm   | 65 pairs   | 90 nm          | <b>A-012258</b> |
| Gold without passivation membrane     | 2 µm  | 2 µm     | 2 mm   | 65 pairs   | 90 nm          | <b>A-012261</b> |
| Platinum without passivation membrane | 2 µm  | 2 µm     | 2 mm   | 65 pairs   | 90 nm          | <b>A-012264</b> |
| Gold                                  | 3 µm  | 3 µm     | 2 mm   | 65 pairs   | 90 nm          | <b>A-012129</b> |
| Platinum                              | 3 µm  | 3 µm     | 2 mm   | 65 pairs   | 90 nm          | <b>A-012130</b> |
| Gold without passivation membrane     | 3 µm  | 3 µm     | 2 mm   | 65 pairs   | 90 nm          | <b>A-012260</b> |
| Platinum without passivation membrane | 3 µm  | 3 µm     | 2 mm   | 65 pairs   | 90 nm          | <b>A-012263</b> |
| Gold                                  | 10 µm | 5 µm     | 2 mm   | 65 pairs   | 90 nm          | <b>A-012125</b> |
| Platinum                              | 10 µm | 5 µm     | 2 mm   | 65 pairs   | 90 nm          | <b>A-012126</b> |
| ITO                                   | 10 µm | 5 µm     | 2 mm   | 65 pairs   | 100 ±20 nm     | <b>A-012128</b> |
| Carbon                                | 10 µm | 5 µm     | 2 mm   | 65 pairs   | 1200 ±100 nm   | <b>A-012127</b> |
| Gold without passivation membrane     | 10 µm | 5 µm     | 2 mm   | 65 pairs   | 90 nm          | <b>A-012259</b> |
| Platinum without passivation membrane | 10 µm | 5 µm     | 2 mm   | 65 pairs   | 90 nm          | <b>A-012262</b> |
| ITO without passivation membrane      | 10 µm | 5 µm     | 2 mm   | 65 pairs   | 100 ±200 nm    | <b>A-012265</b> |
| Carbon without passivation membrane   | 10 µm | 5 µm     | 2 mm   | 65 pairs   | 1200 ±100 nm   | <b>A-012266</b> |

| Options                                      | Catalog n°      |
|--|-----------------|
| Cable kit for IDA electrode                  | <b>A-011066</b> |
| Ag/AgCl ink for reference electrode (2.0 mL) | <b>A-011464</b> |



## Ring disk type electrodes


| Ring disk type electrodes               | Ring OD | Ring ID | Disk | Catalog n°      |
|---|---------|---------|------|-----------------|
| Gold ring disk electrode (3 pieces)     | 6 mm    | 4 mm    | 3 mm | <b>A-002081</b> |
| Platinum ring disk electrode (3 pieces) | 6 mm    | 4 mm    | 3 mm | <b>A-002082</b> |
| Carbon ring disk electrode (3 pieces)   | 6 mm    | 4 mm    | 3 mm | <b>A-002083</b> |

## Screen Printed Electrodes (SPE)

Screen Printed Electrodes (SPE) are low cost and disposable. They include one reference electrode, one counter electrode and one or several working electrodes according to the design needed. The working electrode can be made of glassy carbon, platinum, gold. The SPE allows users to work with a small amount of solution (25-100 µL).

A two-electrode set-up is used for multi working electrode SPE with one shared reference electrode and one shared counter electrode.

Note with Bio-Logic's multichannel potentiostats, it is possible to run simultaneous investigations on each working electrode of the array (up to 16).



| Screen Printed Electrodes | Working electrode | Counter electrode | Ref. electrode | Substrat  | Packaging  | Catalog n°      |
|---------------------------|-------------------|-------------------|----------------|-----------|------------|-----------------|
| 3-electrode               | graphite          | graphite          | Ag/AgCl        | alumina   | 20 pieces  | <b>G-GGA/0</b>  |
| ∅ of WE:                  | 2 mm              |                   |                |           |            |                 |
| size of RE:               | 0.5 x 1 mm        |                   |                |           |            |                 |
| ratio CE/WE surface area: | 4/1               |                   |                |           |            |                 |
| overall dimension:        | 50 x 10 mm        |                   |                |           |            |                 |
|                           | graphite          | graphite          | Ag/AgCl        | alumina   | 100 pieces | <b>G-GGA/1</b>  |
|                           | graphite          | graphite          | Ag/AgCl        | alumina   | 500 pieces | <b>G-GGA/2</b>  |
|                           | graphite          | graphite          | Ag/AgCl        | valox FR1 | 20 pieces  | <b>G-GGP/0</b>  |
|                           | graphite          | graphite          | Ag/AgCl        | valox FR1 | 100 pieces | <b>G-GGP/1</b>  |
|                           | graphite          | graphite          | Ag/AgCl        | valox FR1 | 500 pieces | <b>G-GGP/2</b>  |
|                           | platinum          | graphite          | Ag/AgCl        | alumina   | 100 pieces | <b>G-PTGA/1</b> |
|                           | platinum          | graphite          | Ag/AgCl        | alumina   | 500 pieces | <b>G-PTGA/2</b> |
|                           | platinum          | graphite          | Ag/AgCl        | valox FR1 | 100 pieces | <b>G-PTGP/1</b> |
|                           | platinum          | graphite          | Ag/AgCl        | valox FR1 | 500 pieces | <b>G-PTGP/2</b> |
|                           | gold              | graphite          | Ag/AgCl        | alumina   | 20 pieces  | <b>G-AUGA/0</b> |
|                           | gold              | graphite          | Ag/AgCl        | alumina   | 100 pieces | <b>G-AUGA/1</b> |
|                           | gold              | graphite          | Ag/AgCl        | alumina   | 500 pieces | <b>G-AUGA/2</b> |
|                           | gold              | graphite          | Ag/AgCl        | valox FR1 | 20 pieces  | <b>G-AUGP/0</b> |
|                           | gold              | graphite          | Ag/AgCl        | valox FR1 | 100 pieces | <b>G-AUGP/1</b> |
|                           | gold              | graphite          | Ag/AgCl        | valox FR1 | 500 pieces | <b>G-AUGP/2</b> |

| Option   | Catalog n°   |
|--|--------------|
| Connector for 3-electrode Screen Printed Electrodes (2 mm banana plug) | <b>U-GEM</b> |

| Option             | Working electrode | Counter electrode & ref. electrode | Substrat  | Packaging  | Catalog n°       |
|--------------------|-------------------|------------------------------------|-----------|------------|------------------|
| Array of 16 WE     | graphite          | Ag/AgCl                            | alumina   | 100 pieces | <b>G-A-GA/1</b>  |
| ∅ of WE:           | 0.8 mm            |                                    |           |            |                  |
| size of pseudo-RE: | 20 x 1 mm         |                                    |           |            |                  |
| overall dimension: | 45 x 19 mm        |                                    |           |            |                  |
|                    | graphite          | Ag/AgCl                            | alumina   | 500 pieces | <b>G-A-GA/2</b>  |
|                    | graphite          | Ag/AgCl                            | valox FR1 | 100 pieces | <b>G-A-GP/1</b>  |
|                    | graphite          | Ag/AgCl                            | valox FR1 | 500 pieces | <b>G-A-GP/2</b>  |
|                    | platinum          | Ag/AgCl                            | alumina   | 100 pieces | <b>G-A-PTA/1</b> |
|                    | platinum          | Ag/AgCl                            | alumina   | 500 pieces | <b>G-A-PTA/2</b> |
|                    | platinum          | Ag/AgCl                            | valox FR1 | 100 pieces | <b>G-A-PTP/1</b> |
|                    | platinum          | Ag/AgCl                            | valox FR1 | 500 pieces | <b>G-A-PTP/2</b> |
|                    | gold              | Ag/AgCl                            | alumina   | 100 pieces | <b>G-A-AUA/1</b> |
|                    | gold              | Ag/AgCl                            | alumina   | 500 pieces | <b>G-A-AUA/2</b> |
|                    | gold              | Ag/AgCl                            | valox FR1 | 100 pieces | <b>G-A-AUP/1</b> |
|                    | gold              | Ag/AgCl                            | valox FR1 | 500 pieces | <b>G-A-AUP/2</b> |

## Standard corrosion cells

For corrosion applications, two packages for large volume i.e. ~1 L are available:

- basic corrosion cell kit,
- complete corrosion cell kit (allows temperature control and includes sample holder and cell stand).




|   | Basic corrosion cell kit<br><b>EL-CORR-1</b> | Complete corrosion cell kit<br><b>EL-CORR-2</b> | Catalog n°       |
|---|--|---|------------------|
| Glass cell 1 L  | ■  |   | <b>EL-C-001</b>  |
| Double jacketed glass cell 1 L                                    |  | ■   | <b>EL-C-002</b>  |
| PTFE cap  | ■  | ■   | <b>EL-C-003</b>  |
| PTFE ring, silicon encapsulated, OD 102 mm                        | ■  | ■   | <b>EL-C-004</b>  |
| Cell collar with clamp  | ■  | ■   | <b>EL-C-005</b>  |
| Double purge tube   |  | ■   | <b>EL-C-006</b>  |
| Graphite counter electrode rod (2 pieces)                         | ■  | ■   | <b>EL-C-009</b>  |
| Double nut 25 mm and 12 mm diameter                               |  | ■   | <b>EL-C-011</b>  |
| Telescopic cell stand   |  | ■   | <b>EL-C-012</b>  |
| Sample holder 1 cm <sup>2</sup>                                   |  | ■   | <b>092-C-013</b> |
| Purge tube  | ■  |   | <b>EL-C-016</b>  |
| Bridge tube for 8 mm diameter reference electrode                 | ■  | ■   | <b>EL-C-017</b>  |
| Calomel reference electrode length 80 mm, OD 8 mm screw cap       | ■  | ■   | <b>R-XR110</b>   |
| Cable connection for screw cap 100 mm, 2 mm banana plug           | ■  | ■   | <b>R-A94L111</b> |
| 12 mm OD conical rings for reference electrode of 8 mm (4 pieces) | ■  | ■   | <b>R-X31M012</b> |

### Options

|   |       |                  |
|---|-------|------------------|
| Bridge tube for 6 mm diameter reference electrode                     |       | <b>EL-C-008</b>  |
| PT100 probe (indicate connector type)                                 |       | <b>EL-C-014</b>  |
| Magnetic stirrer & heater, without PT100 probe                        | 220 V | <b>EL-C-015A</b> |
|   | 110 V | <b>EL-C-015B</b> |
| Aluminum base holder for magnetic stirrer                             |       | <b>EL-C-018</b>  |
| Set of 10 porous frit (4 mm CoralPor™) with PTFE heat shrink (200 mm) |       | <b>092-VYC4</b>  |

500 mL or 2 L double jacketed glass cell are also available (EL-C-020 and EL-C-021, respectively).



## Avesta cell

Avesta cell is an electrochemical cell developed for pitting corrosion testing (ASTM G150).

It is designed to avoid microcrevice corrosion formed between the working electrode and the gasket at the bottom aperture of the cell.

A filter paper ring placed between the sample and the gasket is flooded by distilled water in order to eliminate crevice corrosion.

The water flow is controlled by a peristaltic pump (EL-AV-008) which delivers 0.5 to 5 mL/h.

| Avesta cell     | Catalog n°       |
|-----------------|------------------|
| Avesta cell kit | <b>EL-AV-001</b> |



### Content

|  |                  |
|--|------------------|
| Double jacketed cell glass 250 mL          | <b>EL-AV-002</b> |
| PTFE cap 5 holes                           | <b>EL-AV-003</b> |
| O-ring PTFE silicone encapsulated          | <b>EL-AV-004</b> |
| Cell collar with clamp                     | <b>EL-A-005</b>  |
| Double purge tube                          | <b>EL-AV-006</b> |
| Filter paper ring (100 pieces)             | <b>EL-AV-007</b> |
| Graphite counter electrode rods (2 pieces) | <b>EL-C-009</b>  |
| Bridge tube for RE 6 mm                    | <b>EL-AV-010</b> |
| Double nut                                 | <b>EL-C-011</b>  |
| Telescopic cell stand                      | <b>EL-AV-012</b> |
| Skeleton                                   | <b>EL-AV-015</b> |

### Options

|   |                  |
|---|------------------|
| Peristaltic pump for low flow   | <b>EL-AV-008</b> |
| Bridge tube for reference electrode with OD of 8 mm                   | <b>EL-AV-013</b> |
| Single purge tube   | <b>EL-AV-014</b> |
| Temperature probe PT100   | <b>EL-C-014</b>  |
| Set of 10 porous frit (4 mm CoralPor™) with PTFE heat shrink (200 mm) | <b>092-VYC4</b>  |

Principle of Avesta cell





## Flat cells, 1 to 10 cm<sup>2</sup> sample area

This cell with a volume of 250 mL is adapted to perform AC or DC experiments on flat specimen of 1 cm<sup>2</sup> or 10 cm<sup>2</sup> area.

This cell has a double jacket for temperature control and three holes for reference electrode, purge tube and temperature probe with inner diameter of 17.6 mm and the 2 others of 8.3 mm.

Thanks to its design, any sample with a thickness inferior to 10 mm fits in the holder.

| Flat cell 1 cm <sup>2</sup>                                       | Catalog n°        | Flat cell 10 cm <sup>2</sup>                                      | Catalog n°        |
|---|-------------------|---|-------------------|
| Flat cell kit 1 cm <sup>2</sup>                                   | <b>EL-FLAT</b>    | Flat cell kit 10 cm <sup>2</sup>                                  | <b>EL-FLAT-2</b>  |
| <b>Content</b>  |                   | <b>Content</b>  |                   |
| RE-2B calomel reference electrode (length: 90 mm, OD: 6 mm)       | <b>A-013430</b>   | RE-2B calomel reference electrode (length: 90 mm, OD: 6 mm)       | <b>A-013430</b>   |
| Platinum mesh counter electrode (54 mm wire/ 80 mesh), 25 x 35 mm | <b>A-702439</b>   | Platinum mesh counter electrode (54 mm wire/ 80 mesh), 25 x 35 mm | <b>A-702439</b>   |
| Reference electrode bridge tube (6 mm diameter)                   | <b>EL-F-004</b>   | Reference electrode bridge tube (6 mm diameter)                   | <b>EL-F-004</b>   |
| Corrosion flat cell 1cm <sup>2</sup>                              | <b>EL-FLAT-3</b>  | Corrosion flat cell 10 cm <sup>2</sup>                            | <b>EL-FLAT-4</b>  |
| EL-FLAT-3 includes:   |                   | EL-FLAT-4 includes:   |                   |
| - glass part for flat cell (250 mL)                               | <b>EL-F-002</b>   | - glass part for flat cell (250 mL)                               | <b>EL-F-002</b>   |
| - mechanical part for corrosion flat cell 1 cm <sup>2</sup>       | <b>EL-FLAT-3H</b> | - mechanical part for corrosion flat cell 10 cm <sup>2</sup>      | <b>EL-FLAT-4H</b> |

## Investigation in aggressive media

If the experiment is performed in more aggressive media such as fluorhydric acid media, it is possible to get a body of the flat cell in PVDF\* instead of glass.

| Temperature         | Fluorhydric acid 48% |       | Sulfuric acid 98% |       | Phosphoric acid 85% |       | Hydrochloric acid 35% |       | Nitric acid 70% |       | Perchlorique acid |       | Sodium hydroxide 50% |       | Potassium hydroxide concentrated |       |
|---------------------|----------------------|-------|-------------------|-------|---------------------|-------|-----------------------|-------|-----------------|-------|-------------------|-------|----------------------|-------|----------------------------------|-------|
|                     | 20° C                | 50° C | 20° C             | 50° C | 20° C               | 50° C | 20° C                 | 50° C | 20° C           | 50° C | 20° C             | 50° C | 20° C                | 50° C | 20° C                            | 50° C |
| PTFE                | ■                    | ■     | ■                 | ■     | ■                   | ■     | ■                     | ■     | ■               | ■     | ■                 | ■     | ■                    | ■     | ■                                | ■     |
| PVDF*               | ■                    | ■     | ■                 | ■     | ■                   | ■     | ■                     | ■     | ■               | ■     | ■                 | ■     | ■                    | ■     | ■                                | ■     |
| Borosilicated glass | ■                    | ■     | ■                 | ■     | ■                   | ■     | ■                     | ■     | ■               | ■     | ■                 | ■     | ■                    | ■     | ■                                | ■     |

| Options   | Catalog n°          |
|---|---------------------|
| PVDF* body (non double jacketed)                                      | <b>EL-F-PVDF</b>    |
| EPDM O-rings for 1 cm <sup>2</sup> with PEEK ferrule for CE**         | <b>EL-SEAL-1B</b>   |
| EPDM O-rings for 10 cm <sup>2</sup> with PEEK ferrule for CE**        | <b>EL-SEAL-10B</b>  |
| PTFE O-rings for 1 cm <sup>2</sup> with PEEK ferrule for CE**         | <b>EL-SEAL-T1B</b>  |
| PTFE O-rings for 10 cm <sup>2</sup> with PEEK ferrule for CE**        | <b>EL-SEAL-T10B</b> |
| EPDM O-rings for 1 cm <sup>2</sup> **                                 | <b>EL-SEAL-1</b>    |
| EPDM O-rings for 10 cm <sup>2</sup> **                                | <b>EL-SEAL-10</b>   |
| PTFE O-rings for 1 cm <sup>2</sup> **                                 | <b>EL-SEAL-T1</b>   |
| PTFE O-rings for 10 cm <sup>2</sup> **                                | <b>EL-SEAL-T10</b>  |
| Set of 10 porous frit (4 mm CoralPor™) with PTFE heat shrink (200 mm) | <b>092-VYC4</b>     |

\* Polyvinylidene fluoride  
\*\* The O-ring kits include 4 O-rings for the glassware side and 10 O-rings for the sample side

■ excellent  
■ good  
■ not recommended  
■ not compatible  
\* Polyvinylidene fluoride

## Galvanic cells, 1 to 10 cm<sup>2</sup> sample area

Thanks to the modular design of the flat cell, it is possible to place two different materials at each end of the cell.

The surface area may be 1 or 10 cm<sup>2</sup>.

| Galvanic cell 1 cm <sup>2</sup>     | Catalog n°        | Galvanic cell 10 cm <sup>2</sup>          | Catalog n°         |
|-------------------------------------|-------------------|---|--------------------|
| Galvanic cell kit 1 cm <sup>2</sup> | <b>EL-GAL-1</b>   | Galvanic flat cell kit 10 cm <sup>2</sup> | <b>EL-GAL-10</b>   |
| <b>Content</b>                      |                   | <b>Content</b>                            |                    |
| Flat cell kit 1 cm <sup>2</sup>     | <b>EL-FLAT</b>    | Flat cell kit 10 cm <sup>2</sup>          | <b>EL-FLAT-2</b>   |
| Galvanic kit 1 cm <sup>2</sup>      | <b>092-FLAT/1</b> | Galvanic kit 10 cm <sup>2</sup>           | <b>092-FLAT/10</b> |



## Plate material evaluating cell, up to 1 cm<sup>2</sup> sample area

This cell was developed to evaluate a plate material such as a metal, semi-conducting plate...

The sample plate is sandwiched between the two cell blocks.

The required volume of solution is about 1 mL.

| Plate material evaluating cell       | Catalog n°      |
|--------------------------------------|-----------------|
| Plate material evaluating cell       | <b>A-011951</b> |
| <b>Content</b>                       |                 |
| PTFE cell (body & base) (1 piece)    | -               |
| O-ring (1 piece)                     | -               |
| Screw 20 mm (1 piece)                | -               |
| Purging tube, 100 mm                 | -               |
| Platinum counter electrode (1 piece) | <b>A-002222</b> |
| <b>Options</b>                       |                 |
| O-ring (10 pieces)                   | <b>A-012022</b> |

## Coating cell

| Coating cell     | Catalog n°     |
|------------------|----------------|
| Coating cell kit | <b>EL-COAT</b> |

| Content  |                 |
|--|-----------------|
| Glass for coating cell                               | <b>EL-P-002</b> |
| Nylon base with three feet                           | <b>EL-P-003</b> |
| Rubber cup with two holes                            | <b>EL-P-004</b> |
| Metallic clamp                                       | <b>EL-P-005</b> |
| O-ring for coating cell                              | <b>EL-P-006</b> |
| Graphite rod counter electrode (L: 145 mm, OD: 6 mm) | <b>EL-P-009</b> |

| Options                                  |                 |
|--|-----------------|
| Ag/AgCl reference electrode (OD: 6mm)    | <b>A-012167</b> |
| Bridge tube for 6 mm reference electrode | <b>EL-P-008</b> |
| Mask for 1 cm <sup>2</sup> (20 pieces)   | <b>EL-P-011</b> |
| Mask for 5 cm <sup>2</sup> (20 pieces)   | <b>EL-P-012</b> |
| Mask for 10 cm <sup>2</sup> (20 pieces)  | <b>EL-P-013</b> |




## Small volume cell vials, 1 to 100 mL

To complete a kit, the cell vials are also offered separately.

Note that other volumes are available on request.



ID 15.6  
A-011504




A-001209  
cell holder for A-011504




ID 25.6  
A-001056



ID 46.4  
A-012632



A-012672



A-001051



A-012652

## Big volume cell vials, 100 to 2,000 mL



EL-A-001



EL-A-002



EL-C-001



EL-C-002



EL-C-019,  
cell holder for EL-C-001

| Small & big volume cell vials        |                           | Volume   | OD     | ID      | Height | Quantity     | Purpose                    | Catalog n°      |
|--------------------------------------|---------------------------|----------|--------|---------|--------|--------------|----------------------------|-----------------|
| Small volume cell vials              | sample vial               | 5 mL     | 18 mm  | 15.6 mm | 30 mm  | 10           | VC-4                       | <b>A-011504</b> |
|                                      |                           | 20 mL    | 28 mm  | 25.6 mm | 50 mm  | 10           | SVC-2, SVC-3               | <b>A-001056</b> |
|                                      |                           | 100 mL   | 50 mm  | 46.4 mm | 72 mm  | 1            | RRDE-3A, bulk electrolysis | <b>A-012632</b> |
| water jacketed-glass vial            |                           | 5 mL     | 40 mm  | 15.6 mm | 40 mm  | 1            | VC-4                       | <b>A-012672</b> |
|                                      |                           | 20 mL    | 55 mm  | 25.6 mm | 50 mm  | 1            | SVC-2, SVC-3               | <b>A-001051</b> |
|                                      |                           | 100 mL   | 70 mm  | 46.4 mm | 80 mm  | 1            | RRDE-3A, bulk electrolysis | <b>A-012652</b> |
| cell holder for 20 mL vial           |                           |          |        |         | 1      | SVC-2, SVC-3 | <b>A-001209</b>            |                 |
| Big volume cell vials                | sample vial               | 80 mL    | 90 mm  | 62 mm   | 80 mm  | 1            | EL-ELECTRO-1               | <b>EL-A-001</b> |
|                                      |                           | 1,000 mL | 120 mm | 90 mm   | 175 mm | 1            | EL-CORR-1                  | <b>EL-C-001</b> |
|                                      | water jacketed-glass vial | 80 mL    | 90 mm  | 62 mm   | 85 mm  | 1            | EL-ELECTRO-2               | <b>EL-A-002</b> |
|                                      | 1,000 mL                  | 120 mm   | 90 mm  | 200 mm  | 1      | EL-CORR-2    | <b>EL-C-002</b>            |                 |
| cell holder for king size vial (1 L) |                           |          |        |         | 1      | EL-CORR-1    | <b>EL-C-019</b>            |                 |


Dimensions in mm



### Bridge tubes for corrosion cells



### Small size bridge tubes



### Purge tubes for corrosion cells



The tolerance of each dimension is approximately ±0.5 mm.


### Bridge tubes for analytical cells



### Purge tubes for analytical cells



### Bridge tubes with ceramic junction




| Glassware                         | Ref. electrode                    | OD     | Height | Quantity | Purpose      | Vycor type compatibility | Catalog n°      |
|-----------------------------------|-----------------------------------|--------|--------|----------|--------------|--------------------------|-----------------|
| Small size bridge tube            | Ø 6 mm                            | 6 mm   | 68 mm  | 2        | SVC-2, SVC-3 | 092-VYC3 <sup>[1]</sup>  | <b>A-012176</b> |
|                                   | Ø 6 mm                            | 6 mm   | 68 mm  | 22       | SVC-2, SVC-3 |                          | <b>A-012306</b> |
|                                   | Ø 9 mm                            | 9 mm   | 68 mm  | 2        | SVC-2        |                          | <b>A-012177</b> |
|                                   | Ø 9 mm                            | 9 mm   | 68 mm  | 22       | SVC-2        |                          | <b>A-012307</b> |
| Bridge tube for corrosion cells   | Ø 6 mm                            | 10 mm  | 250 mm | 1        | EL-CORR      | 092-VYC4 <sup>[2]</sup>  | <b>EL-C-008</b> |
|                                   | Ø 8 mm                            | 10 mm  | 250 mm | 1        | EL-CORR      |                          | <b>EL-C-017</b> |
|                                   | 10 mm                             | 120 mm | 1      | EL-FLAT  |              | <b>EL-F-004</b>          |                 |
|                                   | 10 mm                             | 140 mm | 1      | EL-FLAT  |              | <b>EL-F-004B</b>         |                 |
| Bridge tube for analytical cells  | Ø 8 mm                            | 10 mm  | 165 mm | 1        | EL-ELECTRO   |                          | <b>EL-A-017</b> |
|                                   | Ø 6 mm                            | 10 mm  | 155 mm | 1        | EL-ELECTRO   |                          | <b>EL-A-008</b> |
| Purge tube for analytical cells   | single                            | 10 mm  | 165 mm | 1        | EL-ELECTRO   |                          | <b>EL-A-016</b> |
|                                   | double                            | 10 mm  | 200 mm | 1        | EL-ELECTRO   |                          | <b>EL-A-006</b> |
| Purge tube for corrosion cells    | single                            | 10 mm  | 200 mm | 1        | EL-CORR      |                          | <b>EL-C-016</b> |
|                                   | double                            | 10 mm  | 220 mm | 1        | EL-CORR      |                          | <b>EL-C-006</b> |
| Bridge tube with ceramic junction | short                             | 12 mm  | 70 mm  | 1        |              |                          | <b>R-AL100</b>  |
|                                   | standard                          | 8 mm   | 250 mm | 1        |              |                          | <b>R-AL110</b>  |
|                                   | with reverse sleeve (non aqueous) | 8 mm   | 138 mm | 1        |              |                          | <b>R-AL210</b>  |
|                                   | standard                          | 8 mm   | 138 mm | 1        |              |                          | <b>R-AL120</b>  |
|                                   |                                   |        |        |          |              |                          |                 |

#### Options

|   |                 |
|---|-----------------|
| [1] Set of 10 porous 2.8 mm glass frit (CoralPor™) with PTFE heat shrink (200 mm) | <b>092-VYC3</b> |
| [2] Set of 10 porous 4 mm glass frit (CoralPor™) with PTFE heat shrink (200 mm)   | <b>092-VYC4</b> |

Dimensions in mm

## Working electrodes



To address every application, a wide range of working electrodes (WE) is offered with diameters ranging from 7  $\mu\text{m}$  up to 6 mm and built with different materials.


## Working electrode maintenance

To refresh the surface of the electrode we recommend polishing the electrode after each measurement.

| Working electrodes maintenance               |  | Catalog n°      |
|--|--|-----------------|
| <b>PK-3 electrode polishing kit</b>          |  |                 |
| <b>Content</b>                               |  |                 |
| 0.05 $\mu\text{m}$ polishing alumina (20 mL) |  | <b>A-001050</b> |
| 1 $\mu\text{m}$ polishing diamond (10 mL)    |  | <b>A-002054</b> |
| Glass plate (1 piece)                        |  | <b>A-002249</b> |
| Alumina polishing pad (10 pieces)            |  | -               |
| Diamond polishing pad (10 pieces)            |  | -               |
| <b>Spare parts</b>                           |  |                 |
| Alumina polishing pad (20 pieces)            |  | <b>A-001040</b> |
| Diamond polishing pad (20 pieces)            |  | <b>A-001041</b> |
| Emery paper UF800 (20 pieces)                |  | <b>A-012611</b> |
| Coarse polishing pad (20 pieces)             |  | <b>A-001042</b> |
| 6 $\mu\text{m}$ polishing diamond (10 mL)    |  | <b>A-002053</b> |




Dimensions in mm



## Working electrodes

|              |   | Isolation                 | OD       | Electrode size $\pm 4\%$ | Catalog n°      |                 |
|--------------|---|---------------------------|----------|--------------------------|-----------------|-----------------|
| Carbon       | long glassy carbon electrode                    | PEEK                      | 6 mm     | 3 mm                     | <b>A-012744</b> |                 |
|              | standard glassy carbon electrode                | PEEK                      | 10 mm    | 5 mm                     | <b>A-002417</b> |                 |
|              |   | PEEK                      | 6 mm     | 3 mm                     | <b>A-002012</b> |                 |
|              |   | PEEK                      | 6 mm     | 1.6 mm                   | <b>A-012297</b> |                 |
|              |   | PEEK                      | 6 mm     | 1 mm                     | <b>A-002411</b> |                 |
|              | small glassy carbon electrode                   | PEEK                      | 3 mm     | 1.6 mm                   | <b>A-012298</b> |                 |
|              |   | PEEK                      | 3 mm     | 1 mm                     | <b>A-002412</b> |                 |
|              | micro carbon fiber electrode                    | glass                     | 4 mm     | 33 $\mu\text{m}$         | <b>A-002002</b> |                 |
|              |   | glass                     | 4 mm     | 7 $\mu\text{m}$          | <b>A-002007</b> |                 |
|              | standard pyrolytic graphite electrode           | basal plane               | PEEK     | 6 mm                     | 3 mm            | <b>A-002252</b> |
|              |   | edge plane                | PEEK     | 6 mm                     | 3 mm            | <b>A-002253</b> |
|              | standard PFCE carbon electrode                  | PEEK                      | 6 mm     | 3 mm                     | <b>A-002408</b> |                 |
|              |   | PEEK                      | 6 mm     | 1 mm                     | <b>A-002409</b> |                 |
|              | small PFCE carbon electrode                     | PEEK                      | 3 mm     | 1 mm                     | <b>A-011854</b> |                 |
| Platinum     | platinum gauze electrode                        | PEEK                      | 80 mesh  | 25 x 35 mm               | <b>A-002250</b> |                 |
|              | long platinum electrode                         | PEEK                      | 6 mm     | 3 mm                     | <b>A-012745</b> |                 |
|              | standard platinum electrode                     | PEEK                      | 10 mm    | 5 mm                     | <b>A-002420</b> |                 |
|              |   | PEEK                      | 6 mm     | 3 mm                     | <b>A-002422</b> |                 |
|              |   | PEEK                      | 6 mm     | 1.6 mm                   | <b>A-002013</b> |                 |
|              | small platinum electrode                        | PEEK                      | 3 mm     | 1.6 mm                   | <b>A-002313</b> |                 |
|              | micro platinum electrode                        | glass                     | 4 mm     | 100 $\mu\text{m}$        | <b>A-002009</b> |                 |
|              |   | glass                     | 4 mm     | 25 $\mu\text{m}$         | <b>A-002003</b> |                 |
|              |   | glass                     | 4 mm     | 15 $\mu\text{m}$         | <b>A-002015</b> |                 |
|              |   | glass                     | 4 mm     | 10 $\mu\text{m}$         | <b>A-002005</b> |                 |
| Gold         | gold gauze electrode                            | PEEK                      | 100 mesh | 25 x 35 mm               | <b>A-002251</b> |                 |
|              | long gold electrode                             | PEEK                      | 6 mm     | 3 mm                     | <b>A-012746</b> |                 |
|              | standard gold electrode                         | PEEK                      | 10 mm    | 5 mm                     | <b>A-002418</b> |                 |
|              |   | PEEK                      | 6 mm     | 3 mm                     | <b>A-002421</b> |                 |
|              |   | PEEK                      | 6 mm     | 1.6 mm                   | <b>A-002014</b> |                 |
|              | small gold electrode                            | PEEK                      | 3 mm     | 1.6 mm                   | <b>A-002314</b> |                 |
|              | micro gold electrode                            | glass                     | 4 mm     | 100 $\mu\text{m}$        | <b>A-002010</b> |                 |
|              |   | glass                     | 4 mm     | 25 $\mu\text{m}$         | <b>A-002004</b> |                 |
|              |   | glass                     | 4 mm     | 10 $\mu\text{m}$         | <b>A-002006</b> |                 |
|              | Silver  | standard silver electrode | PEEK     | 10 mm                    | 5 mm            | <b>A-002416</b> |
|              |   | PEEK                      | 6 mm     | 3 mm                     | <b>A-002419</b> |                 |
|              |   | PEEK                      | 6 mm     | 1.6 mm                   | <b>A-002011</b> |                 |
| Palladium    | small silver electrode                          | PEEK                      | 3 mm     | 1.6 mm                   | <b>A-002315</b> |                 |
|              | standard palladium electrode                    | PEEK                      | 6 mm     | 1.6 mm                   | <b>A-002019</b> |                 |
| Nickel       | small palladium electrode                       | PEEK                      | 3 mm     | 1.6 mm                   | <b>A-002319</b> |                 |
|              | standard nickel electrode                       | PEEK                      | 6 mm     | 1.5 mm                   | <b>A-002016</b> |                 |
| Copper       | micro nickel electrode                          | glass                     | 4 mm     | 100 $\mu\text{m}$        | <b>A-002273</b> |                 |
|              | standard copper electrode                       | PEEK                      | 6 mm     | 1.6 mm                   | <b>A-002017</b> |                 |
| Iron         |   | PEEK                      | 6 mm     | 3 mm                     | <b>A-012584</b> |                 |
|              | micro copper electrode                          | glass                     | 4 mm     | 25 $\mu\text{m}$         | <b>A-002271</b> |                 |
| Carbon paste | standard iron electrode                         | PEEK                      | 6 mm     | 1.5 mm                   | <b>A-002018</b> |                 |
|              |   | PEEK                      | 6 mm     | 3 mm                     | <b>A-012585</b> |                 |
|              | standard carbon paste electrode hole depth 4 mm | PEEK                      | 6 mm     | 3 mm                     | <b>A-002210</b> |                 |
|              | small carbon paste electrode hole depth 4 mm    | PEEK                      | 3 mm     | 1.6 mm                   | <b>A-002223</b> |                 |
|              | CPO carbon paste oil base 1 g                   |                           |          |                          | <b>A-001010</b> |                 |

## Small size reference electrodes for aqueous media (Ag/AgCl, Hg)



Reference electrodes are divided into two groups according to the media in which the electrode is immersed (aqueous or organic media).

### Small size reference electrode for aqueous media (Ag/AgCl, Hg)

|   | Junction | Electrolyte  | Purpose   | Catalog n°      |
|---|----------|--|---|-----------------|
| RE-1B silver, silver chloride reference electrode <sup>(1)</sup>                | IPPG*    | 3 M NaCl   | SVC-2, SVC-3, VC-4, bulk electrolysis, RRDE-3A, flat cell | <b>A-012167</b> |
| RE-1S silver, silver chloride reference electrode <sup>(1)</sup>                | IPPG*    | 3 M NaCl   | SECM  | <b>A-012168</b> |
| RE-1CP silver, silver chloride reference electrode                              | ceramic  | saturated KCl  | SVC-2, SVC-3, VC-4, bulk electrolysis, RRDE-3A, flat cell | <b>A-013429</b> |
| RE-3V silver, silver chloride reference electrode screw type, poly acetal resin | IPPG*    | 3 M NaCl   | for flow cell spectroelectrochemical cell SEC-3F          | <b>A-012169</b> |
| RE-3VP Ag/AgCl screw type electrode, PEEK resin                                 | IPPG*    | 3 M NaCl   | for flow cell spectro electrochemical cell SEC-3F         | <b>A-012170</b> |
| Ag/AgCl ink, 2 mL   |          | - surface resistance: 0.2 Ω/sq/25.4 μm<br>- viscosity: 50,000 ± 10,000 cP @21.1 °C<br>- flash point: 82 °C | for micro CV cell, IDA measurement                        | <b>A-011464</b> |
| RE-2BP calomel reference electrode  | ceramic  | saturated KCl  | SVC-2, SVC-3, VC-4, bulk electrolysis, RRDE-3A, flat cell | <b>A-013430</b> |
| RE-2CP Hg <sub>2</sub> SO <sub>4</sub> reference electrode, free from chloride  | ceramic  | saturated K <sub>2</sub> SO <sub>4</sub>   |   | <b>A-013431</b> |
| RE-61AP Hg/HgO reference electrode main body in polyacetal resin                | ceramic  | 1M NaOH  | for alkaline media  | <b>A-013395</b> |
| RHE Reversible Hydrogen Electrode kit   | IPPG*    | HCl or H <sub>2</sub> SO <sub>4</sub>  | for acidic media [pH < 2]                                 | <b>A-013373</b> |

#### Spare parts

|   |                 |
|---|-----------------|
| (1) Set of 10 porous 2.8 mm glass frit (CoralPor™) with PTFE heat shrink (200 mm) | <b>092-VYC3</b> |
|---|-----------------|

#### Options

|  |                 |
|--|-----------------|
| RE-PV preservative vial for reference electrode, 10 mL | <b>A-012108</b> |
| Bridge tube Ø 9.0 mm [2 pieces]                        | <b>A-012177</b> |
| Bridge tube Ø 9.0 mm [22 pieces]                       | <b>A-012307</b> |
| Double junction chamber kit for RHE                    | <b>A-013375</b> |

### Store your reference electrode immersed in the electrolyte

When it is not in use, the recommended way to maintain the reference electrode capability and lifetime is to keep it in a sealed preservative vial with solution. This storage solution should be identical to the filling solution of the reference electrode.


### Prevent contamination

To prevent contamination of the reference electrode, a sample holder can be used.

### Replace the vycor when needed

Yellowish discoloration indicates contamination. This is caused by the adsorption of organic compounds into the pores.

## Small size reference electrodes for non aqueous media (Ag/Ag<sup>+</sup>)



### Small size reference electrode for non aqueous media (Ag/Ag<sup>+</sup>)

|   | Junction | Electrolyte                                 | Purpose | Catalog n°      |
|---|----------|---|---------|-----------------|
| RE-7 non aqueous reference electrode [Ag/Ag <sup>+</sup> ] <sup>(1)</sup>                 | IPPG*    | Ag/Ag <sup>+</sup> /ACN**/TBAP*** CV        |         | <b>A-012171</b> |
| RE-7S non aqueous reference electrode [Ag/Ag <sup>+</sup> ] <sup>(1)</sup>                | IPPG*    | Ag/Ag <sup>+</sup> /ACN**/TBAP*** SECM      |         | <b>A-012172</b> |
| RE-7V non aqueous reference electrode [Ag/Ag <sup>+</sup> ] screw type, poly acetal resin | IPPG*    | Ag/Ag <sup>+</sup> /ACN**/TBAP*** flow cell |         | <b>A-012173</b> |
| RE-7VP non aqueous reference electrode [Ag/Ag <sup>+</sup> ] screw type, PEEK resin       | IPPG*    | Ag/Ag <sup>+</sup> /ACN**/TBAP*** flow cell |         | <b>A-012174</b> |

#### Spare parts


|   |                 |
|---|-----------------|
| Electrolyte solution (10 mL)  | <b>A-012549</b> |
| PTFE cap with Ag wire (for RE-7)  | <b>A-012057</b> |
| PTFE cap with Ag wire (for RE-7S)   | <b>A-012058</b> |
| Sample holder 6 mm diameter (for RE-7) [2 pieces]                                 | <b>A-012176</b> |
| (1) Set of 10 porous 2.8 mm glass frit (CoralPor™) with PTFE heat shrink (200 mm) | <b>092-VYC3</b> |

#### Options

|  |                 |
|--|-----------------|
| RE-PV preservative vial for reference electrode, 10 mL | <b>A-012108</b> |
| Bridge tube Ø 9.0 mm [2 pieces]                        | <b>A-012177</b> |
| Bridge tube Ø 9.0 mm [22 pieces]                       | <b>A-012307</b> |

\* Ion Permeable Porous Glass  
\*\* Acetonitrile  
\*\*\* Tetra Butyl Ammonium Perchlorate

To avoid electrolyte leakage or concentration due to evaporation during storage or transport, the electrolyte is separated from the body. This gives the possibility to fill it regularly.



|         |   |
|---------|---|
| 0.930 V | Hg/HgO/ NaOH (0.1 M)  |
| 0.650 V | Hg/Hg <sub>2</sub> SO <sub>4</sub> / K <sub>2</sub> SO <sub>4</sub> (sat) |
| 0.624 V | Fc/Fc <sup>+</sup> TBAP (0.1M) ACN  |
| 0.542 V | Ag/Ag <sup>+</sup> TBAP (0.1M) ACN  |
| 0.241 V | Hg/Hg <sub>2</sub> Cl <sub>2</sub> KCl (sat)                              |
| 0.236 V | Hg/Hg <sub>2</sub> Cl <sub>2</sub> NaCl (sat)                             |
| 0.205 V | Ag/AgCl/ KCl (3.5 M)  |
| 0.197 V | Ag/AgCl/ KCl (sat)  |
| 0.194 V | Ag/AgCl/ NaCl (sat)   |
| 0.000 V | NHE Normal Hydrogen Electrode   |

E/V vs NHE at 25 °C

Dimensions in mm

## Maintenance of Vycor-based reference electrode

## King size reference electrodes



| King size reference electrodes  | Connection type          | Junction | Length | OD     | Catalog n°       |
|---|--------------------------|----------|--------|--------|------------------|
| Red rod reference electrode   | banana cable             | ceramic  | 103 mm | 7.5 mm | <b>R-REF201</b>  |
|   | banana cable             | ceramic  | 103 mm | 12 mm  | <b>R-REF251</b>  |
|   | screw cap <sup>(1)</sup> | ceramic  | 103 mm | 8 mm   | <b>R-REF421</b>  |
| Hg/HgSO <sub>4</sub> reference electrode in sat. K <sub>2</sub> SO <sub>4</sub> , chloride free                               | banana cable             | ceramic  | 103 mm | 7.5 mm | <b>R-REF601</b>  |
|   | screw cap <sup>(1)</sup> | ceramic  | 103 mm | 8 mm   | <b>R-REF621</b>  |
| Calomel reference electrode non aqueous (LiCl) for general purpose  | screw cap <sup>(1)</sup> | ceramic  | 103 mm | 8 mm   | <b>R-REF921</b>  |
| Calomel reference electrode in KCl sat for CV, straight   | screw cap <sup>(1)</sup> | ceramic  | 120 mm | 8 mm   | <b>R-XR110</b>   |
| Calomel reference electrode in KCl sat for CV, curved   | screw cap <sup>(1)</sup> | ceramic  | 120 mm | 8 mm   | <b>R-XR150</b>   |
| Hg/HgSO <sub>4</sub> reference electrode in sat K <sub>2</sub> SO <sub>4</sub> for CV   | screw cap <sup>(1)</sup> | ceramic  | 120 mm | 8 mm   | <b>R-XR200</b>   |
| Ag/AgCl reference electrode Sat KCl for CV  | screw cap <sup>(1)</sup> | ceramic  | 120 mm | 8 mm   | <b>R-XR300</b>   |
| Hg/HgO reference electrode in 0.1 M KOH for CV  | screw cap <sup>(1)</sup> | ceramic  | 120 mm | 8 mm   | <b>R-XR400</b>   |
| Hg/HgO reference electrode in 1 M KOH with fiber rod  | screw cap <sup>(1)</sup> | ceramic  | 120 mm | 8 mm   | <b>R-XR440</b>   |
| Ag/AgCl reference electrode for EIS   | screw cap <sup>(1)</sup> | ceramic  | 160 mm | 8 mm   | <b>R-XR820</b>   |
| <b>Options</b>  |                          |          |        |        |                  |
| (1) Cable connection for screw cap electrode, recommended to connect these references electrodes to the Bio-Logic instruments | banana plug of 2 mm      | 1 m      |        |        | <b>R-CRI9439</b> |
|   | banana plug of 2 mm      | 100 mm   |        |        | <b>R-A94L111</b> |
|   | banana plug of 4 mm      | 1 m      |        |        | <b>R-CL111</b>   |
| Salt bridge   | ceramic                  | 70 mm    | 12 mm  |        | <b>R-AL100</b>   |
|   | ceramic                  | 138 mm   | 8 mm   |        | <b>R-AL120</b>   |
|   | ceramic                  | 250 mm   | 8 mm   |        | <b>R-AL110</b>   |
| Salt bridge with reverse sleeve, non aqueous  |                          | 138 mm   | 8 mm   |        | <b>R-AL210</b>   |
| Conical rings for 8 mm OD 12 mm electrodes (4 pieces)   |                          |          |        |        | <b>R-X31M012</b> |


## Counter electrodes



| Counter electrodes                               | Length   | Wire Ø     | Surf. area             | Purpose                                     | Catalog n°      |
|--|----------|------------|------------------------|---|-----------------|
| Platinum   | 57 mm    | 0.5 mm     | ~ 0.7 cm <sup>2</sup>  | SVC-2, VC-4, plate material evaluating cell | <b>A-002222</b> |
|  | 50 mm    | 0.5 mm     | ~ 0.7 cm <sup>2</sup>  | SVC-3                                       | <b>A-002233</b> |
|  | 230 mm   | 0.5 mm     | ~ 3.6 cm <sup>2</sup>  | RRDE-3A, bulk electrolysis, SVC-3           | <b>A-002234</b> |
|  | 230 mm   | 0.5 mm     | ~ 3.6 cm <sup>2</sup>  | RRDE-3A, bulk electrolysis, SVC-3           | <b>A-012638</b> |
| Gold   | 230 mm   | 0.5 mm     | ~ 3.6 cm <sup>2</sup>  | RRDE-3A, bulk electrolysis, SVC-3           | <b>A-012639</b> |
| Nickel   | 230 mm   | 0.5 mm     | ~ 3.6 cm <sup>2</sup>  | RRDE-3A, bulk electrolysis, SVC-3           | <b>A-012639</b> |
| Stainless steel                                  | 50 mm    | 1.5 mm     | ~ 2.35 cm <sup>2</sup> | flow cell                                   | <b>A-012198</b> |
| Platinum gauze electrode, PEEK                   | 80 mesh  | 25 x 35 mm | 0.08 mm                | ~ 1 cm <sup>2</sup>                         | <b>A-002250</b> |
| Platinum gauze electrode, platinum wire of 51 mm |          | 25 x 35 mm | 0.08 mm                | ~ 1 cm <sup>2</sup>                         | <b>A-702439</b> |
| Gold gauze electrode, PEEK                       | 100 mesh | 25 x 35 mm | 0.07 mm                | ~ 1 cm <sup>2</sup>                         | <b>A-002251</b> |
|  | 55 mesh  | 40 x 50 mm | 0.12 mm                | ~ 5.1 cm <sup>2</sup>                       | <b>A-002254</b> |

| Metallic electrodes          | Length | OD    | Wire dimension | Catalog n°     |
|------------------------------|--------|-------|----------------|----------------|
| Platinum <sup>(1)</sup>      | 80 mm  | 8 mm  | Ø 1 mm         | <b>R-XM110</b> |
|                              | 120 mm | 8 mm  | plate 5 x 5 mm | <b>R-XM120</b> |
|                              | 120 mm | 8 mm  | plate 8 x 8 mm | <b>R-XM140</b> |
|                              | 120 mm | 12 mm | disk of 10 mm  | <b>R-XM150</b> |
| Glassy carbon <sup>(1)</sup> | 103 mm | 8 mm  | rod Ø of 3 mm  | <b>R-M291C</b> |

| Options   |                             | Catalog n°       |
|---|-----------------------------|------------------|
| (1) Cable connection for screw cap electrode, recommended to connect these references electrodes to the Bio-Logic instruments | 1 m, banana plug of 2 mm    | <b>R-CRI9439</b> |
|   | 100 mm, banana plug of 2 mm | <b>R-A94L111</b> |
|   | 1 m, banana plug of 4 mm    | <b>R-CL111</b>   |
| Conical rings for 8 mm OD 12 mm electrodes (4 pieces)   |                             | <b>R-X31M012</b> |



Dimensions in mm



## Rotating Ring Disk Electrode (RRDE)

RRDE-3A is a system that can be used to perform hydrodynamic (RDE or RRDE) measurement. Its design (short stainless steel shaft) allows the users to control accurately the electrode rotation and the modulation.

RRDE-3A is electronically controlled by a closed loop circuit driving a DC servo-motor. Electrodes are small and rapidly interchangeable.

The unit also provides an adjustable valve system for inert gas purging inside the cell vial. Additionally, the access to glass cell is easy for rinsing, cleaning, and replacing the electrodes. It is easy to remove and replace the cell vial.

RRDE-3A is able to be operated as a stand-alone unit or directly controlled by Bio-Logic electrochemical workstation.

### Features

- Remote- and manual-controlled rotation
- Small sample volume (up to 10 - 15 mL), gas line connector is available, remote and manual controlled purge lines
- Open architecture for easy access to cell

### Specifications

|                        |   |
|------------------------|---|
| Rotational range       | 100 - 8,000 rpm   |
| Setting resolution     | 1 rpm   |
| Accuracy               | < 0.1%  |
| Rotation control type  | PLL (phase-locked loop)                                       |
| Bandwidth              | 60 Hz at 3,500 rpm base and 1,000 rpm peak-to-peak modulation |
| Inlet gas pressure     | 5 PSI   |
| Temperature            | 10-50 °C  |
| Power                  | 100-240 VAC, 50/60 Hz   |
| Dimensions (W x D x H) | 185 x (base: 230, body: 120) x 400 mm                         |
| Weight                 | 6 kg  |



### Rotating Ring Disk Electrode

RRDE-3A apparatus **A-012180**

Working electrodes (page 25), reference electrodes (page 20) and counter electrodes (page 23) have to be purchased separately.

#### Content

|   |                 |
|---|-----------------|
| RRDE-3A glass cell vial 100 mL              | <b>A-012632</b> |
| Spin coating adapter                        | <b>A-012064</b> |
| Nipple for purge valve                      | <b>A-012065</b> |
| Tygon tubing, OD 6.4 mm x ID 3.2 mm (1.3 m) | <b>A-010058</b> |
| PTFE cap for RRDE-3A                        | <b>A-012631</b> |
| Instruction manual                          | -               |
| Power cable                                 | -               |

#### Options

|   |                  |
|---|------------------|
| Water-jacketed glass cell (100 mL, OD 70 mm, ID 46.4 mm, H 80 mm) | <b>A-012652</b>  |
| DB9 cable to control RRDE-3A                                      | <b>092-22/11</b> |
| Bipotentiostat cable for two standard channels                    | <b>092-22/12</b> |
| Bearing assembly  | <b>A-012625</b>  |

## Working electrode

### Standard tips

|                                    |                                       | Tips length | Tips OD | Ring ID | Ring OD | Disk Ø          | Catalog n°      |
|------------------------------------|---------------------------------------|-------------|---------|---------|---------|-----------------|-----------------|
| RRDE                               | platinum ring/platinum disk           | 25 mm       | 12 mm   | 5.0 mm  | 7.0 mm  | 4.0 mm          | <b>A-011172</b> |
|                                    | platinum ring/glassy carbon disk      | 25 mm       | 12 mm   | 5.0 mm  | 7.0 mm  | 4.0 mm          | <b>A-011162</b> |
|                                    | gold ring/glassy carbon disk          | 25 mm       | 12 mm   | 5.0 mm  | 7.0 mm  | 4.0 mm          | <b>A-011163</b> |
|                                    | platinum ring/gold disk               | 25 mm       | 12 mm   | 5.0 mm  | 7.0 mm  | 4.0 mm          | <b>A-011164</b> |
|                                    | gold ring/platinum disk               | 25 mm       | 12 mm   | 5.0 mm  | 7.0 mm  | 4.0 mm          | <b>A-012617</b> |
| RDE                                | glassy carbon ring/glassy carbon disk | 25 mm       | 12 mm   | 5.0 mm  | 7.0 mm  | 4.0 mm          | <b>A-012618</b> |
|                                    | glassy carbon disk                    | 25 mm       | 12 mm   | -       | -       | 3.0 mm          | <b>A-011169</b> |
|                                    | platinum disk                         | 25 mm       | 12 mm   | -       | -       | 3.0 mm          | <b>A-011170</b> |
|                                    | gold disk                             | 25 mm       | 12 mm   | -       | -       | 3.0 mm          | <b>A-011171</b> |
|                                    | aluminum disk                         | 25 mm       | 12 mm   | -       | -       | 3.0 mm          | <b>A-011966</b> |
|                                    | silver disk                           | 25 mm       | 12 mm   | -       | -       | 3.0 mm          | <b>A-011967</b> |
|                                    | copper disk                           | 25 mm       | 12 mm   | -       | -       | 3.0 mm          | <b>A-011968</b> |
|                                    | nickel disk                           | 25 mm       | 12 mm   | -       | -       | 3.0 mm          | <b>A-011969</b> |
|                                    | tantalum disk                         | 25 mm       | 12 mm   | -       | -       | 3.0 mm          | <b>A-011970</b> |
|                                    | titanium disk                         | 25 mm       | 12 mm   | -       | -       | 3.0 mm          | <b>A-011971</b> |
|                                    | tungsten disk                         | 25 mm       | 12 mm   | -       | -       | 3.0 mm          | <b>A-011972</b> |
| carbon paste disk, hole depth 4 mm | 25 mm                                 | 12 mm       | -       | -       | 3.0 mm  | <b>A-011973</b> |                 |

### Option

|                    |                 |
|--------------------|-----------------|
| PK-3 polishing kit | <b>A-011975</b> |
|--------------------|-----------------|

### Disk Replaceable Electrode tips

|                        |   | Catalog n°      |
|------------------------|---|-----------------|
| RRDE                   | platinum ring/GC disk replaceable electrode kit | <b>A-013336</b> |
|                        | <b>content</b> platinum ring assembly           | <b>A-013337</b> |
|                        | glassy carbon disk                              | <b>A-013338</b> |
| RDE                    | PTFE spacer (3 pieces)                          | <b>A-013339</b> |
|                        | glassy carbon disk replaceable electrode kit    | <b>A-013362</b> |
|                        | <b>content</b> glassy carbon disk               | <b>A-013338</b> |
|                        | PTFE spacer (3 pieces)                          | <b>A-013339</b> |
|                        | disk assembly                                   | <b>A-013361</b> |
|                        | gold disk replaceable electrode kit             | <b>A-013364</b> |
|                        | <b>content</b> gold disk                        | <b>A-013366</b> |
|                        | PTFE spacer (3 pieces)                          | <b>A-013339</b> |
|                        | disk assembly                                   | <b>A-013361</b> |
|                        | platinum disk replaceable electrode kit         | <b>A-013365</b> |
|                        | <b>content</b> platinum disk                    | <b>A-013367</b> |
| PTFE spacer (3 pieces) | <b>A-013339</b>                                 |                 |
| disk assembly          | <b>A-013361</b>                                 |                 |

### Tool kit

|  |                 |
|--|-----------------|
| Disk electrode polishing and exchanging tool kit | <b>A-013340</b> |
|--|-----------------|



RRDE platinum ring/GC disk replaceable electrode kit



Disk electrode polishing and exchanging tool kit

## Rotating Disk Electrode (RDE)

The EDI101 is a versatile and rugged rotating disk electrode ideal for use with any Bio-Logic potentiostat/galvanostat. It is available with a wide choice of quick-fit exchangeable tips in standard or custom materials: platinum, glassy carbon, gold... Thanks to its conic shaft holder, the EDI101 can be used with a wide range of cells.


The CTV101 speed control unit controls the EDI101 rotation speed. It offers accuracy better than ±0.1% for precise and fully reproducible experimental conditions. The measured rotation speed is clearly displayed on a four-digit LCD. The speed control unit can be directly driven from a Bio-Logic potentiostat/galvanostat. You can program the variation of the rotation speed from one method to another method and achieve a "Levich" type experiment all in one go.

This RDE can be used with the EL-ELECTRO-1, EL-ELECTRO-2 or EL-ELECTRO-3 kit.

### Specifications

|                                     |   |
|-------------------------------------|---|
| Rotational range                    | 100 - 5,000 rpm   |
| Accuracy                            | < 0.1%  |
| Eccentricity                        | < ±0.1 mm   |
| Setting resolution                  | 1 rpm   |
| Material of EDI101                  | Kel-F   |
| Gas inlet for brush contact purging | Ø 2 mm  |
| Temperature                         | 5 - 40 °C   |
| Power                               | 230 Vac, 47.5 - 63 Hz   |
| Consumption                         | 12 VA   |
| Dimensions                          | CTV101: 80 x 240 x 230 mm (W x D x H),<br>EDI101: 255 mm (length) |
| Weight                              | CTV101: 1.5 kg<br>EDI101: 0.3 kg                                  |

| Rotating Disk Electrode   |                          | Catalog n°        |
|---|--------------------------|-------------------|
| Rotating Disk Electrode system including control unit and electrode rotator |                          | <b>R-RDE/EDI</b>  |
| <b>Content</b>  |                          |                   |
| CTV101, controlling unit for EDI101   |                          | <b>R-R21V035</b>  |
| EDI101 RDE without tip (10-5,000 tr/min)                                    |                          | <b>R-A35T040</b>  |
| BNC-BNC connector 1 m   |                          | <b>COR28100</b>   |
| <b>Options</b>  |                          |                   |
| BEC/EDI electrochemical cell with glassware and electrodes                  |                          | <b>R-XR51V001</b> |
| EDI tip   | glassy carbon Ø 3 mm     | <b>R-A35T090</b>  |
|   | glassy carbon Ø 5 mm     | <b>R-A35T095</b>  |
|   | platinum Ø 2 mm          | <b>R-A35T100</b>  |
|   | platinum Ø 5 mm          | <b>R-A35T105</b>  |
|   | gold Ø 2 mm              | <b>R-A35T110</b>  |
|   | gold Ø 5 mm              | <b>R-A35T120</b>  |
|   | copper Ø 5 mm            | <b>R-A35T130</b>  |
|   | graphite Ø 4 mm          | <b>R-A35T140</b>  |
|   | silver Ø 5 mm            | <b>R-A35T150</b>  |
|   | nickel Ø 5 mm            | <b>R-A35T160</b>  |
|   | without any material     | <b>R-A35T402</b>  |
|   | stainless steel Ø 5 mm   | <b>R-A35T420</b>  |
|   | sample holder Ø 11 mm    | <b>R-A35T450</b>  |
|   | platinum, L 40 mm Ø 5 mm | <b>R-A35T452</b>  |
|   | gold, L 40 mm Ø 4 mm     | <b>R-A35T456</b>  |
| Adapter EDI (banana 4 mm female, banana 2 mm male)                          |                          | <b>R-X41V001</b>  |
| DB9 cable to control RDE  |                          | <b>092-22/11</b>  |



## RDE

A 20<sup>th</sup> century electro-chemist, Dr. Levich proposed the relationship between the voltammetry limiting current ( $I_L$ , A) of a rotating disk electrode and the physical parameters as follows:

$$I_L = \pm 0.62 n F A D^{2/3} \omega^{1/2} \nu^{-1/6} C$$


Where:

- $\omega$  is the rotational velocity (radian per second),
- $n$  electrons involved within the reaction (dimensionless),
- $F$  is the Faraday constant ( $F = 96,485 \text{ C/mol}$ ),
- $A$  is the disk area ( $\text{cm}^2$ ),
- $D$  is the diffusion coefficient ( $\text{cm}^2/\text{s}$ ),
- $\nu$  is the kinematic viscosity ( $\text{cm}^2/\text{s}$ ),
- $C$  is the electroactive species concentration ( $\text{mol}/\text{cm}^3$ ).

$$D = \left( \frac{m}{0.62 n F A \nu^{-1/6} C} \right)^{3/2}$$


An example of the Levich analysis is performed on a 5.5 mM Potassium Hexacyanoferrate (III), in KCl (0.1 M) solution. The working electrode is a glassy carbon disk electrode.

The sequence comprised of 13 rotational velocities from 10 rpm to 2000 rpm, and 13 twin-cycle cyclic voltammograms. The resulting CVs are overlaid and presented in the image below. The Levich analysis leads to the diffusion coefficient ( $D$ ) determination which is  $D = 6.36 \times 10^{-6} \text{ cm}^2/\text{s}$ .



## RRDE

For RRDE measurement, a bipotentiostat is needed. A bipotentiostat controls the two working electrodes i.e. one channel to control the disk electrode and the other to control the ring. An SP-300 equipped with two channels would be an appropriate instrument for RRDE applications.



Because of the presence of two working electrodes in the same setup, a specific connection mode (to avoid any ground loop trouble) is needed.


This can be achieved in two different ways:


- isolating the two channels (at least one potentiostat in floating mode),
- grounding the counter electrode.

The "CE-to-ground" mode of the Bio-Logic instruments provides this unique capability.

Of these two options, the latter is preferred because the leakage current is less. All of the Bio-Logic multichannel potentiostats (except cyclers) offer such electrode connection type.

At the disk electrode, the electro-active species are oxidized or reduced according to the applied potential. This new species is detected by reduction or oxidation at the ring, respectively. In a typical experiment, a CV is performed on the disk electrode and a constant voltage is applied on the ring electrode. It is the CV-CA technique that is available in EC-Lab® in the "bipotentiostat" techniques folder.





## Spectrometer

Spectroelectrochemistry (SEC) can be useful to elucidate electrochemical reaction mechanisms. The spectroelectrochemical kit is composed of three parts (spectrometer, light source and cuvette holder).

The spectrometer is equipped with a trigger to synchronize electrochemical and spectroscopic measurements.

### Utilization modes

#### Transmittance

- Absorbance/transmittance**
- Concentration of chemicals (solution)
  - Polymer extrusion processes
  - DNA quantification

- Reflectance**
- Freshness testing
  - Film thickness/composition (quality control)
  - Activation energy of photocatalytic species
  - Textile quality control

#### Fluorescence


- Fluorescence**
- Marine organisms
  - Biology (DNA, protein, cell proliferation assay, histamine-analysis, alga monitoring)
  - Environmental fields (waste water analysis, ground water trace studies, hydrocarbon detection, dissolved oxygen)
  - Plant efficiency (plant physiology, plant breeding, horticulture, agronomy, agrochemicals, pollution studies, forestry, ecology)
  - Tissue diagnosis

- Scattering**
- Oil concentrations of oil/water system
  - Raman spectroscopy
  - Physical transition phenomena (e.g. melting point, glasstransition crystallize temperature)

#### Irradiance

- Emission**
- Astronomy (e.g. spectra of Hale-Bopp, plasma monitoring)
  - In situ metal monitoring
  - Luminescence (PL,EL), LED & laser wavelength

|                                    | <b>SEC2000-UV/VIS</b>   | Catalog n°      | <b>SEC2000-VIS/NIR</b>           | Catalog n°       |
|------------------------------------|---|-----------------|----------------------------------|------------------|
| <b>Spectrometer combo kit type</b> | SEC2000-UV/VIS  | <b>A-012838</b> | SEC2000-VIS/NIR                  | <b>A-012973</b>  |
| <b>Content</b>                     |   |                 |                                  |                  |
| Spectrometer                       | SEC2000-UV/VIS  | <b>A-012856</b> | SEC2000-VIS/NIR                  | <b>A-012971</b>  |
| Light source                       | SEC2000-DH UV   | <b>A-012193</b> | SEC2000-TH visible               | <b>A-012194</b>  |
| Cuvette holder                     | SEC2000-CUV   | <b>A-012195</b> | SEC2000-CUV with diffusion plate | <b>A-012244</b>  |
| <b>Option</b>                      | Connecting cable to synchronize the spectrometer and the Bio-Logic potentiostat |                 |                                  | <b>092-22/11</b> |



## Spectrometer

|                        | <b>SEC2000-UV/VIS combo kit</b>                               | <b>SEC2000-VIS/NIR combo kit</b>                              |
|------------------------|---|---|
| Detector               | 2,048 pixels CCD array  | 2,048 pixels CCD array  |
| Wave length range      | 220 - 800 nm  | 500 - 1,000 nm  |
| Grating                | blazed at 400 nm  | blazed at 750 nm  |
| Resolution             | 1.8 ± 0.2 nm: slit 50 µm x 1,000 µm (standard) <sup>(1)</sup> | 1.8 ± 0.2 nm: slit 50 µm x 1,000 µm (standard) <sup>(1)</sup> |
| Accuracy               | < 1% at a Abs   | < 1% at a Abs   |
| Dark noise             | < 2 mAbs  | < 2 mAbs  |
| A/D resolution         | 14 bit  | 14 bit  |
| Optical entrance       | SMA905  | SMA905  |
| Interface              | USB 2.0   | USB 2.0   |
| Operating system       | Windows™ XP   | Windows™ XP   |
| Dimensions (W x D x H) | 98 x 118 x 35 mm  | 98 x 118 x 35 mm  |
| Catalog n°             | <b>A-012856</b>   | <b>A-012971</b>   |

(1). Slit could be selected from: 10, 25, 100, 200 µm.

## Light source

|                               | <b>SEC2000-DH Light source (deuterium &amp; tungsten halogen)</b> | <b>SEC2000-TH Light source (tungsten halogen)</b> |
|-------------------------------|---|---|
| Spectral range (Halogen + D2) | 200 - 1,100 nm  | 360 - 2,000 nm                                    |
| D2 lamp spectral range        | 200 - 400 nm  | -   |
| power consumption (240 nm)    | < 5 x 10 <sup>-8</sup> W/nmsr                                     | -   |
| stability                     | 1 x 10 <sup>-3</sup> AU   | -   |
| drift                         | < 0.25%/h   | -   |
| bulb life <sup>(2)</sup>      | > 1,000 h   | -   |
| Halogen lamp bulb life        | > 2,000 h   | 1,500 h   |
| Connector                     | SMA905 included   | SMA905 not included                               |
| Dimensions (W x D x H)        | 98 x 118 x 35 mm  | 98 x 118 x 35 mm                                  |
| Catalog n° light source       | <b>A-012193</b>   | <b>A-012194</b>                                   |
| <b>Spare parts</b>            |   |   |
| Bulb                          | <b>A-012221</b> (UV bulb)   | <b>A-012222</b> (TV bulb)                         |

(2). Value for less than 50% power consumption of 240 nm.



## Cuvette holder

|            | <b>SEC2000-UV/VIS combo kit</b> | <b>SEC2000-VIS/NIR combo kit</b> |
|------------|---------------------------------|----------------------------------|
| Catalog n° | <b>A-012195</b>                 | <b>A-012244</b>                  |

## Static experiment




| Quartz glass spectroelectrochemical cell kit | Platinum             | Catalog n°      | Gold                  | Catalog n°      |
|--|----------------------|-----------------|-----------------------|-----------------|
| <b>0.5 mm</b>                                |                      | <b>A-012813</b> |                       | <b>A-012814</b> |
| <b>Content</b>                               |                      |                 |                       |                 |
| platinum counter electrode                   |                      | <b>A-012609</b> |                       | <b>A-012609</b> |
| thin layer quartz glass cell                 |                      | <b>A-012815</b> |                       | <b>A-012815</b> |
| PTFE cap                                     |                      | <b>A-011501</b> |                       | <b>A-011501</b> |
| purging tube (ETFE, 100 mm)                  |                      | -               |                       | -               |
| gauze working electrode                      | 80 mesh, height 5 mm | <b>A-012606</b> | 100 mesh, height 5 mm | <b>A-012607</b> |
|  |                      | <b>A-012904</b> |                       | <b>A-012905</b> |
| <b>Content</b>                               |                      |                 |                       |                 |
| platinum counter electrode                   |                      | <b>A-012906</b> |                       | <b>A-012906</b> |
| thin layer quartz glass cell                 |                      | <b>A-012907</b> |                       | <b>A-012907</b> |
| PTFE cap                                     |                      | <b>A-011501</b> |                       | <b>A-011501</b> |
| purging tube (ETFE, 100 mm)                  |                      | -               |                       | -               |
| gauze working electrode                      | 80 mesh, height 6 mm | <b>A-011498</b> | 100 mesh, height 6 mm | <b>A-012017</b> |
| <b>Options</b>                               |                      |                 |                       |                 |
| RE-1B Ag/AgCl reference electrode            |                      |                 |                       | <b>A-012167</b> |
| RE-7 non aqueous reference electrode         |                      |                 |                       | <b>A-012171</b> |
| Purging tube (ETFE), 1 m                     |                      |                 |                       | <b>A-010537</b> |



## Flow experiment

| Flow experiment                                |  | Quantity | Catalog n°      |
|--|--|----------|-----------------|
| <b>SEC-2F spectroelectrochemical flow cell</b> |  |          |                 |
| <b>Content</b>                                 |  |          |                 |
| SEC-2F flow cell                               | base                                       | 1        | -               |
|  | cover                                      | 1        | -               |
|  | block A                                    | 1        | -               |
|  | block B                                    | 1        | -               |
| SEC-2F S500 silicone gasket                    |  | 2        | <b>A-012661</b> |
| Stainless tube OD 1.59 mm (length of 50 mm)    |  | 1        | <b>A-012198</b> |
| Needle adaptor                                 |  | 1        | -               |
| Dynaseal PEEK fingertight                      |  | 2        | -               |
| Silicon tube (300 mm)                          |  | 1        | -               |
| PTFE tube (1 m)                                |  | 1        | -               |
| <b>Options</b>                                 |  |          |                 |
| Reference Ø 10 x 55 mm electrode screw type    | RE-3V aqueous                              |          | <b>A-012169</b> |
|  | RE-3VP aqueous (PEEK)                      |          | <b>A-012170</b> |
|  | RE-7V non aqueous                          |          | <b>A-012173</b> |
|  | RE-7VP non aqueous (PEEK)                  |          | <b>A-012174</b> |
| Working grid electrode for flow cell           | 8 x 27 x 1 mm platinum (1 piece)           |          | <b>A-012655</b> |
|  | gold (1 piece)                             |          | <b>A-012656</b> |
|  | carbon grid electrode (1 piece)            |          | <b>A-012657</b> |
|  | ITO electrode (4 pieces)                   |          | <b>A-012658</b> |
|  | 8 x 27 x 0.5 mm ITO electrode (12 pieces)  |          | <b>A-011465</b> |
|  | 10 x 20 x 0.5 mm ITO electrode (10 pieces) |          | <b>A-010887</b> |
|  | 10 x 10 x 0.5 mm ITO electrode (30 pieces) |          | <b>A-011233</b> |
|  | Ø 4 inch x 0.5 mm ITO electrode (1 piece)  |          | <b>A-011827</b> |
| Gasket   | silicone S500, 500 µm thick (4 pieces)     |          | <b>A-012661</b> |
|  | PTFE T500, 500 µm thick (4 pieces)         |          | <b>A-012664</b> |
|  | T250, 250 µm thick (4 pieces)              |          | <b>A-012665</b> |
|  | T100, 100 µm thick (4 pieces)              |          | <b>A-012666</b> |
| Fibre and lens                                 | 400 m optical fibre SR, 250 mm             |          | <b>A-012667</b> |
|  | 400 µm optical fibre SR, 2 m               |          | <b>A-011522</b> |
|  | UV/VIS collimating lens, 200-2,000 nm      |          | <b>A-012234</b> |



Gold grid electrode **A-012656**

### Solution volumes

| Gasket of |  | Volume  |
|-----------|--|---------|
| 100 µm    |  | 4.6 µl  |
| 250 µm    |  | 11.5 µl |
| 500 µm    |  | 230 µl  |



Thin layer  
spectroelectrochemical  
cell

The UFS set is conceived to perform thin-layer spectroelectrochemistry.




UF-Full

## Spacer




Pseudo-reference electrode  
Working electrode  
Auxiliary electrode



Full kit +  
PTFE adapter for optical fibres + optical fibre

| Thin layer spectroelectrochemical cell kit   | Platinum                   | Catalog n°        | Gold                   | Catalog n°        |
|--|----------------------------|-------------------|------------------------|-------------------|
| <b>Full kit</b>  |                            | <b>UF-Full-PT</b> |                        | <b>UF-Full-AU</b> |
| <b>Content</b>   |                            |                   |                        |                   |
| 1 PTFE windows aligner + 1 spacer (platinum auxiliary electrode/Ag pseudo-reference electrode)                         | platinum working electrode | <b>UF-SEC-PT</b>  | gold working electrode | <b>UF-SEC-AU</b>  |
| Stainless steel body cell  |                            | -                 |                        | -                 |
| XY holder (with magnets) to universally fit the cell in a spectrophotometer  |                            | <b>UF-XY</b>      |                        | <b>UF-XY</b>      |
| Plug   |                            | <b>UF-P</b>       |                        | <b>UF-P</b>       |
| PTFE mask  |                            | <b>UF-M</b>       |                        | <b>UF-M</b>       |
| <b>Options</b>   |                            |                   |                        |                   |
| PTFE adapter for optical fibres with UV/VIS/NIR collimating lenses (200-2500 nm), adjustable focus, SMA-905 connection |                            |                   |                        | <b>UF-OFA</b>     |
| Base stand for optical fibres measurements   |                            |                   |                        | <b>UF-BS</b>      |
| PTFE cuvette holder to be used with UF-OFA and UF-BS   |                            |                   |                        | <b>UF-C</b>       |
| 1 spacer (platinum working electrode/platinum auxiliary electrode/Ag pseudo-reference electrode)                       |                            |                   |                        | <b>UF-SPP</b>     |
| 1 spacer (gold working electrode/platinum auxiliary electrode/Ag pseudo-reference electrode)                           |                            |                   |                        | <b>UF-SPA</b>     |
| Optical fibre: SMA-SMA connector length 1 m  |                            |                   |                        | <b>092-101</b>    |




## Finite diffusion conditions

In such finite diffusion conditions, a thin solution layer ( $\leq 0.2$  mm) adjacent to the electrode is confined by the cell walls, so that the cell thickness is smaller than the diffusion layer and the mass transfer can be ignored. The most significant virtue of thin-layer cells is the absence of the effect of the diffusion process and the rapidity with which the electro-active species can be exhaustively electrolyzed. The dropping to near zero of the current flow following the peak in the current potential plot is a characteristic behavior of thin-layer cells, indicating exhaustive electrolysis of the cell reactant and minimal diffusion effects in thin-layer electrochemical cells<sup>(1) (2)</sup>.

Small potential sweep rates (2-10 mV/s) are necessary both to ensure homogeneity of the reactant/product concentrations in the cell and to control the resistive effects.

A cyclic voltammetry should be registered in your SEC cell to better localize the redox process of interest: in fact, the UF spacer has a silver pseudo-reference electrode, which is sensitive to the solution medium (but is expected to remain constant in time in each given experimental condition). Potential drifts can be observed in the presence of irreversible redox reactions, which may alter the solution/analyte composition.




To maintain the ohmic drop as low as possible, the current should also be kept low which means low scan rates and low sample concentrations may be preferred. However, in some cases, due to the narrowness of the optical path, the use of a relatively high concentration of the sample may be required to study the changes of bands with a low molar extinction coefficient.

A milli-molar concentration appears to be the most optimized concentration, but in some cases, the concentration can be adapted.

Due to the high Infra-Red absorbance of more common solvents and electrolytes, a very carefully measured background should be obtained before each Infra-Red spectroelectrochemical experiment.

For the same reason, it is also important to avoid changing the tightness of the cell screws during the experiment itself to avoid changes of the optical path.



(1): M. Krejčík, M. Daněš and F. Hartl, J. Electroanal. Chem., 1991, 317, 179.

(2): P. Leoni, F. Marchetti, C. Bonaccorsi, F. Fabrizi de Biani, L. Marchetti, P. Zanello, Chem. Eur. J., 2008, 14, 847



092-QCA-FC

SE-CL3

## EQCM


The instrument monitors both the resonant frequency and the resonant resistance which are also reflected on the two analog outputs.

It is possible to connect the Quartz Crystal Analyzer (QCA922) to a Bio-Logic potentiostat using a DB9-8BNC cable (catalog n°: 092-22/1).

| EQCM  |                                      | Catalog n°           |
|---|--------------------------------------|----------------------|
| Quartz crystal microbalance 27 MHz kit            |                                      | <b>SE-QCA922A</b>    |
| <b>Content</b>                                    |                                      |                      |
| EQCM 27 MHz main unit and cable                   |                                      | <b>SE-QCA922A-00</b> |
| Male BNC/BNC cable (2 pieces, length 1 m)         |                                      | <b>COR28100</b>      |
| <b>Options</b>                                    |                                      |                      |
| Holder is needed to get a full QCM or EQCM set-up | static and flow QCA cell             | <b>092-QCA-FC</b>    |
|   | dip cell                             | <b>SE-CL3</b>        |
|   | well cell (PTFE)                     | <b>SE-CL4</b>        |
|   | well cell (PEEK)                     | <b>SE-CL4PK</b>      |
|   | transparent well cell                | <b>SE-CL5</b>        |
|   | flow cell (90 µL) (PTFE)             | <b>SE-CL6</b>        |
|   | flow cell (90 µL) transparent (PTFE) | <b>SE-CL7</b>        |
|   | flow cell (90 µL) (PEEK)             | <b>SE-CL6PK</b>      |
| Connector from QCA to potentiostat                |                                      | <b>092-22/1</b>      |
| Low flow peristaltic pump                         |                                      | <b>EL-AV-008</b>     |

## EQCM specifications

|                           |   |
|---------------------------|---|
| Frequency range           | 1 MHz-30 MHz, resolution 0.1 Hz                                     |
| Resonant resistance range | 10-16 Ω, resolution 0.1 Ω   |
| ΔF output                 | full scale: ±10 V (12 bit)<br>±200 Hz / ±2 kHz / ±20 kHz / ±200 kHz |
| Resistance output         | full scale: 0-10 V (11 bit)<br>1 / 2 / 4 / 8 / 16 kΩ                |
| Gate time                 | variable (0.1, 1.0, 10.0 s)   |
| Display                   | 40 digits x 2 rows  |
| Interface                 | IEEE-488, RS-232C, USB  |
| Dimensions (W x D x H)    | 260 x 230 x 88 mm   |
| Weight                    | 3.3 kg  |
| Temperature               | 0 - 40 °C   |
| Power consumption         | 100-120 V AC, 230-240 V AC, 50-60 Hz, (15 W)                        |



## Quartz resonators

| Quartz resonator (calibration at 5 MHz at factory with no charge when purchased) |  | Catalog n°        |
|--|--|-------------------|
| 5 MHz  | gold electrode (25 pieces)                     | <b>SE-5AU</b>     |
| 9 MHz  | etching finish                                 | <b>SE-9AU-E/1</b> |
|  | gold electrode (100 pieces)                    | <b>SE-9AU-E/3</b> |
| polish finish  | gold electrode (30 pieces)                     | <b>SE-9AU-P/1</b> |
|  | gold electrode (100 pieces)                    | <b>SE-9AU-P/3</b> |
| standard finish  | resonator and lead wire                        | <b>SE-9AL</b>     |
|  | aluminum electrode (25 pieces)                 | <b>SE-9AU</b>     |
|  | gold electrode (25 pieces)                     | <b>SE-9C</b>      |
|  | graphite electrode (25 pieces)                 | <b>SE-9CU</b>     |
|  | copper electrode (25 pieces)                   | <b>SE-9MO</b>     |
|  | molybdenum electrode (25 pieces)               | <b>SE-9NI</b>     |
|  | nickel electrode (25 pieces)                   | <b>SE-9PT</b>     |
|  | platinum electrode (25 pieces)                 | <b>SE-9SS</b>     |
|  | stainless steel (SUS304) electrode (25 pieces) | <b>SE-9AU-S</b>   |
|  | gold electrode (25 pieces)                     | <b>SE-9PT-S</b>   |
| separated lead wire  | platinum electrode (25 pieces)                 | <b>SE-9AL-M</b>   |
| mirror finish  | resonator and lead wire                        | <b>SE-9AU-M</b>   |
|  | aluminum electrode (25 pieces)                 | <b>SE-9AU-M2</b>  |
|  | gold electrode (25 pieces)                     | <b>SE-9C-M</b>    |
|  | gold electrode (500 pieces)                    | <b>SE-9CU-M</b>   |
|  | graphite electrode (25 pieces)                 | <b>SE-9IT-M</b>   |
|  | copper electrode (25 pieces)                   | <b>SE-9NI-M</b>   |
|  | ITO electrode (25 pieces)                      | <b>SE-9PT-M</b>   |
|  | nickel electrode (25 pieces)                   | <b>SE-9SI-M</b>   |
|  | platinum electrode (25 pieces)                 | <b>SE-9SS-M</b>   |
|  | silicon electrode (25 pieces)                  | <b>SE-9TI-M</b>   |
| stainless steel (SUS304) electrode (25 pieces)                                   | <b>SE-9AU-MS</b>                               |                   |
| titanium electrode (25 pieces)   | <b>SE-9IT-MS</b>                               |                   |
| separated lead wire  | gold electrode (25 pieces)                     | <b>SE-9PT-MS</b>  |
|  | ITO electrode (25 pieces)                      | <b>SE-9IT-MS</b>  |
|  | platinum electrode (25 pieces)                 | <b>SE-9PT-MS</b>  |
| <b>Others</b>  |  |                   |
| Resonators lead wire no sputter (50 pieces)                                      |  | <b>SE-LEAD/2</b>  |
| Resonators 9 MHz, no sputter   |  | <b>SE-90/2</b>    |
|  | no lead wire (50 pieces)                       | <b>SE-90-M/2</b>  |
|  | mirror finish no lead wire (50 pieces)         | <b>SE-90-M/2</b>  |

Gold electrode

Platinum electrode

Stainless steel electrode

## Surface finishing

|                 | Roughness | Electrode materials deposition |
|-----------------|-----------|--------------------------------|
| Standard finish | 0.6 µm    | sputtered                      |
| Etching finish  | 0.6 µm    | vacuum deposition              |
| Mirror finish   | 0.06 µm   | sputtered                      |
| Polish finish   | 0.06 µm   | vacuum deposition              |



## Biosensor QCM

The specification (27 MHz) of the Quartz Crystal Microbalance (QCM934) allows users to measure a small mass change. QCM934 is especially suitable in the bio-sensing field.

- Excellent sensitivity and resolution.
- Specially designed oscillation circuit for added stability.
- Simultaneous measurement of resonant frequency, resonant resistance and temperature.
- Expandable up to 4 channels.
- Works with flows cell, well cell, dip cell.
- Computer control and data collection via USB interface.

Resonators of 9 MHz are compatible with QCM934 and the cell (see page 35).

## Specifications


|                           |  |
|---------------------------|--|
| Frequency range           | 25 MHz-27 MHz, resolution 0.01 Hz            |
| Resonant resistance range | 2 Ω, resolution 0.1 Ω                        |
| Max channel               | 4  |
| Gate time                 | variable (0.1, 1.0, 10.0 s)                  |
| Display                   | 20 digits x 2 rows                           |
| Interface                 | USB 2.0                                      |
| Dimensions (W x D x H)    | 320 x 280 x 133 mm                           |
| Weight                    | 5.5 kg                                       |
| Temperature               | 10 - 40 °C                                   |
| Power consumption         | 100-120 V AC, 230-240 V AC, 50-60 Hz, (60 W) |

| Flow mode   | Catalog n°            |
|---|-----------------------|
| Flow cell kit                                     | <b>SE-QCM934-F</b>    |
| <b>Content</b>                                    |                       |
| Main unit   | <b>SE-QCM934-000</b>  |
| WinQCM software                                   | <b>SE-WQCM</b>        |
| Oscillator & temperature control module           | <b>SE-QCM934-200</b>  |
| Cable for oscillator & temperature control module | <b>SE-QCM934-200A</b> |
| Measurement chamber                               | <b>SE-QCM934-500</b>  |
| Flow cell (90 µL)                                 | <b>SE-QCM934-510</b>  |
| <b>Options</b>                                    |                       |
| Additional channel                                | <b>SE-QCM934-1CHF</b> |
| Low flow peristaltic pump                         | <b>EL-AV-008</b>      |



"Flow cell-based system of QCM934", "perista pump" and "notebook computer"

| Batch mode  | Catalog n°            |
|---|-----------------------|
| Batch cell kit                                    | <b>SE-QCM934-B</b>    |
| <b>Content</b>                                    |                       |
| Main unit   | <b>SE-QCM934-000</b>  |
| WinQCM software                                   | <b>SE-WQCM</b>        |
| Oscillator & temperature control module           | <b>SE-QCM934-200</b>  |
| Cable for oscillator & temperature control module | <b>SE-QCM934-200A</b> |
| Cable for OSC module connection                   | <b>SE-QCM934-300</b>  |



"Batch cell-based system of QCM934", "transparent well/dip cell" and "notebook computer"

# Quartz crystal analyzer applications

The Quartz crystal microbalance is a mass-sensitive detector based on frequency changes of an oscillating quartz crystal. The oscillation frequency of the crystal is proportional to the mass of the crystal, as well as solution properties near the surface (including viscosity, density, temperature, and compression waves).

A mass increase results in a frequency decrease. Sauerbrey was the first who provided a description and experimental verification of the mass/frequency relationship between foreign layers firmly attached to the quartz crystal resonator<sup>[1][2]</sup>.

The Sauerbrey equation is defined as:

$$\Delta f = \frac{2f_0^2}{A\sqrt{\rho_q\mu_q}} \Delta m$$

$f_0$ : resonant frequency (Hz)

$\Delta f$ : frequency change (Hz)

$\Delta m$ : mass change (g)


$A$ : piezoelectrically active crystal area (area between electrodes, cm<sup>2</sup>)

$\rho_q$ : density of quartz ( $\rho_q = 2.648 \text{ g/cm}^3$ )

$\mu_q$ : shear modulus of quartz for AT-cut crystal ( $\mu_q = 2.947 \times 10^{11} \text{ g/(cm}\cdot\text{s}^2)$ )

## Electropolymerization of pyrrol

The polypyrrol film was coated on an Au quartz (used as the working electrode) using cyclic voltammetry (20 cycles).



The quartz electrode was immersed in an acetonitrile solution (Bu<sub>4</sub>NPF<sub>6</sub> 0.2 mol.L) containing a solution of 1 methylpyrrol monomer (0.01 mol.L).


## Polypyrrol film growth on the quartz working electrode

Fig. 1 presents the polypyrrol film growth on the quartz electrode during the successive cycles of cyclic voltammetry. The reversibility of the charge transfer in such a polymer film is often dependent on the deposition mode (quasi-reversible in this example). The growth is very regular but tends to slow down in the last cycles. That can be due to an interfacial depletion of the solution in methyl pyrrol species in the layer close to the electrode surface and to a saturation of the working electrode surface area.

## QCM measurements during the film growth


Fig. 1 shows the resonant frequency decrease and the resonant resistance increase while the polymer film is growing. Moreover the variation is dependent on the potential sweep resulting in a pseudo oscillation of frequency and resistance related to the successive cycles. This plot can also be made versus potential (see fig. 1).

Fig. 1: overlaid frequency and current vs. E<sub>we</sub> of the polymer film growth. Scanning at 100 mV/s between 0 and 1.018 V.




The interesting point is to follow the frequency evolution versus potential on one cycle. The figure below shows the time evolution on one cycle.

Fig. 2: graphic zoom on one cycle showing the resonant frequency and the current density versus elapsed time<sup>[3]</sup>.



The mass calculation is done automatically in the process data tool of EC-Lab®. More details can be found in the application note<sup>[3]</sup>.



[1]: G. Sauerbrey, Phys. Verh., 1957, 8, 113-114.

[2]: G. Sauerbrey, Z. Phys., 1959, 155, 206-222.

[3]: Application note #13. Section "Apps & literature of EC-Lab division".

## Test cells



PAT-CORE  
and PAT-TRAY

ECC-PAT-CORE


These cells allow users to characterize the electrode (positive and negative electrode). With the ECC-Ref test cell (E-ECC-REF) a reference electrode (Lithium wire) can be used to follow the behavior of the negative and positive electrode simultaneously. This can be done with a VSP, VMP3 or MPG2 that offer the unique 5-lead type connection. Other cells for specific utilization are also available.

The advantages of the PAT-CORE series over the standard test cells are:

- more stable over time than the previous version
- less cross-contamination.
- higher throughput (Tray and single use consumable)/
- better EIS measurement (no artifact due to the cell)
- better reproducibility
- easy electrolyte filling
- no cleaning effort
- save time and space



3-electrode  
test cell



ECC-AIR



| PAT-CORE Series  |  | Catalog n°        |
|--|--|-------------------|
| <b>Leadless series</b>   |  |                   |
| Leadless 2/3 electrode test cell for battery materials (with PAT-core) - PAT-Cell. It has to be used with the PAT-Tray |  | <b>E-ECC-PATC</b> |
| Docking station for 16 electrochemical test cells (PAT-Cell) incl. data logger - PAT-Tray                              |  | <b>E-ECC-PATT</b> |
| <b>ECC series</b>  |  |                   |
| 2/3 electrode test cell for battery materials (ECC series with PAT-core) - ECC-PAT-Core                                |  | <b>E-ECC-PAT*</b> |
| ECC to PAT Upgrade kit - ECC1-00-0310-A  |  | <b>E-ECC-PATU</b> |

| Standard test cells  |  | Catalog n°        |
|--|--|-------------------|
| For lithium ion battery and other aprotic electrochemical system             | 2-electrode test cell                                | <b>E-ECC-STD</b>  |
|  | 3-electrode test cell                                | <b>E-ECC-REF</b>  |
|  | 2/3-electrode test cell                              | <b>E-ECCCOM</b>   |
|  | 2/3-electrode test cell for gas analysis             | <b>E-ECCDEM</b>   |
| For both aqueous and aprotic electrochemical systems with ref electrode      | 2/3-electrode test cell                              | <b>E-ECCAQU</b>   |
|  |  |                   |
| For measuring gas evolution and for drawing gas samples (pressure test cell) | ECC-press-DL   | <b>E-ECCP-DL</b>  |
|  | ECC-press-Aqu  | <b>E-ECCP-AQ</b>  |
|  | ECC-press-Air  | <b>E-ECCP-A</b>   |
|  | ECC-press  | <b>E-ECCP</b>     |
| For gas diffusion electrodes   | ECC-AIR electrochemical test cell (in aprotic media) | <b>E-ECCAIR</b>   |
|  | ECC-AIR electrochemical test cell                    | <b>E-ECCAIR-N</b> |
| For optical characterization   | ECC-Opto Std face to face                            | <b>E-ECCOPT-S</b> |
|  | ECC SBS side by side, complete                       | <b>E-ECC-SBS</b>  |

## Dilatometers

The dilatometer is a research grade instrument dedicated to the measurement of charge-induced strain (expansion and shrinkage) of electrodes down to the sub-micrometer range.

The dilatometers have been developed for the investigation of Li-ion battery and other insertion-type electrodes. It may be used in organic as well as aqueous electrolyte solutions.



| Dilatometer   | Catalog n°       |
|---|------------------|
| Electro chemical dilatometer ECD-nano-DL (complete)                   | <b>E-ECD-3N</b>  |
| Electro chemical dilatometer ECD-3-DL (complete)                      | <b>E-ECD-3DL</b> |
| <b>Options</b>  |                  |
| Cable to connect Bio-Logic instrument (DB9 connector) to ECD-nano-3DL | <b>E-C-3N</b>    |
| Cable to connect Bio-Logic instrument (DB9 connector) to ECD-3-DL     | <b>E-C-3DL</b>   |
| Upgrade kit from ECD-2N to ECD-3N or ECD-2DL to ECD-3DL               | <b>E-C-UP</b>    |

## Specifications

|                          | ECD-nano-DL   | ECD-3-DL                  |
|--------------------------|---|---------------------------|
| Displacement sensor      | capacitive  | LVDT                      |
| Dilatation range         | 250 μm  | 500 μm                    |
| linearity                | < 0.1% of full range  | < 0.15% of full range     |
| Resolution               | ≤5 nm   | ≤50 nm                    |
| Bandwidth                | ≥10 Hz  | ≥10 Hz                    |
| Drift (sample-free)      | ≤20 nm/h  | ≤100 nm/h                 |
| Load on test specimen    | 1.0 N (fixed)   | 0.3 N or 1.0 N (variable) |
| Chemical compatibility   | aprotic organic electrolytes, optional aqueous electrolytes |                           |
| Cell electrolyte volume  | ≤3 mL (ca. 2 mL electrolyte required)                       |                           |
| Temperature range        | -20 to 70 °C  |                           |
| Electrode configurations | 2-electrode, 3-electrode (reference), auxiliary             |                           |
| DC output voltage        | -10 to 10 V   |                           |
| Dimensions (H x W x D)   | 230 x 100 x 110 mm  |                           |
| Weight                   | 2.5 kg  |                           |

## Tools

### Cutting pliers

Catalog n°

18 mm **E-ECS18**

### Option

Cross hair to reduce cutting scrap, stop rail, assy **E-ECS-R**

### Punching tool for lithium foil

Catalog n°

18 mm **E-ECCL18**



### Electrode alignment and assembly tool

ECC CellLoad

**E-ECCLOAD**

### Options

REF - cleaning tool set **E-ECC-104**

REF - load **E-ECC-103**

Removing tool **E-ECC-102**

### Load the metal lithium

REF-load for test cells (allow us to load lithium metal inot REF bore of the E-ECC-REF)

**E-ECCLOADR**



## Rack

To get a compact setup, it is possible to put five units in a rack with a lateral tablet for a computer.

It is possible to use this rack for VMP3, booster chassis, HCP and CLB only if the sliding option is not purchased.

| Rack for MPG-2xx series up to 5 units                                       | Catalog n°      |
|---|-----------------|
| Including ethernet switch, power plugs, computer tablet                     | <b>092-85/1</b> |
| With sliding tablet including ethernet switch, power plugs, computer tablet | <b>092-85/2</b> |

## Specifications

|   |                     |
|---|---------------------|
| Maximum units                             | 5 units             |
| Dimensions rack                           | 600 x 710 x 1850 mm |
| (W x D x H) shelf with the sliding tablet | 495 x 450 x 257 mm  |
| shelf without the sliding tablet          | 495 x 450 x 295 mm  |
| Temperature range                         | 10 - 40° C          |

## Battery Holder

- 4-point measurement,
- 15 A maximum current,
- adaptable to a wide variety of battery sizes (coin cells, 18 650 or 26 650 cells),
- ability to link together many battery holders.

| Battery Holder                                  | Catalog n°       |
|---|------------------|
| BH-1i, dimensions (W x D x H) 134 x 79 x 185 mm | <b>092-22/15</b> |

## Coin Cell Holder

The coin cell holders are directly plugged to the potentiostat/cycler front panel connectors. No cable is required. The CCH is compatible with VMP3 and MPG2 multichannel potentiostat/galvanostat whereas the CCH-120/CCH-124 are compatible with the BCS-805 and BCS-810 cyclers.

| Coin Cell Holder   | Catalog n°       |
|--|------------------|
| CCH, dimensions (W x D x H) 82 x 83 x 65 mm for VMP3/MPG-2 | <b>092-22/14</b> |
| CCH-120 for BCS-805/810                                    | <b>096-120</b>   |
| CCH-124 for BCS-805/810                                    | <b>096-124</b>   |



## Current Collector

VSP-300/VMP-300 and BCS-815 offer the possibility to connect in parallel several channels and then increase the maximum current that can be passed through the cell. In order to simplify the connections, two current collectors have been designed, CC5 and CC8.


CC5 allows the user to connect up to five current boosters of VSP-300 or VMP-300 and with the CC-8 up to 8 channels of VMP-300 or BCS-815 can be connected in parallel.

| Current Collector for 5 boosters                       | Catalog n°       |
|--|------------------|
| CC5 Current Collector (5 channels) for VSP-300/VMP300  | <b>094-201/1</b> |
| CC8 Current Collector (8 channels) for VMP-300/BCS-815 | <b>096-015</b>   |

## Sense Adapter Module (SAM-50)

To be added to a multichannel system to perform stack measurements up to 50 V for 5 channel boards and a 10-element measurement.

| Sense Adapter Module | Catalog n°    |
|----------------------|---------------|
| SAM-50               | <b>092-26</b> |





## Membrane Electrode Assemblies (MEA)



These MEAs are based on a 50  $\mu\text{m}$ -thick perfluorinated membrane, and a platinum loading of 0.5  $\text{mg}/\text{cm}^2$  ( $\pm 0.05$ ), obtained with the standard catalyst concentration, 70% Pt/C.

- 3 layers MEA corresponds to a Catalyst Covered Membrane (CCM): the proton exchanged membrane is covered with an ink containing the catalyst on the electroactive area. Platinum and platinum-Ruthenium alloys are available for hydrogen and reformat operations respectively.
- 5 layers MEA corresponds to a CCM with Gas Diffusion Layers (GDL) completing the electrodes on the electroactive area on both sides of the membrane.
- 7 layers MEA corresponds to a 5-layer MEA with a Flat Gasket (FG) (fiberglass reinforced silicon) around both electrodes.


A subgasket, protecting the membrane and the edges of the active layer can be added to 5 and 7-layer MEAs.

### MEA selection

Only the 7-layer MEAs is fully active. The other 3-layer and 5-layer MEAs need to be modified.



| Membrane Electrode Assemblies (MEA) |                   |           | Catalog n° |
|-------------------------------------|-------------------|-----------|------------|
| Single cell fixture (without MEA's) | 5 $\text{cm}^2$   |           | 095-08/1   |
|                                     | 16 $\text{cm}^2$  |           | 095-08/2   |
|                                     | 25 $\text{cm}^2$  |           | 095-07     |
|                                     | 50 $\text{cm}^2$  |           | 095-08     |
|                                     | 100 $\text{cm}^2$ |           | 095-08/3   |
| Membrane Electrode Assemblies       | 5 $\text{cm}^2$   | 3 layers  | 095-09/13  |
|                                     |                   | 5 layers  | 095-09/14  |
|                                     |                   | 7 layers  | 095-09/15  |
|                                     | 16 $\text{cm}^2$  | 3 layers  | 095-09/16  |
|                                     |                   | 5 layers  | 095-09/17  |
|                                     |                   | 7 layers  | 095-09/18  |
|                                     | 25 $\text{cm}^2$  | 3 layers  | 095-09/1   |
|                                     |                   | 5 layers  | 095-09/2   |
|                                     |                   | 7 layers  | 095-09/3   |
|                                     | 50 $\text{cm}^2$  | 3 layers  | 095-09/4   |
|                                     |                   | 5 layers  | 095-09/5   |
|                                     |                   | 7 layers  | 095-09/6   |
| 100 $\text{cm}^2$                   | 3 layers          | 095-09/7  |            |
|                                     | 5 layers          | 095-09/8  |            |
|                                     | 7 layers          | 095-09/9  |            |
| 225 $\text{cm}^2$                   | 3 layers          | 095-09/10 |            |
|                                     | 5 layers          | 095-09/11 |            |
|                                     | 7 layers          | 095-09/12 |            |



## Educational kit

The educational fuel cell kit is very simple and easy to use: just plug the hydrogen tank into the cell and let the fuel cell use the oxygen in the surrounding air to produce electricity.

- 1 plastic electrolyser cell (25  $\text{cm}^2$ ) containing a 7-layer membrane electrode assembly. At the anode side of the electrolyser cell there is a water reservoir, the hydrogen gas produced comes out of the cell via one of the gas exits on the top and bottom side of the cell. A plastic screw is included for closing the unused gas exit.
- 1 plastic air-breathing fuel cell (25  $\text{cm}^2$ ).

The cell has one hydrogen gas input and one exit. A porous material is placed inside the kit for the homogeneous diffusion of hydrogen over the whole surface area of the membrane electrode assembly.

| Educational kit                           |   | Catalog n° |
|---|---|------------|
| Educational kit (fuel cell /electrolyzer) | 110 V electrical grid   | 095-20/A   |
|   | 220 V electrical grid   | 095-20/B   |
| <b>Content</b>                            |   |            |
| PEM kit                                   | PEM#1: 25 $\text{cm}^2$ air-breathing cells (MEA S25-7L N50-Pt/C 70%), no possibility to dismantle the cell<br>PEM#2: 25 $\text{cm}^2$ air-breathing cells (MEA S25-7L N125-Pt/C 70%), no possibility to dismantle the cell<br>PEM#3: 25 $\text{cm}^2$ air-breathing cells (MEA S25-7L N212-Pt/C 40%), no possibility to dismantle the cell | 095-20/1   |
|   |   | 095-20/2   |
|   |   | 095-20/3   |
|   | Electrolyzer kit (including gas line)   | 095-20/4   |
|   | Variable resistance   | 095-20/5   |
|   | Data acquisition (including connection, data logger, software)  | 095-20/7   |
|   | Power supply 220 V  | 095-20/9   |
|   | 110 V/220 V adapter for 110 V configuration   | 095-20/10  |
| <b>Options</b>                            |   |            |
|   | 25 $\text{cm}^2$ air-breathing cell (no MEA)  | 095-20/6   |
|   | Power supply by two photovoltaic cells  | 095-20/8   |

## Specific cables

By default, the potentiostat and the booster are provided with a cell cable of 1.5 m long.

The cable connected from the booster to the potentiostat is 0.8 m long for VMP3 based instruments.

For some applications, the user may need cables with different length. That's why, longer cable are offered.

Moreover for applications performed in glove box, hermetic cell cables are also offered.

### Set-up connection

Bad connection could induce some effects on the measurement (effect on the stability of the potentiostat, artifact...).

We recommend the user to optimize the set-up by using the accessories offered in this section.

## Connectors

|   | Catalog n°         | Black colour | Red colour         | Blue colour        | White colour       |
|---|--------------------|--------------|--------------------|--------------------|--------------------|
| Alligator clip with 2 mm connector (pack of 10)                     | <b>092-1001/1</b>  |              | <b>092-1001/2</b>  | <b>092-1001/26</b> | <b>092-1001/22</b> |
| Alligator clip 4 mm (pack of 5)                                     | <b>092-1001/13</b> |              | <b>092-1001/14</b> | <b>092-1001/23</b> |                    |
| Socket 2 mm female (pack of 10)                                     | <b>092-1001/5</b>  |              | <b>092-1001/6</b>  | <b>092-1001/7</b>  | <b>092-1001/8</b>  |
| Socket 4 mm female (pack of 10)                                     | <b>092-1001/25</b> |              | <b>092-1001/4</b>  | <b>092-1001/3</b>  |                    |
| Banana plug 2 mm (pack of 10)                                       | <b>092-1001/9</b>  |              | <b>092-1001/10</b> | <b>092-1001/11</b> | <b>092-1001/12</b> |
| Banana plug 4 mm (pack of 5)  |                    |              | <b>092-1001/16</b> | <b>092-1001/15</b> | <b>092-1001/24</b> |
| Banana adapter 2 mm female socket to 4 mm male plug (pack of 5)     | <b>092-1001/18</b> |              | <b>092-1001/19</b> | <b>092-1001/20</b> | <b>092-1001/21</b> |
| Adapter for banana 2 mm female socket to 4 mm male plug (pack of 5) | <b>092-1001/27</b> |              |                    |                    |                    |
| Banana adapter 4 mm female socket to 2 mm male plug (pack of 5)     | <b>092-1001/28</b> |              |                    |                    |                    |

| Connection kits    |  |  |                    |
|--------------------|--|--|--------------------|
| For standard board | - 4 alligator clips of 2 mm:<br>- 3 sockets of 2 mm:<br>- 3 alligator clips of 2 mm:<br>- 2 alligator clips of 4 mm: | blue, white, red, black,<br>blue, white, red<br>blue, white, red<br>red, black | <b>092-1001/30</b> |
| For booster board  | - 3 sockets of 2 mm:<br>- 2 sockets of 4 mm:   | blue, white, red<br>blue, white  | <b>092-1001/31</b> |



Hermetic cell cable for glove box

| Hermetic cell cable for glove box   | Available length | Catalog n°       |
|---|------------------|------------------|
| SP-50, SP-150, BiStat, VSP, VMP3, MPG2<br>A Jaeger connector 12 pins ensures the hermetic junction (hole diameter to make in the glove box wall: 27 mm):<br>- cable with 2 mm connectors in one side and 12 pins Jaeger connector in the other side (length 1 m),<br>- cable with DB-25 connector in one side and 12 pins Jaeger connector in the other side (length 1.5 m),<br>- one feedthrough seal  |                  | <b>092-23/5</b>  |
| SP-200, SP-240, SP-300, VSP-300, VMP-300<br>A Jaeger connector 25 pins ensures the hermetic junction (hole diameter to make in the glove box wall: 45 mm):<br>- cable with electrometer in one side and 25 pins connector in the other side (length: 1 m)<br>- cable with DB-25 connector in one side and 25 pins connector in the other side (length: 1.5 m)<br>- two feedthrough seals one installed in the glove box wall the other dedicated to the channel board calibration outside the box |                  | <b>094-101/6</b> |

| Connection cable from booster to potentiostat | Available length | Catalog n°      |
|---|------------------|-----------------|
| All boosters                                  | 3 m              | <b>092-33/5</b> |
|   | 5 m              | <b>092-33/6</b> |

| Longer cable                   | Available length | Catalog n°       |
|--------------------------------|------------------|------------------|
| SP-50, SP-150, VSP, VMP3, MPG2 | 2.5 m            | <b>092-23/2</b>  |
|                                | 3 m              | <b>092-23/7</b>  |
|                                | 5 m              | <b>092-23/3</b>  |
|                                | 10 m             | <b>092-23/4</b>  |
| Booster 2 A, 4 A, 5 A          | 2.5 m            | <b>092-33/11</b> |
|                                | 3 m              | <b>092-33/12</b> |
|                                | 5 m              | <b>092-33/13</b> |
| Booster 8 A, 10 A, 20 A        | 2.5 m            | <b>092-33/21</b> |
|                                | 3 m              | <b>092-33/22</b> |
|                                | 5 m              | <b>092-33/23</b> |



Nstat box (8 channels)

Bipot cable: dedicated to RRDE applications

## Multi-electrode investigation cables

| Multi-electrode investigation cables | Channel     | Nb of channels | Length | Catalog n°        |
|--------------------------------------|-------------|----------------|--------|-------------------|
| Nstat box (for VSP, VMP3)            | standard    | 4              | 1.5 m  | <b>092-16</b>     |
|                                      |             | 8              | 1.5 m  | <b>092-22/3</b>   |
|                                      | low current | 8              | 1.5 m  | <b>092-22/4</b>   |
| Bipot cable (for VSP, BiStat, VMP3)  | standard    | 2              | 1.5 m  | <b>092-22/12</b>  |
|                                      |             | 2              | 3 m    | <b>092-22/12A</b> |
|                                      |             | 2              | 4 m    | <b>092-22/12B</b> |

| Option  |  | Catalog n°      |
|---|--|-----------------|
| External power supply for the Nstat box (this option is needed if more than one Nstat box is connected to VMP3 or if the user uses VSP) |  | <b>092-16/1</b> |



Advanced Faraday cage

## Faraday cages

To avoid any external perturbations, especially for low current application, we recommend using the Faraday cage.

Please note that to activate this protection, Faraday cage has to be connected to the ground of the instrument (green plug on the rear panel of the instrument).



FC-45 Faraday cage

### Advanced Faraday cage specifications

|              |                           |
|--------------|---------------------------|
| Temperature  | 0 - 50 °C                 |
| Power supply | 100 VAC-240 VAC, 50/60 Hz |
| Fuse         | 1 A                       |
| Gas pressure | < 34 kPa max.             |
| Interface    | IEEE-488, RS-232C         |
| Weight       | 3.8 kg                    |

| Faraday cages                             | Catalog n°       |
|---|------------------|
| FC-45 Faraday cage, 450 x 450 x 450 mm    | <b>094-084/1</b> |
| Stand for FC-45                           | <b>094-084/2</b> |
| Standard Faraday cage, 400 x 200 x 600 mm | <b>NS-FAR600</b> |
| Advanced Faraday cage, 286 x 230 x 320 mm | <b>A-012033</b>  |



DC2 Dummy cell for booster

## Dummy cells

The dummy cell for booster is especially dedicated to check periodically a booster. It is provided with each booster chassis. The dummy cell for booster and DC2 can be bound together as well.

| Dummy cell             |   | Catalog n°       |
|------------------------|---|------------------|
| DC2                    | 1 R/C circuit.<br>For standard check            | <b>094-111/3</b> |
| Dummy cell for booster | 1 power resistor, 5 mOhm.<br>For standard check | <b>092-32/1</b>  |

## Dummy cell for booster specifications

|                         |            |
|-------------------------|------------|
| Resistance/mΩ           | 5          |
| Standard tolerance      | 1%         |
| Temperature coefficient | ±50 ppm/°C |



Test Box 2 Test Box 3

## Test Boxes

| Test Boxes |  | Catalog n°      |
|------------|--|-----------------|
| Test Box 2 | several circuits with high precision resistors. For calibration and validation | <b>092-22/6</b> |
| Test Box 3 | 3 circuits: linear, exponential, non-linear. For teaching and demonstration    | <b>092-22/7</b> |

## External device connection


| External device connection  | Catalog n°       |
|---|------------------|
| DB9-8 BNC connector for auxiliary I/O   | <b>092-22/1</b>  |
| PT100 temperature probe, to be connected to the auxiliary I/O, temperature range: -50 °C to 250 °C, dimensions: 3 x 20 mm, length of cable: 2.5 m | <b>092-22/13</b> |
| IS1 isolation module for auxiliary I/O for VSP300 based instruments   | <b>094-081/5</b> |



IS1

## Transport cases

| Transport cases | Catalog n°     |
|-----------------|----------------|
| SP-200          | <b>094-082</b> |
| SP-240/SP-300   | <b>094-091</b> |



SECM HR SKP SKP LEIS SVP

## Probes

A range of probes dedicated for use with our SECM, SVP, SKP, HR SKP and LEIS scanning probe applications are available for the M370 and M470 systems.

| Probes                            | Catalog n°         |
|-----------------------------------|--------------------|
| SECM 10 μm diameter Platinum disk | <b>U-23/10</b>     |
| SECM 15 μm diameter Platinum disk | <b>U-23/15</b>     |
| SECM 25 μm diameter Platinum disk | <b>U-23/25</b>     |
| SKP 500 μm diameter               | <b>U-SKP370/1</b>  |
| HR SKP 150 μm diameter            | <b>U-SKP-150</b>   |
| LEIS                              | <b>U-LEIS370/1</b> |
| SVP                               | <b>U-SVP370/1</b>  |



TriCell™



Shallow μTriCell™

## Cells

Three cells are available : The TriCell™ is a large volume, wide scan range cell, dedicated to LEIS, SVP, SKP, SDS techniques. The μTriCell™ and its Shallow version are dedicated to SECM techniques (dc, ac and ic mode). The Shallow μTriCell™ contains a slightly smaller volume of electrolyte than the μTriCell™ is more accessible.

All cells accommodate samples mounted in a 32 mm diameter resin cylinder.

| Cells             | Volume | Catalog n°         |
|-------------------|--------|--------------------|
| TriCell™          | ~1 L   | <b>U-TRICELL</b>   |
| μTriCell™         | 7 mL   | <b>U-μTRICELL</b>  |
| Shallow μTriCell™ | 6 mL   | <b>U-SμTRICELL</b> |

## VCAM3 Video Microscope System

The VCAM3 is a long working distance video microscope which allows users to view the distance between probe tip and sample surface in many scanning probe electrochemistry techniques.


### VCAM3 specifications

|                       |                               |
|-----------------------|-------------------------------|
| Min illumination      | 0.0003 lux                    |
| Field of view         | 1.4 mm (x4.5) to 8.6 mm (0.7) |
| Operation temperature | -30 to +70 °C                 |
| Catalog n°            | <b>U-VCAM3</b>                |





## Sales & technical contacts



Bio-Logic Science Instruments  
is present in most countries through an extensive distribution network

### Headquarters

#### **Bio-Logic SAS**

1, rue de l'Europe - 38 640 Claix

France

Phone: +33 476 98 68 31

Fax: +33 476 98 69 09

### Affiliate offices

#### **Bio-Logic USA, LLC**

P.O.Box 30009 - 37 930 Knoxville

USA

Phone: +1 865 769 3800

Fax: +1 865 769 3801

#### **Bio-Logic Science Instruments Pvt Ltd**

304, Orion Business Park, Next to Cine Wonder,

G. B. Road, Thane(W), 400 607 Mumbai

India

Phone: +91 2225842128



[www.bio-logic.info](http://www.bio-logic.info)