

Characterization of Food and Animal Feed Related Products by the Thermo Scientific FLASH 2000 Elemental Analyzer

Dr. Liliana Krotz, OEA Product Specialist, and Dr. Guido Giuzzi, OEA Product Manager, Thermo Fisher Scientific, Milan, Italy

Key Words

Food supplements, animal feed, meat, combustion, CHNS

Introduction

Food and animal feed is made up of chemical compounds that determine flavor, color, texture and nutritional value, and are carefully regulated by federal authorities and various international organizations to ensure that they are safe to eat and are accurately labelled.

One of the main analyses for quality control and R&D purposes is elemental characterization. The determination of nitrogen, carbon, hydrogen and sulfur, provides useful information on the characterization of these materials. It is therefore very important to have an accurate and precise technique, preferably automatic, that allows fast analysis with excellent reproducibility. The Thermo Scientific™ FLASH 2000 analyzer (Figure 1) copes effortlessly with the wide array of laboratory requirements such as accuracy, precision and day to day reproducibility.



Figure 1. FLASH 2000 Elemental Analyzer



Methods

For CHNS determinations, the elemental analyzer operates according to the dynamic flash combustion of the sample. The sample is weighed in a tin capsule and introduced into the combustion reactor via the Thermo Scientific™ MAS™ 200R autosampler together with a proper amount of oxygen. After combustion, the resultant gases are carried by a helium flow to a layer filled with copper, then swept through a GC column that separates the combustion gases, finally being detected by a thermal conductivity detector (TCD) (Figure 2). Total run time is less than 10 minutes. For NCS or for sulfur only determination, the water produced during combustion is adsorbed through a H₂O trap before entering the GC column.

For NC determination, after combustion, the produced gases are carried by a helium flow to a second reactor filled with copper, then swept through a H₂O trap, a GC column, before finally being detected by a thermal conductivity detector. Total run time is less than five minutes (see Figure 3).

A complete report is automatically generated by the Thermo Scientific™ Eager Xperience dedicated data handling software and displayed at the end of the analysis.

Results

Different food and animal feed related products with a large range of elemental concentrations were analyzed in various configurations to show the performance of the instrument in terms of repeatability.

For the CHNS, NCS and sulfur only determinations, the addition of Vanadium Pentoxide (oxygen donor), was used for a complete conversion of the sulfur, and the instrument was calibrated with the standards BBOT* and Nicotinamide. For NC determination, Acetanilide and Aspartic acid were used as standards to calibrate the instrument. In all cases, K factor was used as the calibration method.

* BBOT: 2,5-Bis (5-*tert*-butyl-benzoxazol-2-yl) thiophene




Figure 2. FLASH 2000 CHNS configuration




Figure 3. FLASH 2000 NC configuration

Table 1 shows the CHNS data obtained from different sample matrices. The weight of sample was 2 - 3 mg for animal gelatines and food supplement A, 3 - 4 mg for starch and food supplement C, 8 - 10 mg for food supplement B. No memory effect was observed when changing the sample nature, indicating complete combustion of all samples with quantitative determination of the elements.

Table 1. CHNS data

Sample	N %	RSD %	C %	RSD %	H %	RSD %	S %	RSD %
Fish gelatine	16.249	0.185	43.023	0.089	6.902	2.632	0.394	2.004
	16.212		43.099		6.608		0.408	
	16.189		43.051		6.586		0.408	
Bovine gelatine	15.796	0.148	44.615	0.037	6.623	0.309	0.531	0.601
	15.835		44.647		6.658		0.536	
	15.838		44.624		6.622		0.537	
Porcine gelatine	16.088	0.226	44.460	0.096	6.631	0.585	0.531	0.970
	16.016		44.397		6.659		0.536	
	16.043		44.379		6.582		0.537	
Starch	2.530	0.329	31.008	0.204	5.396	1.004	0.399	0.902
	2.516		30.850		5.456		0.396	
	2.537		31.000		5.415		0.391	
	2.520		30.956		5.373		0.392	
	2.528		30.967		5.310		0.398	
Food supplement A	13.168	0.137	52.179	0.104	6.665	0.311		
	13.160		52.084		6.626			
	13.194		52.178		6.626			
Food supplement B	0.071	1.109	8.0197	0.036	2.207	0.254		
	0.071		8.0251		2.199			
	0.073		8.0241		2.209			
Food supplement C	0.330	0.533	40.615	0.323	6.274	0.488	0.366	1.581
	0.332		40.396		6.331		0.368	
	0.329		40.630		6.323		0.357	

Table 2 shows the NCS data obtained from different sample matrices. The weight of the sample was 3 - 4 mg. No memory effect was observed when changing the sample nature, indicating complete combustion of all samples with quantitative determination of the elements.

Table 2. NCS data

	N %	RSD %	C %	RSD %	S %	RSD %
Animal feed 1	2.525	1.850	40.432	0.270	0.160	5.723
	2.523		40.342		0.143	
	2.444		40.215		0.155	
Animal feed 2	3.892	1.069	44.752	0.776	0.287	1.773
	3.845		44.893		0.282	
	3.810		44.235		0.277	
Meat 1	13.726	1.199	50.315	0.298	0.776	0.901
	14.057		50.107		0.790	
	13.931		50.026		0.784	
Meat 2	12.939	0.286	51.064	0.206	0.724	0.742
	12.979		50.886		0.732	
	13.013		50.879		0.721	
Meat 3	12.255	0.094	53.381	0.118	0.764	0.798
	12.243		53.486		0.774	
	12.232		53.495		0.764	
Meat 4	12.451	0.399	52.603	0.111	0.785	0.199
	12.369		52.701		0.783	
	12.459		52.597		0.781	
Potato tuber	0.874	1.330	39.542	0.229	0.0709	1.221
	0.893		39.528		0.0698	
	0.896		39.382		0.0716	
	0.868		39.389		0.0721	
	0.883		39.348		0.0709	

Table 3 shows the sulfur data obtained from different soya and maize samples. The weight of the sample was 3 - 4 mg.

Table 3. Sulfur data

Sample	S %	RSD %
Soya 1	0.356 - 0.336 - 0.338 - 0.350 - 0.343 - 0.341	2.214
Soya 2	0.351 - 0.344 - 0.343	1.635
Soya 3	0.372 - 0.373 - 0.363	1.725
Soya 4	0.366 - 0.363 - 0.364	0.388
Maize 1	0.115 - 0.114 - 0.116 - 0.113 - 0.113 - 0.106	3.137
Maize 2	0.119 - 0.111	4.919
Maize 3	0.104 - 0.104	0
Maize 4	0.112 - 0.102	6.608
Maize 5	0.102 - 0.100	1.396

Table 4 shows the NC data of blood flour, meat flour and food supplements. The weight of sample was 9 - 10 mg for blood flour and meat flour while for food supplements the weight was 10 – 20 mg.

Table 4. NC data

Sample	N %	RSD %	C %	RSD %
Blood flour	14.536	0.797	48.791	0.458
	14.307		49.238	
	14.427		49.052	
Meat flour	7.213	1.804	32.008	1.963
	7.474		33.287	
	7.385		32.722	
Food supplement 1	5.739	0.071	14.598	0.084
	5.747		14.623	
	5.744		14.612	
Food supplement 2	11.365	0.123	45.625	0.385
	11.391		45.303	
	11.369		45.343	
Food supplement 3	5.845	0.269	17.127	0.138
	5.831		17.115	
	5.863		17.161	
Food supplement 4	3.630	0.310	37.093	0.201
	3.613		36.975	
	3.634		36.956	

Conclusions

The FLASH 2000 analyzer allows the quantitative recovery of the elements from any matrix with no memory effect observed when changing the sample. The advantage of the FLASH 2000 analyzer lies in its ability to perform NC determination or simultaneous CHNS determination in a single run. By a simple modification of the CHNS configuration, the analysis of NCS or sulfur only can be performed using the same analytical conditions.

www.thermoscientific.com

©2014 Thermo Fisher Scientific Inc. All rights reserved. ISO is a trademark of the International Standards Organization. All other trademarks are the property of Thermo Fisher Scientific and its subsidiaries. This information is presented as an example of the capabilities of Thermo Fisher Scientific products. It is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

Africa +43 1 333 50 34 0
Australia +61 3 9757 4300
Austria +43 810 282 206
Belgium +32 53 73 42 41
Canada +1 800 530 8447
China 800 810 5118 (free call domestic)
 400 650 5118

Denmark +45 70 23 62 60
Europe-Other +43 1 333 50 34 0
Finland +358 9 3291 0200
France +33 1 60 92 48 00
Germany +49 6103 408 1014
India +91 22 6742 9494
Italy +39 02 950 591

Japan +81 45 453 9100
Latin America +1 561 688 8700
Middle East +43 1 333 50 34 0
Netherlands +31 76 579 55 55
New Zealand +64 9 980 6700
Norway +46 8 556 468 00
Russia/CIS +43 1 333 50 34 0

COMPANY WITH QUALITY MANAGEMENT
 SYSTEM CERTIFIED BY DNV
 = ISO 9001:2008 =

Thermo Fisher Scientific S.p.A. Milan, Italy
 is ISO 9001:2008 Certified.

Singapore +65 6289 1190
Spain +34 914 845 965
Sweden +46 8 556 468 00
Switzerland +41 61 716 77 00
UK +44 1442 233555
USA +1 800 532 4752

Thermo
 SCIENTIFIC

Part of Thermo Fisher Scientific