

# Petrochemical Compounds Characterization with the Thermo Scientific FLASH 2000 CHNS/O Elemental Analyzer

Liliana Krotz and Guido Giuzzi  
Thermo Fisher Scientific, Rodano, Milan, Italy

## Overview

**Purpose:** Characterization of petrochemical compounds by Organic Elemental Analysis.

**Methods:** Samples were analyzed using an elemental analyzer with automated autosamplers.

**Results:** CHNS/O data are shown to demonstrate the performance of the FLASH 2000 Organic Elemental Analyzer.

## Introduction

The Thermo Scientific™ FLASH 2000 Organic Elemental Analyzer (Figure 1), which is based on the dynamic combustion of the samples, provides quantitative and automated simultaneous CHNS determination and the oxygen determination by pyrolysis. Solid, viscous, liquid and gas samples can be analyzed without any matrix effect. From the CHNS/O data, the dedicated Thermo Scientific™ Eager Xperience Data Handling Software calculates automatically the GHV and NHV values (Gross Heat and Net Heat, both expressed in kcal/kg) and the CO<sub>2</sub> Emission Trade data.

**FIGURE 1. Thermo Scientific FLASH 2000 with MAS 200R and AS 1310 Autosamplers.**



## Methods


For CHNS determination the Elemental Analyzer operates according to the dynamic flash combustion of the sample. Samples are weighed in tin containers and introduced into the combustion reactor via the Thermo Scientific™ MAS 200R Autosampler or injected directly by a syringe through the Thermo Scientific™ AS 1310 Liquid Autosampler with the proper amount of oxygen. After combustion the resulted gases are carried by a helium flow to a layer filled with copper, then swept through a GC column, providing the separation of the combustion gases. Finally, they are detected by a thermal conductivity detector (TCD). Total run time is less than 10 min. (Figure 2).

For trace sulfur determination, the gases produced by combustion are carried by a helium flow to a layer filled with copper, then swept through a water trap, a short GC column. Finally the sulfur is measured by the flame photometric detector (FPD). Total run time 5 min. (Figure 3).


For oxygen determination, the system operates in pyrolysis mode. Samples are weighed in silver containers and introduced into the pyrolysis chamber via the MAS 200R Autosampler or injected directly by a syringe through the AS 1310 Liquid Autosampler. The reactor contains nickel coated carbon at 1060°C. The oxygen present in the sample, combined with the carbon, forms carbon monoxide which is then chromatographically separated to other products and detected by the TCD detector (Figure 2).

A comprehensive report is automatically generated by the Eager Xperience Data Handling Software and displayed at the end of the analysis.

**FIGURE 2. CHNS/O configuration**



**FIGURE 3. Sulfur configuration by FPD detector**



## Results

Different solids, viscous and liquid petrochemical samples were chosen to show the reproducibility obtained with the system. Coal, coke, graphite, lignite and catalyst samples were homogenized by a ball mill while viscous and liquid samples were analyzed without pre-treatment.

Table 1 and 2 show CHNS/O and CHNS determination of different matrices. Instrument calibration was performed with 2-3 mg of BBOT (2, 5-Bis (5-ter-butyl-benzoxazol-2-yl) thiophene), while sample weight was 2 – 3 mg. No matrix effect was observed when changing the nature of sample. Table 1 indicates also the heat value GHV (Gross Heat Value in kcal/kg) and NHV (Net Heat Value in kcal/kg) calculated automatically by the Eager Xperience Data Handling Software.

**TABLE 1. CHNS/O determination and Heat Value calculation of coal and lignite samples**

Sample	N %	C %	H %	S %	O %	GHV	NHV
Coal	1.740	73.006	5.389	0.838	12.733	7316	7040
	1.732	72.950	5.358	0.857	12.731	7316	7040
	1.747	73.238	5.398	0.828	12.813	7313	7037
RSD%	0.431	0.209	0.390	1.752	0.367	0.024	0.025
Hard Coal	1.287	80.137	4.621	0.488	4.903	7957	7720
	1.288	80.123	4.513	0.486	5.085	7918	7687
	1.329	80.706	4.617	0.497	5.094	8003	7766
RSD%	1.842	0.414	1.336	1.195	2.144	0.535	0.514
Brown Coal	1.988	77.258	3.280	0.386	3.730	7313	7145
	2.006	77.196	3.266	0.386	3.723	7304	7136
	2.028	77.823	3.292	0.397	3.690	7364	7195
RSD%	0.998	0.446	0.397	1.630	0.575	0.442	0.444
Mineralized Lignite	0.275	16.512	2.487	0.222	13.036	1651	1524
	0.273	16.480	2.476	0.226	13.046	1651	1524
	0.274	16.553	2.469	0.223			
RSD %	0.290	0.219	0.366	0.962	0.056	0.019	0.020
Lignite 1	0.758	62.281	4.636	0.377	25.202	5589	5357
	0.750	62.357	4.544	0.373	25.540	5574	5343
	0.757	62.240	4.377	0.371			
RSD %	0.559	0.074	2.899	0.724	0.942	0.184	0.192

**TABLE 2. CHNS determination of petrochemical samples**

Sample	N%	RSD%	C%	RSD%	H%	RSD%	S%	RSD%
Pet coke	1.486	2.266	96.329	0.313	0.210	8.095	0.641	0.384
	1.534		96.756		0.235		0.643	
Coke 1	1.229	0.011	86.195	0.384	4.343	0.346	0.656	4.535
	1.228		86.665		4.322		0.616	
Coke 2	0.377	1.360	98.997	0.172	-	-	0.422	1.017
	0.370		98.917		0.431			
	0.374		98.576		0.429			
	0.363		98.973		0.434			
	0.369		98.836		0.429			
Catalyst	0.007	4.478	0.165	1.716	1.164	0.552	0.408	1.253
	0.007		0.166		1.154		0.409	
	0.006		0.171		1.153		0.417	
Crude oil	0.208	2.802	84.701	0.130	12.368	0.545	2.325	1.513
	0.197		84.563		12.432		2.381	
	0.206		84.781		12.297		2.393	

The performance of the FLASH 2000 OEA was evaluated by comparing the repeatability of the CHN data with the ASTM D5375 requirements showed in Table 3. The method covers the instrumental determination of nitrogen, carbon and hydrogen in coal and coke samples.

**TABLE 3. Concentration Range and Limit of Repeatability accepted by ASTM D5375-02**

Element	Concentration range (%)	Repeatability Limit (r)
Carbon	48.6 to 90.6	0.64
Hydrogen	0.14 to 5.16	0.16
Nitrogen	0.69 to 1.57	0.11

Repeatability Limit (r): the value below which the absolute difference between two test results calculated to a dry basis of separate and consecutive test determinations, carried out on the same sample, in the same laboratory, by the same operator, using the same apparatus.

The precision of the FLASH 2000 OEA for CHN determination was demonstrated performing two series of samples analyzed in duplicate. The calibration of the system was performed with 2-3 mg of acetanilide. Samples were weighted at 2-3 mg. Table 4 shows the CHN data obtained and the difference (Diff.) calculated between both data. All differences fall within or below the Repeatability Limit indicated in the official method, meaning the good homogeneity and the complete combustion of the samples.

**TABLE 4. CHN data of coal and coke samples.**

Serie	Elements	Coal 1		Coal 2		Coke	
		%	Diff.	%	%	%	Diff.
1	N	1.28	0	1.03	1.03	1.03	0
		1.28		1.03	1.03		
	C	86.12 86.59	0.47	87.47 87.61	87.47 87.61	87.47 87.61	0.09
2	H	4.43	0.09	0.30	0.30	0.30	0.01
		4.52		0.29	0.29		
	N	1.31 1.30	0.01	1.04 1.05	1.04 1.05	1.04 1.05	0.01
2	C	86.81	0.61	87.00	87.00	87.00	0.46
		86.20		87.46	87.46		
	H	4.47 4.43	0.04	0.30 0.29	0.30 0.29	0.30 0.29	0.01

Table 5 shows the repeatability of CHNS/O determination carbon black samples while Table 6 shows the oxygen data and the heat values GHV and NHV (Gross Heat Value and Net Heat Value, both expressed in kcal/kg), and the CO<sub>2</sub> Emission Trade data calculated by the Eager Xperience Data Handling Software. For CHNS, the analyzer was calibrated with 2-3 mg of BBOT\* standard using K factor as calibration method while the sample was weighted at 2-2.5 mg. For oxygen determination, 1-2 mg of BBOT was analyzed as standard using K factor while the sample was weighted at 1-2 mg.

**TABLE 5. CHNS data of carbon black samples.**

Sample	N %	RSD%	C %	RSD %	H %	RSD %	S %	RSD %
A	0.188	0.676	96.559	0.111	0.300	1.843	0.829	2.193
	0.185		96.734		0.295			
	0.186		96.817		0.304			
	0.186		96.706		0.308			
B	0.273	0.810	96.178	0.159	0.267	2.086	0.823	0.493
	0.277		95.899		0.280			
	0.272		95.351		0.272			
	0.273		96.079		0.277			
	0.273		96.079		0.277			

**TABLE 6. Oxygen data, Heat values and CO<sub>2</sub> value of carbon black samples**

Sample	O %	RSD %	GHV (kcal/lsg)	RSD %	NHV (kcal/kg)	RSD %	CO <sub>2</sub> E.T.	RSD %
A	1.149	0.397	8005.53	0.003	7990.03	0.003	105.96	0.000
	1.158		8005.13		7989.63		105.96	
	1.152		8005.42		7989.93		105.96	
B	1.661	0.623	7915.43	0.006	7901.37	0.006	106.37	0.005
	1.679		7915.64		7900.58		106.38	
	1.661		7915.44		7901.38		106.37	

Table 7 shows the reproducibility of CH determination in diesel samples while Table 8 shows the reproducibility of oxygen determination of diesel samples, analyzed by automated liquid injection. The volume injected was 2  $\mu$ l.

**TABLE 7. Reproducibility of CH determination of diesel samples.**

Sample	Diesel 1		Diesel 2		Diesel 3	
	C %	H %	C %	H %	C %	H %
%	85.34	14.32	84.67	15.37	86.45	13.87
	85.35	14.31	84.72	15.39	86.51	13.92
	85.64	14.39	84.49	15.32	86.26	13.82
	85.67	14.37	84.62	15.34	86.32	13.86
	85.72	14.38	84.53	15.36	86.41	13.90
	85.73	14.38	84.71	15.36	86.59	13.88
	85.61	14.37	84.60	15.36	86.57	13.90
	85.25	14.34	84.55	15.32	86.52	13.93
	85.26	14.37	84.37	15.35	86.67	13.85
	85.67	14.39	84.43	15.33	86.47	13.88
	RSD %	0.23	0.20	0.14	0.14	0.02

**TABLE 8. Reproducibility of Oxygen determination of diesel samples.**

Sample	O %	RSD %
A	0.0265	8.179
	0.0224	
	0.0214	
	0.0233	
	0.0233	
B	0.0191	3.848
	0.0192	
	0.0181	
	0.0196	
	0.0201	

High density diesel samples were weighed in tin capsules for CHNS analysis and in silver capsule for oxygen determination. Table 9 shows the CHNS/O data and the Heat Values of two diesel samples.

**TABLE 9. CHNS/O and Heat Value calculation of diesel samples.**

Sample	N %	C %	H %	S %	O %	GHV	NHV
1	0.0618	90.54	9.39	0.0587	0.3002	10211	9709
	0.0644	90.20	9.39	0.0585	0.2711	10212	9710
	0.0657	90.15	9.38	0.0596	0.2862	10212	9710
2	0.0579	90.06	9.70	0.0342	0.2259	10289	9770
	0.0681	90.06	9.72	0.0398	0.2560	10288	9769
	0.0614	90.02	9.71	0.0381	0.2431	10288	9769

Table 10 shows the sulfur data of petrochemical samples obtained with the FPD detector (sulfur content is at trace levels). Gasoline and diesel samples were weighed in hard tin capsules closed by a dedicated sealing device. The calibration was performed using reference materials at trace sulfur content.

**TABLE 10. Sulfur determination by FPD detector.**

Sample	S ppm	RSD %
Catalyst 2	13	9.897
	11	
	11	
Coke	398	1.507
	404	
	392	
Crude oil	376	3.619
	371	
	397	
Gasoline	66	2.257
	68	
	69	
Diesel	25	9.03
	22	
	21	
	22	

## Conclusion

The Thermo Scientific FLASH 2000 Elemental Analyzer, based on Dumas method, proved to be a valuable solution for the elemental analysis of petrochemical samples in terms of accuracy, repeatability and sensitivity of results. Its automation, high speed of analysis and the reduction of long sample preparation processes allow efficient analysis and help reduce overall operational costs.

All data were obtained with good reproducibility and no matrix effect was observed when changing the sample indicating complete combustion. CHN data of coal and coke were compliant to the ASTM D5375 method.

The advantage of the FLASH 2000 Analyzer lies in the possibility to perform CHNS determination in a single run, then changing the configuration oxygen determinations without any extra modules.

By introducing minor adaptations to the configuration, trace of sulfur can be analyzed through the FPD detector.

Thanks to the coupling of the liquid Autosampler to the elemental analyzer, liquid samples are analyzed with high precision and accuracy. The coupling proved to be easy to install, with automated syringe alignment, the correction of the volume in weight is automated and performed by the Eager Xperience Data Handling Software.

[www.thermoscientific.com](http://www.thermoscientific.com)

©2015 Thermo Fisher Scientific Inc. All rights reserved. ISO is a trademark of the International Standards Organization. ASTM is a trade mark of American Society for Testing and Materials. All other trademarks are the property of Thermo Fisher Scientific and its subsidiaries. This information is presented as an example of the capabilities of Thermo Fisher Scientific products. It is not intended to encourage use of these products in any manners that might infringe the intellectual property rights of others. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

COMPANY WITH QUALITY MANAGEMENT  
SYSTEM CERTIFIED BY DNV  
= ISO 9001:2008 =

Thermo Fisher Scientific S.p.A. Milan, Italy  
is ISO 9001:2008 Certified.

**Africa** +43 1 333 50 34 0  
**Australia** +61 3 9757 4300  
**Austria** +43 810 282 206  
**Belgium** +32 53 73 42 41  
**Canada** +1 800 530 8447  
**China** 800 810 5118 (free call domestic)  
400 650 5118

**Denmark** +45 70 23 62 60  
**Europe-Other** +43 1 333 50 34 0  
**Finland** +358 10 3292 200  
**France** +33 1 60 92 48 00  
**Germany** +49 6103 408 1014  
**India** +91 22 6742 9494  
**Italy** +39 02 950 591

**Japan** +81 45 453 9100  
**Korea** +82 2 3420 8600  
**Latin America** +1 561 688 8700  
**Middle East** +43 1 333 50 34 0  
**Netherlands** +31 76 579 55 55  
**New Zealand** +64 9 980 6700  
**Norway** +46 8 556 468 00

**Russia/CIS** +43 1 333 50 34 0  
**Singapore** +65 6289 1190  
**Spain** +34 914 845 965  
**Sweden** +46 8 556 468 00  
**Switzerland** +41 61 716 77 00  
**UK** +44 1442 233555  
**USA** +1 800 532 4752

PN42239-EN 1115S

**Thermo**  
SCIENTIFIC

A Thermo Fisher Scientific Brand