



Thermo Scientific
ITQ Series Quadrupole
Ion Trap GC/MS™

Optimized
for real-world sample analysis

Thermo
SCIENTIFIC

Thermo Scientific

ITQ series quadrupole ion traps

the most sensitive ion trap mass spectrometers available today

The Thermo Scientific ITQ Series of external ionization ion trap mass spectrometers is the latest generation in our industry-recognized legacy of ion trap instrumentation – GC/MS[®]. For unmatched results, day in and day out – whether it's in an environmental lab monitoring the water quality near industrial facilities, a state forensic laboratory tracking clues in a criminal investigation, or a food safety laboratory striving to protect food supplies in a global economy, the ITQ™ Series fits your laboratory's needs.

- Significantly less chemical background for enhanced detectivity with the proven high temperature ion source
- Precision quantitation day after day with electronic flow control of CI reagent gas
- Maximum up-time with simple ion volume change for rapid routine maintenance
- Superior manufacturing quality control
- Computer independence and networking with Ethernet communication
- Exceptional MS/MS performance with automatic waveform optimization for more efficient ion storage
- Tuning for USEPA-regulated methods and selection of standard libraries for searching and identification
- Versatile analysis modes including simultaneous full-scan MS and MS/MS, MSⁿ, PPINICI™*, and Data Dependent scanning

* Pulsed Positive Ion/Negative Ion Chemical Ionization

The ITQ Series of GC/MS[®] instruments offers your laboratory choices that allow you to match the GC/MS system you select to your lab's needs.

- **The Thermo Scientific ITQ 700:** The ideal selection for routine, full-scan GC/MS analyses for the cost-conscious laboratory.
- **The Thermo Scientific ITQ 900:** Tap into an expanded full-scan mass range coupled to a powerful GC with a full-featured touch screen interface.
- **The Thermo Scientific ITQ 1100:** The ultimate choice for applications from research to routine, with powerful new tools that expand your lab's capabilities.

The ITQ Series is designed with the ability to upgrade in mind, protecting your investment by adapting to your lab's changing workflows and needs over time. If your needs change, upgrade your instrument to gain access to new features, greater flexibility, more power. Better yet, regardless of your choice, you will have the most sensitive GC-ion trap mass spectrometer available, giving you lower detection limits, even in matrix.

The ITQ Series offers a range of operating modes, from full-scan MS and MS/MS (MSⁿ), to positive and negative chemical ionization. Dual modes for sequential full scan and MS/MS or positive ion/negative ion chemical ionization allow you to acquire both types of data in a single injection. Smart Data Dependent scanning allows you to quickly collect data, confirm the identity of compounds, and further reduce sample cleanup costs. Variable damping gas, an option available exclusively on the ITQ Series, further improves GC/MS sensitivity up to 5X or more across a broad range of real-world samples.



Unsurpassed External Ionization Ion Trap Technology

Every ITQ system features a unique, high temperature external ionization ion source to effectively handle dirty samples and significantly improve response stability. The enhanced source incorporates an external passive collector and utilizes construction materials that are able to withstand higher operating temperatures. With an upper temperature limit of 300 °C, the ITQ source stays cleaner longer, and hundreds of samples can be analyzed before routine source maintenance is required.


Unmatched Full-Scan Sensitivity for Routine GC/MS Applications

Each ITQ instrument offers the best full-scan electron ionization (EI) sensitivity available today, meaning your laboratory can easily achieve low detection limits, even in matrix. External ionization and advanced tuning algorithms ensure library-searchable mass spectral data that allow for confidence in identification and quantification. Whether it's used for routine environmental analysis, industrial quality control, spectral confirmation, or as a university training instrument, the ITQ Series offers the performance you need for routine methods.

Upgrade to MSⁿ for Selective Matrix Elimination

Adding the power of mass spectrometry/mass spectrometry (MSⁿ) to your ITQ system allows your lab to routinely test for compounds that cannot typically be detected with other GC/MS techniques. MSⁿ provides the ultimate in selectivity for target analytes, eliminating false positives and false negatives. This powerful tool comes standard on the ITQ 1100 system, and is an available option for the ITQ 700 and ITQ 900 instruments. MSⁿ is applicable to most compounds, and is as simple to use as selected ion monitoring (SIM). The superb selectivity allows detection and quantitation at sub-picogram levels for target analytes in matrices such as plant and animal tissues, soils and sludge, biological fluids, and milk.

- Unequaled performance in difficult matrices
- Setup as simple as selected ion monitoring
- Far higher degree of confidence – eliminates false positives and false negatives



Innovative technology powers the ITQ series, and your productivity

the right system, for right now and into the future

Wider Dynamic Range

Quantitative performance is a vitally important aspect of any GC/MS system. Our engineers have taken the strength of our system and made it even better. Improvements to the algorithms used to trap ions have enhanced linear dynamic range and spectral quality. The quantitation and spectral consistency range extends from femtograms (fg) up to nanograms (ng) and beyond.

Exchange Ion Volumes in Minutes

Removable ion volumes provide added versatility and convenience for all types of analysis. For example, you can easily switch from GC/MS to probe analysis in under three minutes. With the vacuum interlock option, you can exchange the ion volume without breaking vacuum. This allows you to quickly replace a dirty ion volume with a clean one or change to an ion volume optimized for electron ionization (EI), chemical ionization (CI), or the combined EI/CI ion volume, depending on your application. With the versatile EI/CI volume, use digital CI with Data Dependent acquisition to obtain the molecular ion, the MS/MS spectra, and the EI library-searchable spectra – in one injection. Now that's a way to increase your sample throughput – fast!



Removable ion volumes allow quick switching from one technique to another, as well as easy changeover to solids probe analysis.

Easy Positive and Negative Chemical Ionization

The reagent gas pressure inside the ion source determines the quantitative reproducibility. These variations are minimized in ITQ systems equipped with chemical ionization (CI) by incorporating electronic flow control of the reagent gas. The method software controls this flow so that any time the data system activates the analytical method, a precise flow of the CI reagent gas is delivered to the mass spectrometer. Calibration curves are stable from day to day and quantitative precision is improved. Positive and negative chemical ionization are easy to set up, and a wide variety of reagent gases are supported through the software. In addition, proprietary Pulsed Positive Ion Negative Ion/Chemical Ionization (PPINICI) allows for the acquisition of both positive and negative chemical ionization data in a single injection.

- Digital electronic flow control of reagent gas
- Improved quantitative accuracy and precision
- Reproducible ion ratios

ITQ Series Features


- External ionization ion source – standard across the ITQ Series
- Wide spectral/quantitative dynamic range
- Full-scan EI mode for identifying unknowns – library search information
- Removable ion volumes
- ± 10 kV conversion dynode
- Powerful Thermo Scientific Xcalibur data system

Analytical Versatility

- Three mass range options – match MS capabilities to your analytical needs
- MSⁿ mode ($n \leq 5$) for targeting difficult compounds
- PCI for molecular weight confirmation
- NCI (ECD-MS) for ultimate sensitivity and selectivity
- Simultaneous full-scan MS and MS/MS in a single analysis
- Variable damping gas option for enhanced sensitivity
- PPINICI – hardware/software option with alternating PCI/NCI scans
- Data Dependent scanning option for advanced spectral pattern recognition

Expanded Capabilities

- Vacuum interlock
- Liquid and headspace autosamplers
- DEP/DIP solids probes
- High-capacity, 250 L/s turbopump system
- Two positions available for traditional GC detectors
- Turn-key applications packages



New advances in ion trap performance and ease of use

easily optimize and improve your MSⁿ experiments

Continuing with our legacy of innovation, the ITQ Series comes with two new advances in technology, designed to improve data quality and system ease of use. Now, the advantages of MS/MS are more readily achieved.

Automated Collision Energy

For the ITQ, we have introduced new Automated Collision Energy (ACE), which greatly simplifies MSⁿ optimization, making this powerful tool easier and more routine. Determining the collision energy setting required to ensure maximal product ion generation typically means performing a number of separate experiments using different collision energy settings. The results of these experiments would then have to be manually evaluated to determine which setting resulted in the greatest product ion intensities. Any changes – such as to damping gas settings to get more sensitivity – requires a repeat of the process. Plus, the more compounds, the more complex and time-consuming the determination becomes.


ACE automatically calculates an estimated optimal collision energy based on the operating parameters of the scan. Because of differences in the ease of which ions are fragmented, this calculation is an estimated value. ACE runs three energies in a single scan to ensure optimal fragmentation and product ion production. The optimal energy no longer needs to be determined – ACE covers all of the bases every time you run MSⁿ. If you want to run at a fixed collision energy, the estimated optimal collision energy calculated by ACE gives you a good starting point.

- Three collision energies in each scan – the estimated optimal plus a lower and higher value
- Easily apply the power of MSⁿ for routine use
- Simplified starting point for manual methods

Advanced Pulsed Q Dissociation Mode – For Enhanced MS/MS Quality

Our engineers originally developed the PQD technique for use in our line of liquid chromatography-linear trap mass spectrometers. Now, for the first time, this advanced technology is available on the ITQ 1100 GC/MSⁿ system. PQD is used to generate spectral data that are qualitatively similar to data produced using standard collision induced dissociation (CID). The key difference is that PQD increases the overall spectral quality, particularly in low mass range, which lets you see low *m/z* fragments that are usually excluded from CID spectra. PQD also helps you access higher energy dissociation channels. PQD is a novel fragmentation mechanism that involves precursor ion activation at high *q*, a time delay to allow the precursor to fragment, then a rapid pulse to low *q* where all fragment ions are trapped. The product ions can then be scanned out of the ion trap and detected.

- Generate information-rich mass spectral data, even in the low mass range
- Eliminates low-mass cutoff
- Produces precise, reproducible fragmentation




Thermo Scientific ITQ 700 GC/MS

the ideal choice for small budgets, tight spaces

The Thermo Scientific ITQ 700 is fully automated and designed for the analyst who performs routine, high volume GC/MS applications. The optimized external ion source design meets all standard tuning requirements and allows quantitation of samples from the low picogram to the mid-nanogram range. Whether it's used for routine environmental analysis, industrial quality control, spectral confirmation, or as a university training instrument, the ITQ 700 fulfills your lab's needs at an attractive entry-level price.

- Value-packed system for any budget
- Superior software from the MS leader
- Enhanced low mass spectral quality

The ITQ 700 was designed for the laboratory seeking a basic GC/MS system for routine full scan quantitation or basic teaching applications. Combining the sensitivity of full scan ion trap MS with the user-exchangeable plug-in injector and detector technology of the Thermo Scientific TRACE 1300 GC, this package is perfect for labs with limited budgets. A mass range to 700 u covers most general GC/MS applications including environmental, QA/QC, forensic, and general teaching and instructional settings.



Easily meet BFB and DFTPP tune requirements for EPA methods

ITQ 700 Highlights

Standard Features


- External ionization ion source
- Scan range of 10–700 u
- Most sensitive EI full scan GC/ion trap MS
- Reliable 70 L/s turbomolecular pump
- Powerful Xcalibur™ data system – for instrument operation, data analysis and reporting
- TRACE™ 1300 GC with user-exchangeable plug-in injector and detector technology

Upgrade Options

- 250 L/s turbomolecular pump
- Vacuum Probe Interlock – for easy source maintenance without venting
- MSⁿ with ACE and PQD
- Chemical Ionization – expanded analytical capabilities
- Direct sample probe – analyze samples without chromatographic separation and without venting or moving the GC

Added Value Thermo Scientific Software Options

- Thermo Scientific TraceFinder software follows your lab's workflows



Thermo Scientific ITQ 900 GC/MS

enhanced analytical flexibility with
added standard features

The Thermo Scientific ITQ 900 combines the power of ion trap mass spectrometry with the increased productivity and ease of use of the technology-leading Thermo Scientific TRACE 1310 gas chromatograph.

Larger routine QA/QC laboratories will benefit from the TRACE™ 1310 GC, which features an icon-based touch screen control. While retaining the unique plug-in flexibility of the TRACE 1300 GC, the touch screen interface allows for temperature status updates, maintenance commands, run logs, multiple language capabilities and video tutorials that drive simple instrument interaction.

The ITQ 900 also features a broader mass range, up to 900 u, allowing it to detect a wider range of compounds. Tap into even more analytical capability with an upgrade to MSⁿ, chemical ionization and solids probes.

- Analyze a broader range of targets and unknowns with expanded detector and injector options
- Upgradeable MS options for increased versatility
- On-instrument GC control through easy-to-use touch screen interface

Larger routine QA/QC laboratories will benefit from the TRACE 1310 GC which features a complete touch screen and an icon-driven easy-to-use local interface which is ideal for direct instrument control when method development and storage are required.

While retaining all of the capabilities and performance of the TRACE 1300 GC model, local status update of the oven, injectors and detectors, maintenance commands, run log, multiple language capabilities and video tutorials are additional tools that drive simple instrument interaction.

ITQ 900 Highlights

Standard Features

- External ionization ion source
- Scan range of 10–900 u
- Most sensitive EI full scan GC/ion trap MS
- Reliable 70 L/s turbomolecular pump
- Xcalibur data system – for instrument operation, data analysis and reporting
- TRACE 1310 GC – high-performance GC with user-exchangeable plug-in injector and detector technology, and an easy-to-use touch screen interface

Upgrade Options

- Vacuum Probe Interlock – for easy source maintenance without venting
- Add GC detectors, plus select from a broad range of injection ports
- 250 L/s turbomolecular pump
- MSⁿ with ACE and PQD
- Chemical Ionization – expanded analytical capabilities
- Direct sample probe – analyze samples without chromatographic separation and without venting or moving the GC

Added Value Thermo Scientific Software Options

- TraceFinder software follows your lab's workflows



Thermo Scientific ITQ 1100 GC/MSⁿ

the ultimate performer, for applications
from research to routine

The Thermo Scientific ITQ 1100 is the ultimate ion trap-based GC/MS system, perfect for the laboratory seeking the most powerful and flexible GC/MS platform. This system provides an extended mass range – up to 1100 u – to increase the number of compounds which can be detected and identified. New advanced MSⁿ functions – Automated Collision Energy (ACE) and our proprietary Pulsed Q Dissociation Mode (PQD) are standard on the ITQ 1100. The ITQ 1100 comes standard with a full suite of popular options, including the 250 L/s turbomolecular pump, vacuum probe interlock, MSⁿ, sequential full scan/MS-MS and Data Dependent Scanning. Available options include the direct insertion probe and chemical ionization.

New Advances Offer Enhanced Performance and Ease of Use


New to the ITQ 1100 system are two advanced modes of operation, each designed to support your lab's needs for advanced benchtop GC/MS capabilities. Automated collision energy (ACE) simplifies the development of MSⁿ methodologies by automatically determining the proper energy set up for the collision-induced dissociation step. MS/MS experiments are easier to develop, which makes it easier for your lab to apply this powerful technique, even for routine applications. The Thermo Scientific Pulsed Q Dissociation Mode (PQD) increases the number of product ions formed during the CID process, yielding greater information for the qualitative MSⁿ experiment. Both of these advanced functions come with the ITQ 1100 to let your laboratory unlock the true potential of MS/MS, and expand into new markets and application areas.

Restore Performance or Change Modes, Quickly

The ITQ 1100 comes standard with a vacuum probe interlock system. This system allows for insertion and removal of the "ion volume" in the source without venting the analyzer. Removable ion volumes provide added versatility and convenience for all types of analysis. With the vacuum interlock option, you can exchange the ion volume without breaking vacuum. Replace dirty ion volumes quickly to restore response factors. With optional chemical ionization, you can switch between EI and CI and back again, without venting the analyzer to make this change. Optimizing data quality for EI and CI work is further enhanced through the use of specially designed ion volumes, each developed to provide optimal performance depending on the ionization mode.


Data Dependent Scanning Facilitates Method Development

Data dependency is a powerful and time-saving tool for dynamic and complex qualitative analyses. Data dependent algorithms probe mass spectra for ions matching given mass selection settings. The list of ions is compiled in an intensity-sorted dependent mass list that is used for MS/MS scans in subsequent events. Data dependency allows MSⁿ methods to be written without knowing in advance what specific precursor ions may be found during an analysis. Additional data are obtained without spending time running multiple methods. You benefit from library-searchable full-scan MS spectra and additional structural information from the MSⁿ data.




Flame Retardant in Pond Water


4 ppt 2,2',4,4',6-Pentabromodiphenyl ether in pond water




Full-scan spectrum of 4 fg/μL 2,2',4,4',6-PBDE spiked into pond water.



MS/MS product ion of 2,2',4,4',6-PBDE spiked into pond water.



MS/MS product ion of 2,2',4,4',6-PBDE spiked into benzene.



Theoretical product ion cluster of a pentabromodiphenyl ether.

403.78
405.78
407.78
408.78
409.78
410.78
411.78
412.78
413.78
414.78
415.78
416.78
417.78
418.78
419.78
420.78
421.78
422.78
423.78
424.78
425.78
426.78
427.78
428.78
429.78
430.78
431.78
432.78
433.78
434.78
435.78
436.78
437.78
438.78
439.78
440.78
441.78
442.78
443.78
444.78
445.78
446.78
447.78
448.78
449.78
450.78
451.78
452.78
453.78
454.78
455.78
456.78
457.78
458.78
459.78
460.78
461.78
462.78
463.78
464.78
465.78
466.78
467.78
468.78
469.78
470.78
471.78
472.78
473.78
474.78
475.78
476.78
477.78
478.78
479.78
480.78
481.78
482.78
483.78
484.78
485.78
486.78
487.78
488.78
489.78
490.78
491.78
492.78
493.78
494.78
495.78
496.78
497.78
498.78
499.78
500.78
501.78
502.78
503.78
504.78
505.78
506.78
507.78
508.78
509.78
510.78
511.78
512.78
513.78
514.78
515.78
516.78
517.78
518.78
519.78
520.78
521.78
522.78
523.78
524.78
525.78
526.78
527.78
528.78
529.78
530.78
531.78
532.78
533.78
534.78
535.78
536.78
537.78
538.78
539.78
540.78
541.78
542.78
543.78
544.78
545.78
546.78
547.78
548.78
549.78
550.78
551.78
552.78
553.78
554.78
555.78
556.78
557.78
558.78
559.78
560.78
561.78
562.78
563.78
564.78
565.78
566.78
567.78
568.78
569.78
570.78
571.78
572.78
573.78
574.78
575.78
576.78
577.78
578.78
579.78
580.78
581.78
582.78
583.78
584.78
585.78
586.78
587.78
588.78
589.78
590.78
591.78
592.78
593.78
594.78
595.78
596.78
597.78
598.78
599.78
600.78
601.78
602.78
603.78
604.78
605.78
606.78
607.78
608.78
609.78
610.78
611.78
612.78
613.78
614.78
615.78
616.78
617.78
618.78
619.78
620.78
621.78
622.78
623.78
624.78
625.78
626.78
627.78
628.78
629.78
630.78
631.78
632.78
633.78
634.78
635.78
636.78
637.78
638.78
639.78
640.78
641.78
642.78
643.78
644.78
645.78
646.78
647.78
648.78
649.78
650.78
651.78
652.78
653.78
654.78
655.78
656.78
657.78
658.78
659.78
660.78
661.78
662.78
663.78
664.78
665.78
666.78
667.78
668.78
669.78
670.78
671.78
672.78
673.78
674.78
675.78
676.78
677.78
678.78
679.78
680.78
681.78
682.78
683.78
684.78
685.78
686.78
687.78
688.78
689.78
690.78
691.78
692.78
693.78
694.78
695.78
696.78
697.78
698.78
699.78
700.78
701.78
702.78
703.78
704.78
705.78
706.78
707.78
708.78
709.78
710.78
711.78
712.78
713.78
714.78
715.78
716.78
717.78
718.78
719.78
720.78
721.78
722.78
723.78
724.78
725.78
726.78
727.78
728.78
729.78
730.78
731.78
732.78
733.78
734.78
735.78
736.78
737.78
738.78
739.78
740.78
741.78
742.78
743.78
744.78
745.78
746.78
747.78
748.78
749.78
750.78
751.78
752.78
753.78
754.78
755.78
756.78
757.78
758.78
759.78
760.78
761.78
762.78
763.78
764.78
765.78
766.78
767.78
768.78
769.78
770.78
771.78
772.78
773.78
774.78
775.78
776.78
777.78
778.78
779.78
780.78
781.78
782.78
783.78
784.78
785.78
786.78
787.78
788.78
789.78
790.78
791.78
792.78
793.78
794.78
795.78
796.78
797.78
798.78
799.78
800.78
801.78
802.78
803.78
804.78
805.78
806.78
807.78
808.78
809.78
810.78
811.78
812.78
813.78
814.78
815.78
816.78
817.78
818.78
819.78
820.78
821.78
822.78
823.78
824.78
825.78
826.78
827.78
828.78
829.78
830.78
831.78
832.78
833.78
834.78
835.78
836.78
837.78
838.78
839.78
840.78
841.78
842.78
843.78
844.78
845.78
846.78
847.78
848.78
849.78
850.78
851.78
852.78
853.78
854.78
855.78
856.78
857.78
858.78
859.78
860.78
861.78
862.78
863.78
864.78
865.78
866.78
867.78
868.78
869.78
870.78
871.78
872.78
873.78
874.78
875.78
876.78
877.78
878.78
879.78
880.78
881.78
882.78
883.78
884.78
885.78
886.78
887.78
888.78
889.78
890.78
891.78
892.78
893.78
894.78
895.78
896.78
897.78
898.78
899.78
900.78
901.78
902.78
903.78
904.78
905.78
906.78
907.78
908.78
909.78
910.78
911.78
912.78
913.78
914.78
915.78
916.78
917.78
918.78
919.78
920.78
921.78
922.78
923.78
924.78
925.78
926.78
927.78
928.78
929.78
930.78
931.78
932.78
933.78
934.78
935.78
936.78
937.78
938.78
939.78
940.78
941.78
942.78
943.78
944.78
945.78
946.78
947.78
948.78
949.78
950.78
951.78
952.78
953.78
954.78
955.78
956.78
957.78
958.78
959.78
960.78
961.78
962.78
963.78
964.78
965.78
966.78
967.78
968.78
969.78
970.78
971.78
972.78
973.78
974.78
975.78
976.78
977.78
978.78
979.78
980.78
981.78
982.78
983.78
984.78
985.78
986.78
987.78
988.78
989.78
990.78
991.78
992.78
993.78
994.78
995.78
996.78
997.78
998.78
999.78
1000.78

ITQ 1100 Highlights

Standard Features

- External ionization ion source
- Scan range of 10–1100 u
- Most sensitive EI full scan GC/ion trap MS
- 250 L/s turbomolecular pump
- MSⁿ, for MS/MS experiments
- Automated Collision Energy
- Pulsed Q Dissociation Mode
- Vacuum Probe Interlock – for easy source maintenance without venting
- Intelligent Data Dependent Scanning
- Powerful Xcalibur data system – for instrument operation, data analysis and reporting
- TRACE 1310 GC – high-performance GC with user-exchangeable plug-in injector and detector technology and an easy-to-use touch screen interface

Upgrade Options


- Chemical Ionization, including PPNICI – expand your analytical capabilities
- Direct sample probe – analyze samples without chromatographic separation and without venting or moving the GC

Added Value Software Options

- TraceFinder software follows your lab's workflows
- Mass Frontier™ software – for structural elucidation and fragmentation studies



Advanced options boost performance and productivity



Interchangeable solids probe tips:
(1) Direct Exposure-DEP
(2) Direct Insertion-DIP

Solids Probe

Solids probes are most often used for qualitative or semi-quantitative analysis of materials that are difficult, if not impossible to elute chromatographically.

With a GC injection, polar high-molecular weight compounds will often decompose during the volatilization process leaving the analyst with a spectrum containing the desired compound as well as the decomposition products. Our Direct Exposure Probe (DEP) and Direct Insertion Probe (DIP) are especially suited to analysis of thermally labile compounds. Both are controlled by a single, easy-to-use controller.

Most customers want a simple interface that allows analysis of a solid sample without disturbing the GC setup. With an ITQ system, probe analysis can be performed without disconnecting any part of the gas chromatograph or breaking vacuum.

- Ideal for thermally labile samples
- Convenient, interchangeable probe types




ITQ 1100 GC/MS with direct sample probe

Variable Damping Gas Control

Only the ITQ Series includes a variable damping gas option. This allows you to increase the amount of helium pressure in the ion trap independent of the column carrier gas. The increased helium pressure traps more ions per unit time thereby increasing detectivity, especially in MSⁿ mode. This option can increase the detected area of compounds 500% or more.

Thermo Scientific TriPlus RSH Autosampler – Flexible Sampling Solutions

The TriPlus™ RSH autosampler provides unrivaled performance in liquid sampling automation, with tremendous productivity for high-throughput labs.



Unique, robotic tool change capabilities enable seamless automated switching from liquid to headspace to SPME injection – all in the same sequence.

Thermo Scientific AI/AS 1310 Series Autosampler – Affordable, Powerful Simplicity


The AI/AS 1310 Series liquid autosampler combines utmost precision in liquid sampling with simple, self-configuring setup that ensures confident performance in high-throughput and research environments.

AI/AS 1310 Series Options:

- AI 1310 Autoinjector – Plug & Inject 8-position modular sampling system
- AS 1310 Autosampler – 105-position sample tray with integrated control
- AI/AS 1310 Gemini Kit – Inject into two injection ports on the TRACE GC Ultra simultaneously



DEP Probe Analysis of TNT using PPINICI



Full-scan Direct Exposure Probe-PPINICI analysis of 2,4,6-TNT


Comprehensive software solutions

intuitive platforms for GC/MS, LC/MS and advanced MS instruments provide confident control from method development to reporting

Maximum Versatility

Thermo Scientific Xcalibur Data System

The Xcalibur™ data system is the core 'operating system' for each of the broad range of Thermo Scientific mass spectrometry systems. Offering a common set of easy-to-use, yet powerful, tools for GC/MS and LC/MS alike, Xcalibur software provides a unified experience for every user of every system. It offers instrument control, sample sequencing, and a set of programs for both qualitative and quantitative applications. Xcalibur software is also compatible with commercially-available mass spectral libraries, such as NIST, Wiley and Maurer-Pfleger-Weber, and others. It also includes tools for generating and maintaining your own spectral libraries. Supporting



each type of autosampler and GC available, Xcalibur software allows each mass spectrometer to be utilized to its fullest.


Optimized Productivity

Thermo Scientific TraceFinder Software

The TraceFinder™ software application provides a streamlined workflow for the needs of a wide variety of high-throughput quantitative applications. Customized versions target the critical needs for key application areas, such as environmental, food safety, clinical research and forensic toxicology.

Beginning at the dashboard and working through each of the specialized modules, TraceFinder software guides the user through a highly productive workflow that is focused on the fastest route from sample to report.

Automated processing and QA/QC review, coupled to a reporting system capable of both standard and customized hardcopy and electronic deliverables, this software provides a comprehensive system for high-productivity quantitation.



Easy Interpretation of Mass Spectral Data

Thermo Scientific Mass Frontier software provides a suite of options that allows the MS operator to interpret mass spectral data. This optional software package assists the operator in structural identification, isotope pattern determination, and spectral classifications.



Thermo Scientific chromatography columns and consumables

the perfect partner for optimal analytical performance

Thermo Scientific chromatography columns and consumables are designed to complement our innovative range of GC and GC-MS systems. Get the most out of the TSQ 8000 GC-MS system by pairing it with advanced, high-performance Thermo Scientific products. The wide range of consumables and accessories offer customers applications-focused solutions in the environmental, food analysis, forensics/toxicology, petrochemical, pharmaceutical and general analytical industries.



Sample Handling

- Comprehensive range of vials and closures
- Mass Spec Certified vials, the only pre-cleaned chromatography vial available

Sample Preparation

- Thermo Scientific HyperSep Retain polymeric solid phase extraction products
- QuEChERS solutions for efficient sample preparation and clean-up
- GC derivatization reagents
- Reacti-Therm sample derivatization system

GC Columns and Accessories

- Thermo Scientific TraceGOLD GC columns
- Thermo Scientific TracePLOT columns – for permanent gas and light hydrocarbon
- GC Consumables
- Essential tools and equipment, including the following:
 - GLD-Pro Leak Detector
 - GFM-Pro Electronic Flowmeter



www.thermoscientific.com

©2012 Thermo Fisher Scientific Inc. All rights reserved. ISO is a trademark of the International Standards Organization. iPad is a registered trademark of Apple Inc. Mass Frontier is a trademark of HighChem Ltd. Microsoft and Windows are registered trademarks of Microsoft Corporation. Supelco is a registered trademark of Sigma-Aldrich Inc. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.



Africa-Other +27 11 570 1840	Europe-Other +43 1 333 50 34 0	Japan +81 45 453 9100	Spain +34 914 845 965
Australia +61 3 9757 4300	Finland/Norway/Sweden +46 8 556 468 00	Latin America +1 561 688 8700	Switzerland +41 61 716 77 00
Austria +43 1 333 50 34 0	France +33 1 60 92 48 00	Middle East +43 1 333 50 34 0	UK +44 1442 233555
Belgium +32 53 73 42 41	Germany +49 6103 408 1014	Netherlands +31 76 579 55 55	USA +1 800 532 4752
Canada +1 800 530 8447	India +91 22 6742 9434	New Zealand +64 9 980 6700	
China +86 10 8419 3588	Italy +39 02 950 591	Russia/CIS +43 1 333 50 34 0	
Denmark +45 70 23 62 60		South Africa +27 11 570 1840	

BR10234_E 07/12M

Thermo
SCIENTIFIC
Part of Thermo Fisher Scientific