



Catalog

Electrochemistry accessories.



A one-stop-shop to suit your every need.

As one of the world's first designers and manufacturers of high-performance electrochemical measurement instruments, BioLogic has forged its place in the international market.

Our comprehensive product portfolio covers cutting-edge scientific products in such diverse applications as electrochemistry, battery testing and fuel cell/material testing.

And supporting these high-precision measurement tools is an extensive line of product accessories ranging from sophisticated quartz crystal microbalances to electrodes.

So whatever your research field, we can provide you with high-quality, hard-wearing equipment.

BioLogic: a one-stop-shop to suit your every need.



Analytical Cells	Small Volume Cells	4
	Large Volume Cells	6
	Flow Cells	7
Corrosion Cells	Standard Corrosion Cells	8
	Avesta Cell	9
	Flat Cells	10
	Galvanic Cells	11
	Plate Material Evaluating Cell	11
	Coating Cell	11
Glassware	Small & Large Volume Cell Vials	12
	Bridge & Purge Tubes for Corrosion Cells	14
	Small-Size Bridge Tubes	14
	Bridge & Purge Tubes for Analytical Cells	15
	Bridge Tubes for King-Size Reference Electrodes	15
Electrodes	Working Electrodes	16
	King-Size Reference Electrodes	18
	Counter Electrodes	19
	Metallic Electrodes	19
	Small-Size Reference Electrodes	20
	Maintenance of Reference Electrodes	22
	Screen Printed Electrodes	23
	Inter Digital Array (IDA) Electrodes	23
	Ring-Disk Type Electrodes	23
Rotating Electrodes	BluRev Rotating Disk Electrode (RDE)	24
	Rotating Ring-Disk Electrode (RRDE)	26
Spectroelectrochemistry	Spectrometer System	28
	Static Cell	30
	Flow Cell	31
	Thin Layer Cells	32
Quartz Crystal Microbalance	BluQCM QSD	34
	QCA922A	36
Battery Accessories	Coin Cell Holders	38
	Pouch Cell Holders	38
	Cylindrical Cell Holders	39
	Prismatic & Pouch Cell Holders	39
	Current Collectors	40
	Sense Adapter Module (SAM-50)	40
	Redox Flow Battery Cells	41
Connection Accessories	External Device Connection	42
	Test Boxes & Dummy cells	42
	Glove Box & Long Cables	43
	Connectors	43
	Multi-Electrode Investigation Cables	44
	Faraday Cages	44
Scanning Products Accessories	Probes	45
	Video Microscope System	45
	Cells	46
	M470 Glove Box Cables	47
	USB-PIO	47
Material Testing Accessories	Temperature Control Units	48
	Sample Holders	49

Analytical Cells.



SVC-2



SVC-3

Small Volume Cells

Each voltammetry cell is designed for a specific application (specific working electrode, volume of solution, oxygen-free condition, etc.).


For example, for voltammetry investigations using standard working electrodes with an outer diameter (OD) of 6 mm, the fixed configuration of SVC-3 kit is recommended.

For applications requiring other working electrode shapes, SVC-2 is more suitable.



If only a small amount of the electroactive compound is available, SVC-2 in the microvolume mode is recommended. Here is a list of the available cells:

- SVC-2, modular
- SVC-3, for a volume of 5 to 20 mL, only for working electrode with OD of 6 mm
- VC-4, for a volume of 1 to 3 mL, only for working electrode with OD of 6 mm
- Bulk electrolysis cell, for a volume of 100 mL

SVC-2 modularity:
microvolume mode



oxygen free mode



Please note that a full, purpose-built, analytical kit is also available SK-2 (A-012763) but the reference electrode must be purchased separately.

This kit includes:

- SVC-3 kit (A-012669)
- PK-3 polishing kit (A-011975) see page 16
- one glassy carbon electrode, OD 6.0 mm, ID 3.0 mm (A-002012) see page 17
- one platinum electrode, OD 6.0 mm, ID 1.6 mm (A-002013) see page 17



VC-4






Bulk electrolysis cell

Small Volume Cells

Products	SVC-2	A-012668	SVC-3	A-012669	VC-4	A-011224	Bulk electrolysis cell	A-001197
Content								
Sample vial/mL	20 (7 pieces)	A-001056	20 (7 pieces)	A-001056	5 (7 pieces)	A-011504	100 (1 piece)	A-012632
Counter electrode (CE)/mm	57	A-002222	50	A-002222	57	A-002222	230	A-002234
PTFE cap		A-012670		A-012671		A-011226		A-O 12551
Purge tube (ETFE), 100 mm		-		-		-		-
Additional items	Adapter 10 to 6 mm	-			Cell holder	A-011227	Porous carbon electrode	A-010530
							Lid for CE	A-001198
							Chamber for CE	A-001196
							O-ring	A-001236
							Port plug	A-009131
							Stirrer bar	A-000178
Options								
Sample holder/mm	9.0 (2 pieces)	A-012177	6.0 (2 pieces)	A-012176				
Cell holder	for 20 mL	A-001209	for 20 mL	A-001209				
Purge tube (ETFE)/m	1	A-010537	1	A-010537	1	A-010537	1	A-010537
Working electrodes	See page 17		See page 17		See page 17		See page 18	
Reference electrodes	See page 18		See page 18		See page 18			

PTFE cap



Cell geometry

The geometry of the cell should be optimized to reduce the ohmic drop. Working and reference electrodes must be positioned close to one another. The counter electrode should not limit the transfer of electrons, so its contact surface should be larger than the contact surface of the working electrode.

Large Volume Cells



An analytical cell available in sizes ranging from 50 to 250 mL that is particularly well-suited to larger volumes of electrolyte.

Two types of packages are available:

- A standard analytical cell kit
- A full analytical cell kit (which allows temperature control and gas purging)

	Standard analytical cell kit (80 mL) EL-ELECTRO-1	Complete analytical cell kit (80 mL) EL-ELECTRO-2	Complete analytical cell kit (150 mL) EL-ELECTRO-3	Catalog n°
Glass cell 80 mL	■			EL-A-001
Double jacketed glass cell 80 mL		■		EL-A-002
Double jacketed glass cell 150 mL			■	EL-A-020
PTFE cap 5 holes	■	■	■	EL-A-003
PTFE ring, silicon encapsulated, OD 102 mm	■	■	■	EL-A-004
Cell collar with clamp	■	■	■	EL-A-005
Double purge tube		■	■	EL-A-006
Bridge tube for reference electrode, OD 6 mm	■	■	■	EL-A-008
Platinum counter electrode	■	■	■	EL-A-009
Purge tube	■			EL-A-016
Reference electrode RE-2BP Hg/Hg ₂ Cl ₂ (Cells kits following RoHS recommendations are also provided, see below)	■	■	■	A-013430
Double nut 25 mm and 12 mm diameter		■	■	EL-A-011
Telescopic cell stand		■	■	EL-A-012
Options				
Electrode bridge extension for electroanalytical cell				EL-A-022
Bridge tube for reference electrode of OD 8 mm				EL-A-017
PT100 probe, indicate connector type				EL-C-014
Magnetic stirrer & heater, 220 V				EL-C-015A
without PT100 probe 110 V				EL-C-015B
Aluminum base holder for magnetic stirrer				EL-C-018
Set of 10 porous 4 mm glass frits (CoralPor™) with PTFE heat shrink (200 mm)				092-VYC4
Optional cells kits following RoHS recommendations				
EL-ELECTRO-1 cell with Ag/AgCl reference electrode				EL-ELECTRO-1A
EL-ELECTRO-2 cell with Ag/AgCl reference electrode				EL-ELECTRO-2A
EL-ELECTRO-3 cell with Ag/AgCl reference electrode				EL-ELECTRO-3A
EL-ELECTRO-1A compatible with BluRev RDE (comes with 094-A-CAP an additional PTFE cap compatible with BluRev)				EL-BLUREV



EL-A-012



EL-A-002

EL-A-003


EL-A-004

EL-A-005

EL-A-008

EL-A-009

EL-A-022



Cross flow cell



Radial flow cell

Flow Cells

	Cross Flow Cell	Catalog n°	Radial Flow Cell	Catalog n°
Product	Cross flow cell*	A-012798	Radial flow cell*	A-012799
Schematic diagram				

* Gasket, working, reference electrodes and connectors are sold separately

Options

Single/radial flow	Electrode	Glassy carbon	Ø 3 mm x 2 (dual type), size 25 x 25 mm	A-001000
		Gold	Ø 3 mm x 2 (dual type), size 25 x 25 mm	A-001002
Cross flow	Electrode	Platinum	Ø 3 mm x 2 (dual type), size 25 x 25 mm	A-001012
		Silver	Ø 3 mm x 2 (dual type), size 25 x 25 mm	A-001008
Dual, series/cross flow	Electrode	Carbon paste	Ø 3 mm x 2 (dual type), size 25 x 25 mm	A-001004
		Nickel	Ø 3 mm x 2 (dual type), size 25 x 25 mm	A-001009
Dual, parallel/cross flow	Electrode	Glassy carbon/gold	Ø 3 mm x 2 (dual type), size 25 x 25 mm	A-001006
		Glassy carbon/platinum	Ø 3 mm x 2 (dual type), size 25 x 25 mm	A-012583
Cross flow	Gasket	TG-2M PTFE (4 pieces)	Film thickness 12 µm	A-001046
		TG-5M PTFE (4 pieces)	Film thickness 25 µm	A-001047
Radial flow	Gasket	TG-6M PTFE (4 pieces)	Film thickness 50 µm	A-001048
		TG-8M PTFE (4 pieces)	Film thickness 100 µm	A-001049
Dual, parallel/cross flow	Gasket	TG-11M silicon (4 pieces)	Film thickness 500 µm	A-001092
		TG-12M silicon (4 pieces)	Film thickness 1,000 µm	A-001093
Cross flow	Cross flow cell block		Material PEEK	A-001032
Radial flow	Electrode	Glassy carbon	Ø 6 mm (single type), size 25 x 25 mm	A-001016
		Glassy carbon	Ø 3 mm (single type), size 25 x 25 mm	A-012124
Dual, parallel/cross flow	Electrode	Platinum	Ø 3 mm (single type), size 25 x 25 mm	A-009908
		Gold	Ø 3 mm (single type), size 25 x 25 mm	A-011155
Radial flow	Electrode	PFCE (Plastic Formed Carbon Electrode)	Ø 3 mm (single type), size 25 x 25 mm	A-000999
		Carbon paste	Ø 3 mm (single type), size 25 x 25 mm	A-010251
Radial flow	Gasket	TG-2MR PTFE (4 pieces)	Film thickness 12 µm	A-001146
		TG-5MR PTFE (4 pieces)	Film thickness 25 µm	A-001147
Dual, parallel/cross flow	Gasket	TG-6MR PTFE (4 pieces)	Film thickness 50 µm	A-001148
		TG-8MR PTFE (4 pieces)	Film thickness 100 µm	A-012802
Radial flow	Gasket	TG-11MR PTFE (4 pieces)	Film thickness 500 µm	A-001192
		TG-12MR PTFE (4 pieces)	Film thickness 1,000 µm	A-001193
Cross/radial	Radial flow cell block		Material PEEK	A-001031
Cross/radial	RE-3VT Ag/AgCl reference electrode screw type		Size Ø 10 x 48 mm	A-013488
		RE-7VT Non-aqueous Ag/Ag ⁺ reference electrode screw type	Size Ø 10 x 48 mm	A-013489
Cross/radial	0.04" Single lead connector (2 pieces)			A-012912
		0.04" MM connector (10 pieces)		A-013273
Cross/radial	1/16" Peek tube		ID 0.25 mm, length 3.0 m	A-001531
		Dynasael peek fingertight (10 pieces)	Screw for pipe connecting integrated 1/16"	A-004130

Corrosion Cells.



Standard Corrosion Cells

For corrosion applications, two packages for large volumes i.e. ~1 L are available:

- A basic corrosion cell kit
- A complete corrosion cell kit (allows temperature control and includes sample holder and cell stand)



	Basic corrosion cell kit EL-CORR-1	Complete corrosion cell kit EL-CORR-2	Catalog n°
Glass cell 1 L	■		EL-C-001
Double jacketed glass cell 1 L		■	EL-C-002
PTFE cap	■	■	EL-C-003
PTFE ring, silicon encapsulated, OD 102 mm	■	■	EL-C-004
Cell collar with clamp	■	■	EL-C-005
Double purge tube		■	EL-C-006
Graphite counter electrode rod (2 pieces) $\rho = 1.070 \mu\Omega \text{ cm}$	■	■	EL-C-009
Double nut 25 mm and 12 mm diameter		■	EL-C-011
Telescopic cell stand		■	EL-C-012
Sample holder 1 cm ² (max sample thickness 3.4 mm and max diameter 14.6 mm)		■	092-C-016
Purge tube	■		EL-C-016
Bridge tube for 8 mm diameter reference electrode	■	■	EL-C-017
Calomel reference electrode length 80 mm, OD 8 mm screw cap	■	■	R-XR300
Cable connection for screw cap 100 mm, 2 mm banana plug	■	■	R-A94L111
12 mm OD conical rings for reference electrode of 8 mm (4 pieces)	■	■	R-X31M012
Options			
Bridge tube for 6 mm diameter reference electrode			EL-C-008
PT100 probe (indicate connector type)			EL-C-014
Magnetic stirrer & heater, without PT100 probe		220 V	EL-C-015A
		110 V	EL-C-015B
Aluminum base holder for magnetic stirrer and 1 L cell vial			EL-C-018
Set of 10 porous frits (4 mm CoralPor™) with PTFE heat shrink (200 mm)			092-VYC4
500 mL double jacketed glass cell is also available (EL-C-020).			



Avesta Cell

The Avesta Cell is an electrochemical cell developed for pitting corrosion testing (ASTM G150).



It is designed to avoid microcrevice corrosion formed between the working electrode and the gasket at the bottom aperture of the cell.

A filter paper ring placed between the sample and the gasket is flooded by distilled water in order to eliminate crevice corrosion.

The water flow is controlled by a peristaltic pump (EL-AV-008). The flow rate depends on the chosen tu- bings, contact your local retailer for more information.

Avesta Cell	Catalog n°
Avesta cell kit	EL-AV-001
Content	
Double jacketed cell glass 250 mL	EL-AV-002
PTFE cap 5 holes	EL-AV-003
O-ring PTFE silicone encapsulated	EL-AV-004
Cell collar with clamp	EL-A-005
Double purge tube	EL-AV-006
Filter paper ring (100 pieces)	EL-AV-007
Graphite counter electrode rods (2 pieces) $\rho = 1.070 \mu\Omega \text{ cm}$	EL-C-009
Bridge tube for RE 6 mm	EL-AV-010
Double nut	EL-C-011
Telescopic cell stand	EL-AV-012
Skeleton	EL-AV-015
Options	
Peristaltic pump for low flow	EL-AV-008
Bridge tube for reference electrode with OD of 8 mm	EL-AV-013
Single purge tube	EL-AV-014
Temperature probe PT100	EL-C-014
Set of 10 porous frits (4 mm CoralPor™) with PTFE heat shrink (200 mm)	092-VYC4

Principle of Avesta cell



Corrosion Cells.



Flat Cells, 1 to 10 cm² sample area

Sample size:
30.5 mm min. diameter (EL-FLAT)
40.5 mm min. diameter (EL-FLAT-2)
6 mm max. thickness

This cell with a volume of 250 mL is suited to perform experiments on flat specimens of 1 or 10 cm² surface area.

This cell has a double jacket for temperature control and three holes for reference electrodes, purge tubes and temperature probes with an inner diameter of 17.6 mm and the two others of 8.3 mm.

Flat Cell 1 cm ²	Catalog n°
Flat cell kit 1 cm ²	EL-FLAT
Content	
RE-2BP calomel reference electrode (length: 90 mm, OD: 6 mm)	A-013430
Platinum mesh counter electrode (54 mm wire/80 mesh), 25 x 35 mm	A-702439
Reference electrode bridge tube (6 mm diameter)	EL-F-004
Corrosion flat cell 1cm ²	EL-FLAT-3
EL-FLAT-3 includes:	
- glass part for flat cell (250 mL)	EL-F-002
- mechanical part for corrosion flat cell 1 cm ²	EL-FLAT-3H

Flat Cell 10 cm ²	Catalog n°
Flat cell kit 10 cm ²	EL-FLAT-2
Content	
RE-2BP calomel reference electrode (length: 90 mm, OD: 6 mm)	A-013430
Platinum mesh counter electrode (54 mm wire/ 80 mesh), 25 x 35 mm	A-702439
Reference electrode bridge tube (6 mm diameter)	EL-F-004
Corrosion flat cell 10 cm ²	EL-FLAT-4
EL-FLAT-4 includes:	
- glass part for flat cell (250 mL)	EL-F-002
- mechanical part for corrosion flat cell 10 cm ²	EL-FLAT-4H

Investigations in aggressive media

If the experiment is performed in more aggressive media such as fluorhydric acid, it is possible to obtain the body of the flat cell in PVDF* instead of glass.

	Fluorhydric acid 48%		Sulfuric acid 98%		Phosphoric acid 85%		Hydrochloric acid 35%		Nitric acid 70%		Perchloric acid		Sodium hydroxide 50%		Potassium hydroxide concentr.	
Temperature	20° C	50° C	20° C	50° C	20° C	50° C	20° C	50° C	20° C	50° C	20° C	50° C	20° C	50° C	20° C	50° C
PTFE	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
PVDF*	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
Borosilicated glass	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■

Options	Catalog n°
PVDF* body (single-jacketed)	EL-F-PVDF
EPDM O-rings for 1 cm ² with PEEK ferrule for CE**	EL-SEAL-1B
EPDM O-rings for 10 cm ² with PEEK ferrule for CE**	EL-SEAL-10B
PTFE O-rings for 1 cm ² with PEEK ferrule for CE**	EL-SEAL-T1B
PTFE O-rings for 10 cm ² with PEEK ferrule for CE**	EL-SEAL-T10B
EPDM O-rings for 1 cm ² **	EL-SEAL-1
EPDM O-rings for 10 cm ² **	EL-SEAL-10
PTFE O-rings for 1 cm ² **	EL-SEAL-T1
PTFE O-rings for 10 cm ² **	EL-SEAL-T10
Set of 10 porous frits (4 mm CoralPor™) with PTFE heat shrink (200 mm)	092-VYC4

* Polyvinylidene fluoride

** The O-ring kits include 4 O-rings for the glassware side and 10 O-rings for the sample side

■ excellent
■ good
■ not recommended
■ not compatible

Sample size:
30.5 mm min. diameter (EL-GAL-1)
40.5 mm min. diameter (EL-GAL-10)
6 mm max. thickness



Galvanic Cells, 1 to 10 cm² sample area

Thanks to the modular design of the flat cell, it is possible to place two different materials at each end of the cell. The surface area may be 1 or 10 cm².

Galvanic Cell 1 cm²	Catalog n°
Galvanic cell kit 1 cm ²	EL-GAL-1
Content	
Flat cell kit 1 cm ²	EL-FLAT
Galvanic kit 1 cm ²	092-FLAT/1

Galvanic Cell 10 cm²	Catalog n°
Galvanic flat cell kit 10 cm ²	EL-GAL-10
Content	
Flat cell kit 10 cm ²	EL-FLAT-2
Galvanic kit 10 cm ²	092-FLAT/10



Plate Material Evaluating Cell, up to 1 cm² sample area

This cell was developed to evaluate a plate material such as metals, semi-conducting plates, etc. The sample plate is sandwiched between the two cell blocks. The required volume of solution is about 1 mL.

Plate Material Evaluating Cell	Catalog n°
Plate material evaluating cell	A-011951
Content	
PTFE cell [body & base] (1 piece)	-
O-ring (1 piece)	-
Screw 20 mm (1 piece)	-
Purging tube, 100 mm	-
Platinum counter electrode (1 piece)	A-002222
Options	
O-ring (10 pieces)	A-012022



Coating Cell

Coating Cell	Catalog n°
Coating cell kit	EL-COAT
Content	
Glass for coating cell	EL-P-002
Nylon base with three feet	EL-P-003
Rubber cup with two holes	EL-P-004
Metallic clamp	EL-P-005
O-ring for coating cell	EL-P-006
Graphite rod counter electrode (L: 145 mm, OD: 6 mm, ρ = 1.070 μΩ cm)	EL-P-009
Options	
RE-1B Ag/AgCl reference electrode (OD: 6 mm)	A-012167
Bridge tube for 6 mm reference electrode	EL-P-008
Mask for 1 cm ² (20 pieces)	EL-P-011
Mask for 3 cm ² (20 pieces)	EL-P-012
Mask for 10 cm ² (20 pieces)	EL-P-013

Small Volume Cell Vials up to 200 mL

To complete a kit, the cell vials are also available separately.
Please note that other volumes are available on request.



Small & Large Volume Cell Vials		Volume/mL	OD/mm	ID/mm	Height/mm	Quantity	Purpose	Catalog n°
Small Volume Cell Vials	Simple vial	5	18	15.6	30	10	VC-4	A-011504
		20	28	25.6	50	10	SVC-2, SVC-3	A-001056
		100	50	46.4	72	1	RRDE-3A, bulk electrolysis	A-012632
	Water jacketed-glass vial	5	40	15.6	40	1	VC-4	A-012672
		20	55	25.6	50	1	SVC-2, SVC-3	A-001051
		100	70	46.4	80	1	RRDE-3A, bulk electrolysis	A-012652
	Simple vial for alkaline solution	100	51.5	46.5	72	10	RRDE-3A, bulk electrolysis	A-013580
		200	67	62	72	8	RRDE-3A, bulk electrolysis	A-013581
	Cell holder for 20 mL vial	-	-	-	-	1	SVC-2, SVC-3	A-001209
	PTFE cap for 100 mL cells	-	-	-	-	1	For A-012632, A-012652, A-013580	A-012631
Large Volume Cell Vials	Simple vial	-	-	-	-	1	For A-013581	A-013582
		80	90	62	80	1	EL-ELECTRO-1	EL-A-001
		250	-	-	-	1	EL-FLAT, EL-FLAT-2, EL-GAL-1, EL-GAL-10	EL-F-002
		1,000	120	90	175	1	EL-CORR-1	EL-C-001
	Water jacketed-glass vial	80	90	62	85	1	EL-ELECTRO-2	EL-A-002
		150	-	-	-	1	EL-ELECTRO-3	EL-A-020
		1,000	120	90	200	1	EL-CORR-2	EL-C-002
		2,000	-	-	-	1	EL-CORR-2	EL-C-021
	Cell holder for king-size vial (1 L)	-	-	-	-	1	EL-CORR-1	EL-C-018

Large Volume Cell Vials up to 2,000 mL



EL-A-001



EL-A-002



EL-C-001

EL-C-018,
cell holder for EL-C-001




EL-C-002

Dimensions in mm

Bridge & Purge Tubes
for corrosion cells

The tolerance of each dimension is approximately ± 0.5 mm.



Dimensions in mm

Small-Size Bridge Tubes



Bridge & Purge Tubes for analytical cells

Bridge Tubes for king-size reference electrodes



Dimensions in mm

Glassware		OD/ mm	Height/ mm	Quan- tity	Purpose	Vycor type compati- bility	Catalog n°
Small-size bridge tube	Ø 6 mm	6	68	2	SVC-2, SVC-3	092-VYC3⁽¹⁾	A-012176
	Ø 6 mm	6	68	22	SVC-2, SVC-3		A-012306
	Ø 9 mm	9	68	2	SVC-2		A-012177
	Ø 9 mm	9	68	22	SVC-2		A-012307
Bridge tube for corrosion cells	Ø 6 mm	10	250	1	EL-CORR	092-VYC4⁽²⁾	EL-C-008
	Ø 8 mm	10	250	1	EL-CORR		EL-C-017
	Ø 6 mm	10	120	1	EL-FLAT		EL-F-004
	Ø 8 mm	10	140	1	EL-FLAT		EL-F-004B
Bridge tube for analytical cells	Ø 8 mm	10	165	1	EL-ELECTRO		EL-A-017
	Ø 6 mm	10	155	1	EL-ELECTRO		EL-A-008
Purge tube for analytical cells	Single	10	165	1	EL-ELECTRO		EL-A-016
	Double	10	200	1	EL-ELECTRO		EL-A-006
Purge tube for corrosion cells	Single	10	200	1	EL-CORR		EL-C-016
	Double	10	220	1	EL-CORR		EL-C-006
Bridge tube with ceramic junction	Short	12	70	1			R-AL100
	Standard	8	250	1			R-AL110
	With reverse sleeve (non aqueous)	8	138	1			R-AL210
	Standard	8	138	1			R-AL120
Options							
(1) Set of 10 porous 2.8 mm glass frits (CoralPor™) with PTFE heat shrink (200 mm)							092-VYC3
(2) Set of 10 porous 4 mm glass frits (CoralPor™) with PTFE heat shrink (200 mm)							092-VYC4


Working Electrodes

To address every application, a wide range of working electrodes (WE) is available with diameters ranging from 7 μm up to 6 mm and built with different materials.

BioLogic exclusive



M-BDD-3: Boron-doped diamond 3 mm diameter disk




Dimensions in mm

Maintenance of working electrodes

To refresh the electrode surface, we recommend polishing it before each measurement.

	Catalog n°
PK-3 electrode polishing kit	A-011975
Content	
0.05 μm polishing alumina (20 mL)	A-001050
1 μm polishing diamond (10 mL)	A-002054
Glass plate (1 piece)	A-002249
Alumina polishing pad (10 pieces)	-
Diamond polishing pad (10 pieces)	-
Spare parts	
Alumina polishing pad (20 pieces)	A-001040
Diamond polishing pad (20 pieces)	A-001041
Emery paper UF800 (20 pieces)	A-012611
Coarse polishing pad (20 pieces)	A-001042
6 μm polishing diamond (10 mL)	A-002053



Working Electrodes		Isolation	OD/mm	Electrode Ø ±4%	Catalog n°	
Carbon	Long glassy carbon electrode	PEEK	6	3 mm	A-012744	
	Standard glassy carbon electrode	PEEK	10	5 mm	A-002417	
		PEEK	6	3 mm	A-002012	
		PEEK	6	1.6 mm	A-012297	
		PEEK	6	1 mm	A-002411	
	Small glassy carbon electrode	PEEK	3	1.6 mm	A-012298	
		PEEK	3	1 mm	A-002412	
	Micro carbon fiber electrode	glass	4	33 µm	A-002002	
		glass	4	7 µm	A-002007	
	Standard pyrolytic graphite electrode	Basal plane	PEEK	6	3 mm	A-002252
		Edge plane	PEEK	6	3 mm	A-002253
	Standard pfce carbon electrode	PEEK	6	3 mm	A-002408	
		PEEK	6	1 mm	A-002409	
		PEEK	3	1 mm	A-011854	
Platinum (99.95% purity)	Platinum gauze electrode 80 mesh	PEEK	-	25 x 35 mm	A-002250	
	Long platinum electrode	PEEK	6	3 mm	A-012745	
	Standard platinum electrode	PEEK	10	5 mm	A-002420	
		PEEK	6	3 mm	A-002422	
		PEEK	6	1.6 mm	A-002013	
	Small platinum electrode	PEEK	3	1.6 mm	A-002313	
	Micro platinum electrode	glass	4	100 µm	A-002009	
		glass	4	25 µm	A-002003	
		glass	4	15 µm	A-002015	
		glass	4	10 µm	A-002005	
	Gold	Gold gauze electrode 100 mesh	PEEK	-	25 x 35 mm	A-002251
		Long gold electrode	PEEK	6	3 mm	A-012746
		Standard gold electrode	PEEK	10	5 mm	A-002418
PEEK			6	3 mm	A-002421	
PEEK			6	1.6 mm	A-002014	
Small gold electrode		PEEK	3	1.6 mm	A-002314	
Micro gold electrode		glass	4	100 µm	A-002010	
		glass	4	25 µm	A-002004	
		glass	4	10 µm	A-002006	
Silver	Standard silver electrode	PEEK	10	5 mm	A-002416	
		PEEK	6	3 mm	A-002419	
		PEEK	6	1.6 mm	A-002011	
	Small silver electrode	PEEK	3	1.6 mm	A-002315	
Palladium	Standard palladium electrode	PEEK	6	1.6 mm	A-002019	
	Small palladium electrode	PEEK	3	1.6 mm	A-002319	
Nickel	Standard nickel electrode	PEEK	6	1.5 mm	A-002016	
	Micro nickel electrode	glass	4	100 µm	A-002273	
Copper	Standard copper electrode	PEEK	6	1.6 mm	A-002017	
		PEEK	6	3 mm	A-012584	
	Micro copper electrode	glass	4	25 µm	A-002271	
Iron (99.65% purity)	Standard iron electrode	PEEK	6	1.5 mm	A-002018	
		PEEK	6	3 mm	A-012585	
Carbon paste	Standard carbon paste electrode hole depth 4 mm	PEEK	6	3 mm	A-002210	
	Small carbon paste electrode hole depth 4 mm	PEEK	3	1.6 mm	A-002223	
	Cpo carbon paste oil base 1 g				A-001010	
Boron-doped diamond	Doping level between 500 and 1000 ppm. The electrode is 500 µm thick disk attached to a conductive rod in brass. It is polished with an Ra<10 nm.	PEEK	7	3 mm	M-BDD-3	

King-Size Reference Electrodes



King-Size Reference Electrodes		Connection type	Junction	Length	OD/mm	Catalog n°
Red rod reference electrode		Banana cable	Ceramic	103 mm	7.5	R-REF201
Ag/AgCl reference electrode Sat KCl for CV		Screw cap ⁽¹⁾	Ceramic	120 mm	8	R-XR300
Red rod reference electrode in 1 M KNO ₃		Screw cap ⁽¹⁾	Fiber rod	120 mm	8	R-XR440
Ag/AgCl reference electrode for EIS		Screw cap ⁽¹⁾	Ceramic	160 mm	8	R-XR820*
Options						
(1) Cable connection for screw cap electrode. We recommend you use this cable to connect these reference electrodes to BioLogic instruments.		Banana plug of 2 mm		100 mm		R-A94L111
		Banana plug of 4 mm		1 m		R-CL111
Salt bridges			Ceramic	70 mm	12	R-AL100
			Ceramic	138 mm	8	R-AL120
			Ceramic	250 mm	8	R-AL110
Salt bridge with reverse sleeve, non aqueous				138 mm	8	R-AL210
Conical rings for 8 mm OD 12 mm electrodes (4 pieces)						R-X31M012

*Ground joint Rin = 14.5/23

How to check your reference electrode

To check your reference electrode, perform an EIS measurement with a two-electrode connection. Choose Galvano Mode EIS (GEIS) to avoid polarization of your reference electrode.

The reference electrode is used as a working electrode and a platinum or gold electrode can be used as counter electrode.

The impedance of the reference electrode should be below 1 kΩ. If it is higher, the junction needs to be replaced.

Counter Electrodes



Dimensions in mm

Metallic Electrodes



Counter Electrodes	Size/mm	Wire Ø/mm	Surf. area/ cm ²	Purpose	Catalog n°
Platinum*	57	0.5	~ 0.7	SVC-2, VC-4, plate material evaluating cell	A-002222
	50	0.5	~ 0.7	SVC-3	A-002233
	230	0.5	~ 3.6	RRDE-3A, bulk electrolysis, SVC-3	A-002234
Gold	230	0.5	~ 3.6	RRDE-3A, bulk electrolysis, SVC-3	A-012638
Nickel	230	0.5	~ 3.6	RRDE-3A, bulk electrolysis, SVC-3	A-012639
Stainless steel	50	1.5	~ 2.35	Flow cell	A-012198
Platinum* gauze electrode, PEEK body 80 mesh	25x35	0.08	~ 22.9		A-002250
Platinum* gauze electrode, 54 mm wire 80 mesh	25x35	0.08	~ 22.9	Flat cell	A-702439
Platinum* electrode for bulk electrolysis cell 80 mesh	40x50	0.5	~ 47.4	Bulk electrolysis cell	A-013612
Gold gauze electrode, PEEK body 100 mesh	25x35	0.07	~ 29		A-002251

*99.95 % purity


Metallic Electrodes	Length/mm	OD/mm	Wire dimension/ mm	Catalog n°
Platinum ⁽¹⁾	80	8	Ø 1	R-XM110
	120	8	Plate 5x5	R-XM120
	120	8	Plate 8x8	R-XM140
	120	12	Disk of 10	R-XM150

Options

(1) Cable connection for screw cap electrode. We recommend you use this cable to connect these refer- ences electrodes to BioLogic instruments.	100 mm, banana plug of 2 mm	R-A94L111
	1 m, banana plug of 4 mm	R-CL111
Conical rings for 8 mm OD 12 mm electrodes (4 pieces)		R-X31M012

Small-Size Reference Electrodes for aqueous media


Reference electrodes are divided into two groups according to the media in which the electrode is immersed (aqueous or organic media).




Small-Size Reference Electrode for Aqueous Media (Ag/AgCl, Hg)

	Junction	Electrolyte	Purpose	Catalog n°
RE-1B Ag/AgCl reference electrode ⁽¹⁾	IPPG*	3 M NaCl	SVC-2, SVC-3, VC-4, bulk electrolysis, RRDE-3A, flat cell	A-012167
RE-1S Ag/AgCl reference electrode ⁽¹⁾	IPPG*	3 M NaCl	SECM	A-012168
RE-1BP reference electrode (Ag/AgCl)	Ceramic	3 M NaCl	SVC-2, SVC-3, VC-4, bulk electrolysis, RRDE-3A, EQCM, flat cell	A-013613
RE-1CP Ag/AgCl reference electrode	Ceramic	Saturated KCl	SVC-2, SVC-3, VC-4, bulk electrolysis, RRDE-3A, flat cell	A-013429
RE-3VT Ag/AgCl reference electrode screw type	Ceramic	3 M NaCl	For flow spectroelectrochemical cell SEC-3F	A-013488
Ag/AgCl ink, 2 mL	Surface resistance: 0.2 Ω /25.4 μ m ² Viscosity: 50,000 \pm 10,000 CP @ 21.1 $^{\circ}$ C Flash point: 82 $^{\circ}$ C		For micro CV cell, IDA measurement	A-011464
RE-2BP Hg/Hg ₂ Cl ₂ reference electrode	Ceramic	Saturated KCl	SVC-2, SVC-3, VC-4, bulk electrolysis, RRDE-3A, flat cell	A-013430
RE-2CP Hg/Hg ₂ SO ₄ reference electrode, free from chloride	Ceramic	Saturated K ₂ SO ₄	SVC-2, SVC-3, VC-4, bulk electrolysis, RRDE-3A, flat cell	A-013431
RE-61AP Hg/HgO reference electrode main body in polyacetal resin	Ceramic	1 M NaOH	For alkaline media	A-013395
RHE Reversible Hydrogen Electrode kit	IPPG*	HCl or H ₂ SO ₄	For acidic media (pH < 2)	A-013373
Spare parts				
Set of 10 porous 2.8 mm glass frits (CoralPor™) with PTFE heat shrink (200 mm)				092-VYC3
Options				
RE-PV preservative vial for reference electrode, 10 mL				A-012108
Bridge tube \varnothing 9.0 mm (2 pieces)				A-012177
Bridge tube \varnothing 9.0 mm (22 pieces)				A-012307
Double junction chamber kit for RHE				A-013375

* Ion Permeable Porous Glass



Small-Size Reference Electrodes for non aqueous media



Small-Size Reference Electrodes for Non Aqueous Media (Ag/Ag⁺)

	Junction	Electrolyte	Purpose	Catalog n°
RE-7 non aqueous reference electrode (Ag/Ag ⁺)	IPPG*	Ag/Ag ⁺ /ACN**/TBAP***	CV	A-012171
RE-7S non aqueous reference electrode (Ag/Ag ⁺)	IPPG*	Ag/Ag ⁺ /ACN**/TBAP***	SECM	A-012172
RE-7VT non aqueous reference electrode (Ag/Ag ⁺) with poly-methyl pentene body	Ceramic	Ag/Ag ⁺ /ACN**/TBAP***	For flow cell (LC, EQCM, SEC-2F).	A-013489

Spare parts

Electrolyte solution (10 mL)	A-012549
PTFE cap with Ag wire (for RE-7)	A-012057
PTFE cap with Ag wire (for RE-7S)	A-012058
Sample holder 6 mm diameter (for RE-7) (2 pieces)	A-012176
Set of 10 porous 2.8 mm glass frits (CoralPor) with PTFE heat shrink (200 mm)	092-VYC3

Options

RE-PV preservative vial for reference electrode, 10 mL	A-012108
Bridge tube Ø 9.0 mm (2 pieces)	A-012177
Bridge tube Ø 9.0 mm (22 pieces)	A-012307

* Ion Permeable Porous Glass
 ** Acetonitrile
 *** Tetra Butyl Ammonium Perchlorate



Maintenance of Reference Electrodes.

Maintenance of reference electrodes

Store your reference electrode immersed in the electrolyte

When not in use, we recommend that you keep reference electrodes in sealed, air-tight vials in order to prolong their life. The storage solution should be identical to the filling solution of the reference electrode. Prefer a cold and dark place.

* Ion Permeable Porous Glass

Prevent contamination

To prevent contamination of the reference electrode, a sample holder can be used (see p. 14).

Replace the junction when needed

If you are using IPPG junctions, yellowish discoloration indicates contamination. This is caused by the absorption of organic compounds into the pores. The average pore diameter of IPPG is about 40 – 200 Å.

If you are using CoralPor™ junctions, you might want to use the replacement kit 092-VYC3 (see page 20).

The average pore diameter of CoralPor™ is about 4 – 10 nm.

Storage of reference electrodes

To avoid electrolyte leakage or concentration due to evaporation during storage or transport, the electrolyte is separated from the body. This enables it to be regularly filled up.

Potentials of common reference electrodes


E/V vs. NHE at 25 °C

0.930	Hg/HgO/ NaOH (0.1 M)
0.650	Hg/Hg ₂ SO ₄ / K ₂ SO ₄ (sat)
0.624	Fc/Fc ⁺ TBAP (0.1M) ACN
0.542	Ag/Ag ⁺ TBAP (0.1M) ACN
0.241	Hg/Hg ₂ Cl ₂ KCl (sat)*
0.236	Hg/Hg ₂ Cl ₂ NaCl (sat)*
0.205	Ag/AgCl/ KCl (3.5 M)
0.197	Ag/AgCl/ KCl (sat)
0.194	Ag/AgCl/ NaCl (sat)
0.000	NHE Normal Hydrogen Electrode

* Hg/Hg₂Cl₂ : Calomel

Other Electrodes.


Screen printed electrodes



Description	Catalog n°
3 electrode SPEs: WE and RE: graphite ; CE: silver printed on alumina (48 pcs)	L-CCAG-AL203
3 electrode SPEs: WE and RE: graphite ; CE: silver printed on PET (48 pcs)	L-CCAG-PET
Connector for L-CCAG-AL203 and L-CCAG-PET with 2 mm receptacles	L-CONNECTOR

InterDigitated Array (IDA) electrodes


The passivation membrane is a Novolac resin + naphtoquinone-diazo compounds.



IDA Electrode	Width/ μ m	Interval/ μ m	Length/mm	N° of feet	Film thickness/nm	Catalog n°
Gold	2	2	2	65 pairs	90	A-012257
Platinum	2	2	2	65 pairs	90	A-012258
Gold without passivation membrane	2	2	2	65 pairs	90	A-012261
Platinum without passivation membrane	2	2	2	65 pairs	90	A-012264
Gold	3	3	2	65 pairs	90	A-012129
Platinum	3	3	2	65 pairs	90	A-012130
Gold without passivation membrane	3	3	2	65 pairs	90	A-012260
Platinum without passivation membrane	3	3	2	65 pairs	90	A-012263
Gold	10	5	2	65 pairs	90	A-012125
Platinum	10	5	2	65 pairs	90	A-012126
ITO	10	5	2	65 pairs	100 \pm 20	A-012128
Carbon	10	5	2	65 pairs	1200 \pm 100	A-012127
Gold without passivation membrane	10	5	2	65 pairs	90	A-012259
Platinum without passivation membrane	10	5	2	65 pairs	90	A-012262
ITO without passivation membrane	10	5	2	65 pairs	100 \pm 200	A-012265
Carbon without passivation membrane	10	5	2	65 pairs	1200 \pm 100	A-012266

Options	
Ag/AgCl ink for reference electrode (2.0 mL)	A-011464
Cable kit for IDA electrode	A-011066

Content	
Electrode fixer	-
Mini vice	-
Connecting cable	-



12.5 x 22 x 0.5 mm

Ring-disk type electrodes

Ring-disk type electrodes	Ring OD/mm	Ring ID/mm	Disk/mm	Catalog n°
Gold ring-disk electrode (3 pieces)	6	4	3	A-002081
Platinum ring-disk electrode (3 pieces)	6	4	3	A-002082
Carbon ring-disk electrode (3 pieces)	6	4	3	A-002083

Rotating Electrodes.

BluRev Rotating Disk Electrode (RDE)


The **BluRev RDE** is a versatile and rugged rotating disk electrode ideal for use with any BioLogic potentiostat/galvanostat. It is available with a wide choice of quick-fit exchangeable tips.

The **BluRev RDE** has been designed to be used with EL-ELECTRO cells series. Please note that if you already have an EL-ELECTRO cell an additional PTFE cap is needed: **094-A-CAP**. Otherwise, **EL-BLUREV** is the cell compatible with **BluRev RDE**.

Furthermore, a special enclosure has been designed for an easy operation and set-up of the **BluRev RDE**: the **BluRev Enclosure**. This makes it easy to hold the RDE body as well as the BioLogic potentiostat cables.

The RC-10k unit controls the RDE rotation speed. It offers an accuracy of ± 1 rpm over the whole rotational range for precise and fully reproducible experimental conditions. The measured rotation speed is clearly displayed on a five-digit LCD screen. The speed can be set manually or remotely by using the analog output of a BioLogic instrument. Any other device with an analog output can be used.

Using EC-Lab[®], you can easily program and perform a Levich experiment all in one go, and also use the new EIS element W_{inf} that will directly give you the diffusion coefficient.



BluRev RDE		Catalog n°
BluRev RDE kit (tips & cell are not included)		094-RC/RDE
Content		
RC-10k Rotation controller		094-RC
Rotating disk electrode (motor, shaft, electrode body with contact)		094-RDE
DB9 to BNC connector for external control of RC-10k		092-22/1
1 m BNC/BNC cable		COR28100
1 replacement contact		094-RDE-BRUSH
1 transport case		-
Tips		
2 mm diameter 999% Platinum disk electrode		094-Pt/2
2 mm diameter 999% Gold disk electrode		094-Au/2
3 mm diameter Glassy Carbon disk electrode		094-GC/3
3 mm diameter 999% Silver disk electrode		094-Ag/3
3 mm diameter 999% Copper disk electrode		094-Cu/3
3 mm diameter 999% Nickel disk electrode		094-Ni/3
3 mm diameter 999% Aluminium disk electrode		094-Al/3
3 mm diameter 316L Stainless steel disk electrode		094-316L/3
5 mm diameter Glassy Carbon disk electrode		094-GC/5
Cells		
EL-ELECTRO-1A cell compatible with BluRev RDE		EL-BLUREV
PTFE 5 holes cap compatible with BluRev RDE (needed if you already have an EL-ELECTRO cell)		094-A-CAP

Specifications

Rotational range/rpm	100 - 10,000
Accuracy/rpm	1 typical over the whole rotational range
Precision/rpm	3 typical over the whole rotational range
Eccentricity/mm	<0.1
Setting resolution/rpm	10 (Manual control) or 25 (Remote control)
Materials of RDE	PEEK, Al alloy
Gas inlet for brush contact purging/mm	Ø 2
Operating temperature/°C	10 - 40
Power	24 Vdc, 1 A max
Max consumption/W	24
Dimensions/mm	RC-10k: 95x227x178 (HxWxD), RDE: 233.6 (length with tip)
Weight/kg	RC-10k: 1.00 RDE: 0.36 (without tip)



BluRev Enclosure		Catalog n°
Protective housing and stand kit for the BluRev		094-ENCL
Content		
1 stainless steel plate with M6 threaded holes to fix the support poles		-
2 support poles to hold the BluRev RDE and the cell		-
(The clamp and the clamping nut are not included)		
1 clamping flange for the BluRev		-
3 half-clamps for the potentiostat cables (all BioLogic cables are supported)		-
1 plexiglas protective housing		-
Specifications		
Dimensions with protective housing (HxWxD)/mm		287x318x308
Weight (with protective housing)/kg		5.3

RDE


W_{inf} : direct access to the diffusion coefficient.

Levich and Koutecký-Levich methods¹ are powerful analysis tools used to obtain kinetic electrochemical parameters such as the diffusion coefficient of a redox species in a given medium and the reaction constant. These analyses require potentiodynamic curves at various rotation speeds.


However, fitting impedance measurements made on a redox reaction occurring at a rotating electrode, at only one rotation speed, also enables the direct measurement of the diffusion coefficient.

A PEIS measurement was performed on a 2 mm Pt electrode, using an equimolar solution of $K_3Fe(CN)_6$ and $K_4Fe(CN)_6$ with concentrations of 5 mM in 0.1 M KCl, a BluRev RDE rotating at 2000 rpm.


The considered electrochemical reaction is the following:



The obtained impedance graph is shown below:



Using ZFit and the equivalent circuit shown below that contains W_{inf} as a diffusion impedance, we directly obtain the diffusion coefficient of the species of interest, in our case $7.8 \times 10^{-6} \text{ cm}^2/\text{s}$, which is in agreement with the data found in the literature^{2,3}. For more detailed information please see the EC-Lab application note #66⁴.



RRDE

For RRDE measurements, a bipotentiostat is needed. A bipotentiostat controls the two working electrodes i.e. one channel to control the disk electrode and the other to control the ring. An SP-300 equipped with two channels would be an appropriate instrument for RRDE applications.

Because of the presence of two working electrodes in the same setup, a specific connection mode (to avoid any ground loop trouble) is needed. This can be achieved in two different ways:

- isolating the two channels (at least one potentiostat in floating mode),
- grounding the counter electrode. The "CE-to-ground" mode of the BioLogic instruments provides this unique capability.

Of these two options, the latter is preferred because there is less leakage of current. All of the BioLogic multichannel potentiostats (except cyclers) offer such types of electrode connection.

At the disk electrode, the electro-active species are oxidized or reduced according to the applied potential. This new species is detected by reduction or oxidation at the ring, respectively.

In a typical experiment, a CV is performed on the disk electrode and a constant voltage is applied on the ring electrode. This is the CV-CA technique that is available in EC-Lab® in the "bipotentiostat" techniques folder.

(1) Application Note #56 "Electrochemical reaction kinetics measurement: the Levich and Koutecký-Levich analysis tools"

(2) A. J. Bard, W. Faulkner, in: Electrochemical Methods, Fundamentals and Applications, 2nd Ed., Wiley, New York (2001) 381.

(3) D. R. Lide, H. V. Kehiaian, in: CRC handbook of thermophysical and thermochemical data, CRC Press, Inc., Boca Raton, (1994)

(4) Application Note 66 "EIS measurements on a Rotating Disk Electrode (RDE) Part I: Determination of a diffusion coefficient using the new element W_{inf} "



Rotating Electrodes.

Rotating Ring-Disk Electrode (RRDE)

RRDE-3A is a system that can be used to perform hydrodynamic (RDE or RRDE) measurement. Its design (short stainless steel shaft) allows users to accurately control electrode rotation and modulation.

RRDE-3A is electronically controlled by a closed loop circuit driving a DC servo-motor. Electrodes are small and rapidly interchangeable.

The unit also provides an adjustable valve system for inert gas purging inside the cell vial. Additionally, access to the glass cell facilitates rinsing, cleaning, and replacing the electrodes. It is also easy to remove and replace the cell vial.

RRDE-3A can be operated as a stand-alone unit or directly controlled by a BioLogic electrochemical workstation.

Features

- Remote and manual-controlled rotation
- Small sample volume (up to 10 – 15 mL), gas line connector is available, remote and manual controlled purge lines
- Open architecture for an easy access to the cell

Specifications

Rotational range/rpm	100 – 8,000
Setting resolution/rpm	1
Accuracy	< 0.1%
Rotation control type	PWM (Pulse Width Modulation)
Bandwidth	60 Hz at 3,500 rpm Base and 1,000 rpm Peak-to-peak modulation
Inlet gas pressure	5 PSI
Temperature/°C	10–50
Power	100–240 VAC, 50/60 Hz
Dimensions (HxWxD)/mm	400x185x(base: 230, body: 120)
Weight/kg	6



Rotating Ring-Disk Electrode

Catalog n°

RRDE-3A apparatus V 2.0

A-012180


Working electrodes (page 27), reference electrodes (page 20) and counter electrodes (page 19) must be purchased separately.

Content

RRDE-3A glass cell vial 100 mL	A-013580
Spin coating adapter	A-012064
Nipple for purge valve	A-012065
Tygon tubing, OD 6.4 mm x ID 3.2 mm (1.3 m)	A-010058
PTFE cap for RRDE-3A	A-012631
Instruction manual	-
Power cable	-



Options

Sample vials for alkaline solution (100mL) (10 pcs)	A-013580
Sample vials for alkaline solution (200mL) (8 pcs)	A-013581
Water-jacketed glass cell (100 mL, OD 70 mm, ID 46.4 mm, H 80 mm)	A-012652
DB9 cable to control RRDE-3A	092-22/11
Bipotentiostat cable for two standard channels	092-22/12
Corrosion resistant bearing assembly	A-013605



A-011172 A-011162 A-011163 A-011164 A-011170 A-011171 A-011169

Working Electrode

Standard tips		Tips length/mm	Tips OD/mm	Ring ID/mm	Ring OD/mm	Disk Ø/mm	Catalog n°
 RRDE	Platinum ring/platinum disk	25	12	5.0	7.0	4.0	A-011172
	Platinum ring/glassy carbon disk	25	12	5.0	7.0	4.0	A-011162
	Gold ring/glassy carbon disk	25	12	5.0	7.0	4.0	A-011163
	Platinum ring/gold disk	25	12	5.0	7.0	4.0	A-011164
	Gold ring/platinum disk	25	12	5.0	7.0	4.0	A-012617
	Glassy carbon ring/glassy carbon disk	25	12	5.0	7.0	4.0	A-012618
 RDE	Glassy carbon disk	25	12	-	-	3.0	A-011169
	Glassy carbon disk	25	12	-	-	5.0	A-013482
	Glassy carbon disk with PTFE spacer	25	12	-	-	3.0	A-013490
	Glassy carbon disk with PTFE spacer	25	12	-	-	5.0	A-013491
	Platinum disk	25	12	-	-	3.0	A-011170
	Gold disk	25	12	-	-	3.0	A-011171
	Aluminum disk	25	12	-	-	3.0	A-011966
	Silver disk	25	12	-	-	3.0	A-011967
	Copper disk	25	12	-	-	3.0	A-011968
	Nickel disk	25	12	-	-	3.0	A-011969
	Tantalum disk	25	12	-	-	3.0	A-011970
	Titanium disk	25	12	-	-	3.0	A-011971
	Tungsten disk	25	12	-	-	3.0	A-011972
	Carbon paste disk, hole depth 4 mm	25	12	-	-	3.0	A-011973

Option

PK-3 polishing kit

A-011975

Disk Replaceable Electrode tips		Catalog n°
RRDE	Platinum ring/GC disk replaceable electrode kit	A-013336
	Content Platinum ring assembly	A-013337
	Glassy carbon disk	A-013338
	PTFE spacer (3 pieces)	A-013339
RDE	Glassy carbon disk replaceable electrode kit	A-013362
	Content Glassy carbon disk	A-013338
	PTFE spacer (3 pieces)	A-013339
	Disk assembly	A-013361
	Gold disk replaceable electrode kit	A-013364
	Content Gold disk	A-013366
	PTFE spacer (3 pieces)	A-013339
	Disk assembly	A-013361
	Platinum disk replaceable electrode kit	A-013365
	Content Platinum disk	A-013367
	PTFE spacer (3 pieces)	A-013339
	Disk assembly	A-013361

Tool kit

Disk electrode polishing/exchange kit **A-013340**




RRDE platinum ring/GC disk replaceable electrode kit



Disk electrode polishing/exchange kit

Spectroelectrochemistry.



Spectrometer System

Spectroelectrochemistry (SEC) can be useful to elucidate electrochemical reaction mechanisms.

The spectroelectrochemical kit is made up of three parts (spectrometer, light source and cuvette holder).

The spectrometer is equipped with a trigger to synchronize electrochemical and spectroscopic measurements.

SEC2020 Spectrometer system

Catalog n°

SEC2020 spectrometer kit

A-013609

Content

SEC2021 Spectrometer (x 1)	-
SEC2022 Deuterium halogen light source (x 1)	-
SEC2023 Cuvette holder (x 1)	-
SEC2024 Platform (x 1)	-
AC adaptor (x 1)	-
Power cable (x 1)	-
USB cable (x 1)	-
Collimator (x 2)	-
Fiber collimator (x 1)	-
Platform screw (x 7)	-
External device connection trigger cable (x 1)	-
Light source control trigger cable (x 1)	-
Plastic cuvette (x 1)	-
SMA905 adaptor for light source (x 3)	-
SMA905 adaptor for light shielding (x 2)	-
Hexagon wrench 0.89 mm (x 1)	-
Hexagon wrench 1.50 mm (x 1)	-
Software (USB memory) (x 1)	-
Waterproof box (x 1)	-

Quick manual, wavelength calibration data sheet, linearity test data sheet and warranty certificate are also included.


Option

Connecting cable to synchronize the SEC2020 with BioLogic instrument.

092-22/11

Spectrometer specifications

Detector	2048 element linear silicon CCD array
Full description	SEC2021-025-DUVN
Detector range/nm	200 - 1025
Grating	Blaze wavelength (300 nm)
Slit/ μm	25
Wavelength resolution/nm	1.3
Fiber connector	SMA905 Core diameter: 600 μm NA=0.22
Interface	USB2.0
Operating system	WindowsTM 7 / 8.1/ 10 (32bit/ 64bit)
Dimensions (HxWxD)/mm	32x86x110




Light source specifications

Light type	Deuterium halogen light source
Wavelength range/nm	200 - 1700
Stability	<0.1%
Drift/h	0.25%
Bulb life/h	>1000 (D2 lamp) >2000 (halogen lamp)
Fiber connector	SMA905
Size (HxWxD)/mm	46x100x165

The SEC2020 spectrometer system uses the Czerny-Turner optical mount. This system is an M-shaped structure symmetrical to the grating (4) and is an optical system with extremely small aberration.

Light source structure



1. SMA905 Connector
2. Slit
3. Collimating mirror
4. Grating
5. Focus mirror
6. 2048 element CCD array

Utilization modes

Transmittance

- Absorbance/transmittance**
- Concentration of chemicals (solution)
 - Polymer extrusion processes
 - DNA quantification

- Reflectance**
- Freshness testing
 - Film thickness/composition (quality control)
 - Activation energy of photocatalytic species
 - Textile quality control


Fluorescence

- Fluorescence**
- Marine organisms
 - Biology (DNA, protein, cell proliferation assay, histamine-analysis, alga monitoring)
 - Environmental fields (waste water analysis, ground water trace studies, hydrocarbon detection, dissolved oxygen)
 - Plant efficiency (plant physiology, plant breeding, horticulture, agronomy, agrochemicals, pollution studies, forestry, ecology)
 - Tissue diagnosis


- Scattering**
- Oil concentrations of oil/water system
 - Raman spectroscopy
 - Physical transition phenomena (e.g. melting point, glass transition crystallize temperature)

Irradiance



- Emission**
- Astronomy (e.g. spectra of Hale-Bopp, plasma monitoring)
 - In situ metal monitoring
 - Luminescence (PL,EL), LED & laser wavelength



Static Cell




Quartz glass spectroelectrochemical cell kit	Platinum	Catalog n°	Gold	Catalog n°
0.5 mm		A-012813		A-012814
Content				
Platinum counter electrode		A-012609		A-012609
Thin layer quartz glass cell		A-012815		A-012815
PTFE cap		A-011501		A-011501
Purging tube (ETFE, 100 mm)		-		-
Gauze working electrode	80 mesh, height 5 mm	A-012606	100 mesh, height 5 mm	A-012607
		A-013510		A-013511
1 mm				
Content				
Platinum counter electrode		A-012906		A-012906
Thin layer quartz glass cell		A-012907		A-012907
PTFE cap		A-011501		A-011501
Purging tube (ETFE, 100 mm)		-		-
Gauze working electrode	80 mesh, height 6 mm	A-011498	100 mesh, height 6 mm	A-012017
Options				
RE-1BP Ag/AgCl reference electrode				A-013613
RE-7 non aqueous reference electrode				A-012171
Purging tube (ETFE), 1 m				A-010537



Flow Cell

Flow Cell		Quantity	Catalog n°
SEC-2F spectroelectrochemical flow cell			A-012660
Content			
SEC-2F flow cell	Base	1	-
	Cover	1	-
	Block a	1	-
	Block b	1	-
SEC-2F S500 silicone gasket		2	A-012661
Stainless tube OD 1.59 mm (length of 50 mm)		1	A-012198
Needle adaptor		1	-
Dynaseal PEEK fingertight		2	-
Silicon tube (300 mm)		1	-
PTFE tube (1 m)		1	-
Options			
Reference electrode screw type	Ø 10x55 mm	RE-3V aqueous	A-012169
		RE-3VP aqueous (PEEK)	A-012170
		RE-7VT non aqueous	A-013489
Working grid electrode for flow cell	8x27x1 mm	Platinum (1 piece)	A-012655
		Gold (1 piece)	A-012656
		Carbon grid electrode (1 piece)	A-012657
		ITO electrode (4 pieces)	A-012658
	8x27x0.5 mm	ITO electrode (12 pieces)	A-011465
	10x20x0.5 mm	ITO electrode (10 pieces)	A-010887
Gasket	Silicone	10x10x0.5 mm	ITO electrode (30 pieces)
		Ø 4 inchx0.5 mm	ITO electrode (1 piece)
			A-011233
			A-011827
Fibre and lens	PTFE	S500, 500 µm thick (4 pieces)	A-012661
		T500, 500 µm thick (4 pieces)	A-012664
		T250, 250 µm thick (4 pieces)	A-012665
		T100, 100 µm thick (4 pieces)	A-012666
		400 µm optical fibre SR, 250 mm	A-012667
		400 µm optical fibre SR, 2 m	A-011522
		UV/VIS collimating lens, 200-2,000 nm	A-012234




Gold grid electrode **A-012656**

Solution volumes

		Volume
Gasket of	100 µm	4.6 µl
	250 µm	11.5 µl
	500 µm	230 µl


Thin Layer Cells

The UFS set has been designed to perform thin-layer (optical path ~ 0.2 mm) spectroelectrochemistry experiments.



UF-Full

Spacer



Pseudo-reference electrode

Working electrode

Auxiliary electrode

Embedded in LDPE

Thin layer spectroelectrochemical cell kit

Full kit

Content

1 PTFE windows aligner + 1 spacer (platinum auxiliary electrode/Ag pseudo-reference electrode)
Stainless steel body cell
XY holder (with magnets) to universally fit the cell in a spectrophotometer
Plug
PTFE mask

Platinum

Catalog n°

UF-Full-PT

Gold

Catalog n°

UF-Full-AU

Platinum working electrode

UF-SEC-PT

Gold working electrode

UF-SEC-AU

-

UF-XY

-

UF-P

UF-M

-

UF-XY

-

UF-P

UF-M

Options

PTFE adapter for optical fibres with UV/VIS/NIR collimating lenses (200-2500 nm), adjustable focus, SMA-905 connection

UF-OFA

Base stand for optical fibre measurements

UF-BS

PTFE cuvette holder to be used with UF-OFA and UF-BS

UF-C

1 spacer (platinum working electrode/platinum auxiliary electrode/Ag pseudo-reference electrode)

UF-SPP

1 spacer (gold working electrode/platinum auxiliary electrode/Ag pseudo-reference electrode)

UF-SPA

Optical fibre: SMA-SMA connector length 1 m

092-101



Full kit + PTFE adapter for optical fibres + optical fibres



Finite Diffusion Conditions

In such finite diffusion conditions, a thin solution layer (≤ 0.2 mm) adjacent to the electrode is confined by the cell walls, so that the cell thickness is smaller than the diffusion layer and the mass transfer can be ignored. The most significant virtue of thin-layer cells is the absence of the effect of the diffusion process and the rapidity with which the electro-active species can be completely electrolyzed. The dropping to near zero of the current flow following the peak in the current potential plot is a characteristic behavior of thin-layer cells, indicating exhaustive electrolysis of the cell reactant and minimal diffusion effects in thin-layer electrochemical cells^{(1) (2)}.

Small potential sweep rates (2-10 mV/s) are necessary both to ensure homogeneity of the reactant/product concentrations in the cell and to control resistive effects.

Cyclic voltammetry should be used in your SEC cell to better identify the redox process of interest: in fact, the UF spacer has a silver pseudo-reference electrode, which is sensitive to the solution medium (but is expected to remain constant in time in each given experimental

condition). Potential drifts can be observed in the presence of irreversible redox reactions, which may alter the solution/analyte composition.

To keep the ohmic drop as low as possible, the current should also be kept low which means low scan rates and low sample concentrations may be preferred. However, in some cases, due to the narrowness of the optical path, the use of a relatively high concentration of the sample may be required to study the changes of bands with a low molar extinction coefficient.



A milli-molar concentration appears to be the most optimized concentration, but in some cases, the concentration can be adapted.

Due to the high Infra-Red absorbance of more common solvents and electrolytes, a very carefully measured background should be obtained before each Infra-Red spectroelectrochemical experiment.

For the same reason, it is also important to avoid changing the tightness of the cell screws during the experiment itself to avoid changes of the optical path.

(1): M. Krejčík, M. Daněk and F. Hartl, J. Electroanal. Chem., 1991, 317, 179.

(2): P. Leoni, F. Marchetti, C. Bonaccorsi, F. Fabrizi de Biani, L. Marchetti, P. Zanello, Chem. Eur. J., 2008, 14, 847



Quartz Crystal Microbalance.

BluQCM QSD

Acoustic sensing principle

The acoustic sensing principle is based on the precise detection of changes on the properties of an acoustic (mechanical) wave traveling through the bulk (QCM-AWS, HFF-AWS) or the surface (LOVE-AWS) of the acoustic wave sensor.

QCM-AWS sensors consist of a thin piece of quartz confined between a pair of metal-based electrodes. An alternating current applied to the quartz crystal induces mechanical oscillations on the quartz due to the piezoelectric effect. A wave is generated and propagated through the sensor and the films attached to it.

The resonance frequency of this wave depends on the oscillating mass of the sensor and its adhering layers. When a thin film is attached to the sensor, the properties of the wave change as well, modifying the resonance frequency and amplitude. If the film is thin and rigid, the decrease in frequency is proportional to the mass of the film.

In this way, the QCM works in the gravimetric regime and the mass of the film can be calculated using the well-known Sauerbrey equation. If the film is not rigid, the measurement of the damped resonance enables the measurement of the viscoelastic changes of the film, allowing the characterization of the QCM working in the non-gravimetric regime.

A mass increase results in a frequency decrease. Sauerbrey was the first to provide a description and experimental verification of the mass/frequency relationship between foreign layers firmly attached to the quartz crystal resonator⁽¹⁾⁽²⁾.

The Sauerbrey equation is defined as:

$$\Delta f = \frac{2f_0^2}{A\sqrt{\rho_q\mu_q}} \Delta m$$

f_0 : resonant frequency (Hz)

Δf : frequency change (Hz)

Δm : mass change (g)

A: piezoelectrically active crystal area (area between electrodes, cm²)


ρ_q : density of quartz ($\rho_q = 2.648 \text{ g/cm}^3$)

μ_q : shear modulus of quartz for AT-cut crystal ($\mu_q = 2.947 \times 10^{11} \text{ g/(cm.s}^2\text{)}$)



Special features

- High accuracy temperature control
- Dissipation measurement
- Measurement in air
- Measurement at overtones



(1): G. Sauerbrey, Phys. Verh., 1957, 8, 113-114.

(2): G. Sauerbrey, Z. Phys., 1959, 155, 206-222.

Sensor	14 mm wrapped	1" wrapped	HFF
Cell			
In-batch eQCM 	AW-BEQ01Q	AW-BEQ02Q	AW-BEQ01HQ
Flow eQCM 	AW-FEQ01Q	-	AW-FEQ01HQ
Flow QCM 	AW-FQ01Q	-	AW-FQ01HQ
In-batch QCM 	AW-BQ01Q	AW-BQ02Q	AW-BQ01HQ
Hermetic Li research in batch eQCM 	AW-BEQLIQ	-	-
In batch probe eQCM 	AW-PEQ11Q	-	-

All these cells feature the patented Quick-Lock design that facilitates operation and ensures an excellent reproducibility.

Most of these cells are available with SMA connectors to be used with standard impedance and network analyzers. Please contact your local reseller for further details.

Sensor type	14 mm	1"	HFF
Resonant frequency/MHz	5, 10	5, 9	50, 100, 150
Available finish	Rough (Au/Ti only), polished		N/A
Available material	Au/Cr Au/Ti C, Pt (5 MHz only)		Au/Cr

To obtain the catalog n° and if a different material is needed, please ask your local distributor.

Reference and counter electrodes have to be purchased separately (except for eQCM flow cells where the Pt plate counter electrode is integrated in the lid of the cell). Compatible counter and reference electrodes are shown in the following table:

	Reference electrode		Counter electrode
	Aqueous	Non-aqueous	
In-batch eQCM cells	RE-1B A-012167	RE-7 A-012171	Pt wire 23 cm coiled A-002234
Flow eQCM cells	RE-1S A-012168	RE-7S A-012172	Pt disk integrated in the cell lid

QSD-300

General function	Tracking* and high resolution** at single and multiple overtones
Operation modes	QCM, HFF-QCM, LOVE-SAW
Sensors technologies	See p.3
Cells available	3 – 45 (depending on sensor and cell type)
Liquid volume above sensor/μL	Optional (QSD-TCU)
Temperature control	Yes
Measurement in air	Up to 7 (up to 13th)
Simultaneous overtones measurements	90x220x260
Dimensions (HxWxD)/mm	3
Weight/ kg	
Sensor	
Frequency range/MHz	4 – 160
Best frequency resolution/Hz	0.1
Best frequency accuracy/Hz	± 0.5
Max. acquisition rate***/points/s	250
Best mass sensitivity in liquid****/pg/cm²	50
Best dissipation sensitivity	10 ⁻⁷
Catalog n°	
BluQCM QSD-300	AW-QSD-300

* Tracking mode provides the full impedance spectrum of the sensor around resonance frequency

** Patented fast and high-resolution single frequency point measurement

*** High-resolution mode at single frequency

**** For HFF-QCM

QSD-TCU

General function	
Temperature control range/°C	15 – 45
Temperature stability/°C	±0.05
Dimensions (H x W x D)/mm	60x220x260
Weight/kg	4.5
Catalog n°	
BluQCM QSD-TCU	AW-QSD-TCU

QSD-FCU

General function	
Syringe volume/μL	250 (default)*
Flow rate range for a 250 μL syringe/μL/min**	12.5-14500 (Standard) 0.625 – 1062.5 (Smooth)
Dimensions (H x W x D)/mm	195x70x250
Weight/kg	0.75
Catalog n°	
Standard flow control unit	AW-QSD-FCU
Smooth flow control unit	AW-QSD-FCUS

*Other syringe volumes are available upon request, from 12.5 μL to 5000 μL

** Flow rates depend on the syringe volume. For the standard flow unit, the flow rate change is 0.6250 – 290000 μL/min. For the smooth flow unit, it is 0.0313– 21250 μL/min. For more information, please contact your local reseller.

Quartz Crystal Microbalance.

QCA922A

This instrument monitors both the resonant frequency and the resonant resistance which are also reflected on the two analog outputs.


It is possible to connect the QCM device to a BioLogic potentiostat using a DB9-8BNC cable (catalog n°: 092-22/1).



092-QCA-FC



SE-CL3



eQCM			Catalog n°			
Quartz crystal microbalance 27 MHz kit			SE-QCA922A			
Content						
EQCM 27 MHz main unit and cable			SE-QCA922A-00			
Male BNC/BNC cable (2 pieces, length 1 m)			COR28100			
Options						
Connector from QCA to potentiostat			092-22/1			
Low flow peristaltic pump			EL-AV-008			
Cells						
Holder is needed to get a full QCM or EQCM set-up						
Static and flow QCA cell			092-QCA-FC			
Dip cell			SE-CL3			
Well cell (PTFE)			SE-CL4			
Well cell (PEEK)			SE-CL4PK			
Transparent well cell			SE-CL5			
Flow cell (90 µl) (PTFE)			SE-CL6			
Flow cell (90 µl) transparent (PTFE)			SE-CL7			
Flow cell (90 µl) (PEEK)			SE-CL6PK			
Resonators						
5 MHz						
Standard finish	Resonator and lead wire	Gold electrode (25 pieces)	SE-5AU			
		Aluminum electrode (25 pieces)	SE-9AL			
		Gold electrode (25 pieces)	SE-9AU			
		Copper electrode (25 pieces)	SE-9CU			
		Molybdenum electrode (25 pieces)	SE-9MO			
		Nickel electrode (25 pieces)	SE-9NI			
		Platinum electrode (25 pieces)	SE-9PT			
		Stainless steel (SS304) electrode (25 pieces)	SE-9SS			
		Separated lead wire	Gold electrode (25 pieces)	SE-9AU-S		
			Platinum electrode (25 pieces)	SE-9PT-S		
	Mirror finish	Resonator and lead wire	Aluminum electrode (25 pieces)	SE-9AL-M		
			Gold electrode (25 pieces)	SE-9AU-M		
			Gold electrode (500 pieces)	SE-9AU-M2		
			Copper electrode (25 pieces)	SE-9CU-M		
			ITO electrode (25 pieces)	SE-9IT-M		
			Nickel electrode (25 pieces)	SE-9NI-M		
			Platinum electrode (25 pieces)	SE-9PT-M		
			Silicon electrode (25 pieces)	SE-9SI-M		
			Stainless steel (SS304) electrode (25 pieces)	SE-9SS-M		
			Titanium electrode (25 pieces)	SE-9TI-M		
		Separated lead wire	Gold electrode (25 pieces)	SE-9AU-MS		
			ITO electrode (25 pieces)	SE-9IT-MS		
			Platinum electrode (25 pieces)	SE-9PT-MS		
			Others			
			Resonators lead wire no sputter (50 pieces)			SE-LEAD/2
			Resonators 9 MHz, no sputter			
No lead wire (50 pieces)			SE-9Q/2			
Mirror finish no lead wire (50 pieces)			SE-9Q-M/2			
Control software for QCM922, QCM922A (no potentiostat control)			SE-WQCM			



Gold electrode



Platinum electrode




Stainless steel electrode

Measurement examples

Electropolymerization of pyrrol

The polypyrrol film was deposited on an Au-coated quartz using cyclic voltammetry (twenty cycles).



The quartz electrode was immersed in an acetonitrile solution (Bu_4NPF_6 0.2 mol/L) containing a solution of 1 methylpyrrol monomer (0.01 mol/L).

Polypyrrol film growth on the quartz working electrode

Fig. 1 represents polypyrrol film growth on the quartz electrode during successive cycles of cyclic voltammetry. The reversibility of the charge transfer in such a polymer film is often dependent on the deposition mode (quasi-reversible in this example). This growth is very regular but tends to slow down during the last cycles. This can be due to an interfacial depletion of the solution in methyl pyrrol species in the layer close to the electrode surface and to a saturation of the working electrode surface area.

(1): G. Sauerbrey, Phys. Verh., 1957, 8, 113–114.

(2): G. Sauerbrey, Z. Phys., 1959, 155, 206–222.

(3): Application note #13. Section "Apps & literature of EC-Lab division".


QCA922A specifications

Frequency range	5 MHz–30 MHz, resolution 0.01 Hz
Resonant resistance range	10–16 Ω , resolution 0.1 Ω
ΔF output	Full scale: ± 10 V (14 bit) ± 100 Hz / ± 500 kHz
Mass range	0.1 ng / Hz // 3.5 ng / Hz
Resistance output	Full scale: 0–10 V (14 bit) 10 Ω to 10 k Ω
Gate time	Variable (10 ms / 20 ms / 100 ms / 1 s / 10 s)
Interface	USB

QCM measurements during the film growth


Fig. 1 shows the resonant frequency decrease and the resonant resistance increase while the polymer film is growing. Moreover, the variation is dependent on the potential sweep resulting in a pseudo oscillation of frequency and resistance related to successive cycles. This plot can also be made versus potential (see fig. 1).

Fig. 1: overlaid frequency and current vs. E_{we} of the polymer film growth. Scanning at 100 mV/s between 0 and 1.018 V.




Here, it is interesting to note the frequency evolution versus potential on one cycle. The figure below shows the time evolution on one cycle.

Fig. 2: graphic zoom on one cycle showing the resonant frequency and the current density versus elapsed time⁽³⁾.



The mass calculation is carried automatically by the EC-Lab[®] process data tool. More details can be found in the application note⁽³⁾.




Surface finishing

	Roughness	Electrode materials deposition
Standard finish	0.6 μm	Sputtered
Mirror finish	0.06 μm	Sputtered

A high roughness means a large surface area. A low roughness means an exact surface area.

Battery Accessories.



Coin Cell Holders

	CCH	CCH-120	CCH-124	CCH-8	CCH-8
Cell max diameter/mm	24	20	24	24	
Cell height/mm	3	3.2	3.2	1.6 - 3.2	
Number of channels	4	1	1	8	
Measurement type	2 point	2 point	2 point	4 point	
To be used with	MPG2	BCS-805	BCS-805	VSP	BCS-805
	VMP3	BCS-810	BCS-810	VMP3 VSP-300 VMP-300 MPG2	BCS-810
Climatic chamber compatibility	No	No	No	Yes (-30 to 80°C)	
Catalog n°	092-22/14	096-120	096-124	092-22/24	096-128



Pouch Cell Holders

	PBH-125	PBH-150	PBH-4	PBH-8
Min leads separation distance/mm	0		12	
Max leads separation distance/mm	110*		44	
Number of channels	1		4	8
Max current/A	25	50	32	
Measurement type	4 point			
Receptacles diameter/mm	4 (power)	6 (power)**	4	
	2 (voltage)	4 (voltage)		
To be used with	All instruments			
Max operating T°/°C	80	100	80	
Size : HxWxD/mm	40x50x210***		135x325x180	135x650x180
Weight/kg	0.2***		1.9	3.8
Catalog n°	092-P25/1	092-P50/1	092-P32/4	092-P32/8

*Measured using the guide rail and the middle of the clamp.

Eyelet ring (The connection kit **094-110/CNT can be used for an easier connection to 6 mm diameter cables).


***Measured with the two clamps mounted on the guide rail.





BH-1i

2 and 4 mm receptacles are available for the current (power) cables. For voltage (sense) cables, only 2 mm receptacles are available.




CBH-8

4 mm receptacles are available for the current (power) and voltage (sense) cables.



PBH-4

4 mm receptacles are available for the current (power) and voltage (sense) cables.



PPBH-132

4 mm receptacles are available for the current (power) and voltage (sense) cables.



PPBH-1100

4 mm receptacles can be used for currents up to 32 A. For higher currents, the 6 mm flush mounting plugs should be used. These are compatible with FlexP 0160, HCV-3048, CC4-60A and CC8 cables (See page 40).

Cylindrical Cell Holders

	BH-1i	CBH-4	CBH-8
Cell max diameter/mm	26	60	
Cell min height/mm	0	30	
Cell max height/mm	76	100	
Number of channels	4	4	8
Max current/A	15	32	
Measurement type	4 point	4 point	
Receptacles diameter/mm	2 and 4	4	
To be used with	All instruments		
Max operating T°/°C	60	80	
Size : HxWxD/mm	205x150x95	335x260x150	335x520x150
Weight/kg	0.6	1.9	3.8
Catalog n°	092-22/15	092-C32/4	092-C32/8

Prismatic and Pouch Cell Holders

	PPBH-132	PPBH-1100
Cell min height/mm		0
Cell max height/mm		139
Min leads separation distance/mm		66
Max leads separation distance/mm		155
Number of channel		1
Max current/A	32	100
Measurement type	4 point	
Receptacles diameter/mm	4 (power and sense)	4 (power and sense) and 6 (power)
To be used with	All instruments	
Max operating T°/°C	80	
Size : HxWxD/mm	265x320x300	320x320x360
Weight/kg	3	5.1
Catalog n°	092-PC32/1	092-PC100/1

Battery Accessories.



CC4-200A



CC8

Current Collectors

BioLogic's current collectors offer the possibility to connect in parallel several channels and increase the maximum current that can be passed through the cell, in order to simplify and reduce the footprint of your setup.



CC4-60A

	CC4-60A	CC8	CC4-200A
Connection details			
Input			
Power cables/receptacles diameter/mm	4		6 (IP2x)
Voltage sense receptacles diameter/mm	2		4
Number of input channels	4	8	4
Max current/channel/A	15		50
Output			
Power receptacles diameter/mm	6 (IP2x)		8 (Amphenol, IP2x)
Voltage sense receptacles diameter/mm	2 (IP2x)		4 (IP2x)
Max output current/A	60	120	200
Cables details			
Output power cables	1 pair of 2 m power cables with 6 mm receptacles		1 pair of 2.5 m power cables with 8 mm receptacles and M8 threads
Output voltage cables	1 pair of 2 m sense cables with 2 mm banana plugs		1 pair of 2.5 m sense cables with 4 mm banana plugs
Instrument compatibility	BCS-815 VSP-300 VMP-300	BCS-815*	FlexP0160 FlexP0060 HCV-3048
Measurement type	4 point		
Max operating T°/°C	80		
Size (with feet) : HxWxD/mm	70x170x88	70x300x88	120x248x169
Weight/kg	3.8		
Catalog n°	096-022	096-015	093-100/CC4

*The CC8 comes with BCS tablets and cables. It is also compatible with the VSP-300 and VMP-300 and can be provided without cables and tablets using the following part number: 096-015/1.



SAM-50

Sense Adapter Module (SAM-50)

This can be added to a multichannel system to perform stack measurements up to 60 V for 5 channel boards and a 10-element measurement. 3 SAM-50s can be linked to follow-up 30 elements.

Sense Adapter Module	Catalog n°
SAM-50	092-26


Redox Flow Battery Cells.

This range of redox flow battery cells are manufactured by **Pinflow energy storage**.

The lab cells are specially designed to control the pressure applied on the carbon felts that are used as electrodes. Using rigid components and non-flat bipolar plates, you can not only perform reproducible experiments but you can also use various electrode thicknesses.

Cells stacks, complete turnkey set-ups with climatic chambers and flow control are also available.

Please ask your local reseller for more information.




Description	Catalog n°
Lab cells	
Redox Flow Lab Cell 20 cm ² and 3.75 mm spacer thickness	P-LABCELL/20
Redox Flow Lab Cell 50 cm ² and 3.75 mm spacer thickness	P-LABCELL/50
Spares	
Set of sealings for <25 cm ²	P-SSEAL/20
Pack of sealings for <25 cm ² cell, contains 10 sets	P-PSEAL/20
Set of bipolar plates for different compression for <25 cm ² cell	P-SBIPOL/20
Set of sealings for 50 cm ²	P-SSEAL/50
Pack of sealings for 50 cm ² cell, contains 10 sets	P-PSEAL/50
Set of bipolar plates for different compression for 50 cm ² cell	P-SBIPOL/50
Accessories	
Mounting stand	P-STAND
Felt cutter	P-CUTTER
Torque wrench	P-WRENCH
Set of fittings with new FPM O-rings	P-FPM

Pinflow cells were used and characterized in the following papers:

<https://doi.org/10.1016/j.memsci.2018.02.011>

<https://doi.org/10.1016/j.jpowsour.2018.01.079>

Connection Accessories.



DC3



Dummy cell for booster



Dummy cell for BCS

Dummy Cells

A specific dummy cell is available for boosters. This cell and DC3s have been designed to clamp onto one another, for ease of use, and to help you better manage your experiments.

Dummy Cell		Catalog n°
DC3	1 R+R/C+R/C circuit.	094-111/3
Dummy cell for booster	1 power resistor, 5 mΩ. Precision: 1% Temperature coefficient: ±50 ppm/°C	092-32/1
Dummy cell for BCS series	Three test resistors with 1% precision.	096-016/1

Test Boxes

Test Boxes		Catalog n°
Test Box 2	Several circuits with high precision resistors, for calibration and validation	092-22/6
Test Box 3	Three circuits: linear, two non-linear systems (Tafel & passivating) for teaching and demonstration	092-22/7

External Device Connection

External device connection	Catalog n°
DB9-8 BNC connector for auxiliary I/O	092-22/1
IS1 isolation module for auxiliary I/O for VMP-300 based instruments	094-081/5



IS1




DB9-8 BNC

Pt Probe	Catalog n°
PT100 temperature probe, to be connected to the auxiliary I/O, temperature range: -50 °C to 250 °C, Dimensions: 3 x 20 mm, Length of cable: 2.5 m, Accuracy: ±1 °C	092-22/13

Test Boxes

Test Box 2



Test Box 3



Transport Cases

Transport Cases	Catalog n°
SP-200	094-082
SP-240/SP-300	094-091

Specific Cables

By default, the potentiostat and the booster are provided with a 1.5 m long cell cable.

The cable connected from the booster to the potentiostat is 0.8 m long for VMP3 based instruments.

For some applications, the user may need different length cables. For this reason, longer cables are available.

For applications carried out in glove boxes, cell cables are also available.



Hermetic cell cable for glove box

Hermetic cell cable for glove box

	VMP3	VMP-300
Catalog n°	092-23/5	094-101/6
Content:		
Feedthrough type	12 pins	25 pins*
Inside glovebox (length: 1 m)	Cable with 2 mm connectors on one side and 12-pin Jaeger connector on the other side	Cable with electrometer on one side and 25-pin connector on the other side
Outside glovebox (length: 1.5 m)	Cable with SubD25 connector on one side and 12-pin Jaeger connector on the other side (length 1.5 m)	Cable with SubD25 connector on one side and 12-pin Jaeger connector on the other side (length 1.5 m)
Requirement:		
Hole to make in the glove box/mm	27	45

*Two feedthrough seals one installed in the glove box wall the other dedicated to the channel board calibration outside the box

Set-up connection

Bad connections can affect measurements (stability of potentiostat, artifacts etc).

In order to optimise your set-up, we recommend you use the accessories described in this section.

Connection cable from booster to potentiostat	Available length/m	Catalog n°
All boosters	3	092-33/5
	5	092-33/6

Longer cable	Available length/m	Catalog n°
SP-50, SP-150, VSP, VMP3	2.5	092-23/2
	3	092-23/7
	5	092-23/3
	10	092-23/4
Booster 2 A, 5 A for essential range (VMP3)	2.5	092-33/11
	3	092-33/12
	5	092-33/13
Booster 10 A, 20 A for essential range (VMP3)	2.5	092-33/21
	3	092-33/22
	5	092-33/23
BCS-815	0.25	096-011/1
	2.5	096-011/2
	5	096-011/3
	10	096-011/4

Connectors

2 mm alligator clip (pack of 10)

4 mm alligator clip (pack of 5)

2 mm receptacle (pack of 10)

4 mm receptacle (pack of 10)

2 mm banana plug (pack of 10)

4 mm banana plug (pack of 5)

2 mm receptacle to 4 mm plug adapter (pack of 5)

2 mm receptacle to 4 mm plug adapter for banana plug (pack of 5)


4 mm receptacle to 2 mm plug adapter for banana plug (pack of 5)

Connection kits

For standard board	- 4 alligator clips of 2 mm: blue, white, red, black - 3 receptacles of 2 mm: blue, white, red	092-1001/30
For booster board	- 3 alligator clips of 2 mm: blue, white, red - 2 alligator clips of 4 mm: red, black - 3 receptacles of 2 mm: blue, white, red - 2 receptacles of 4 mm: blue, white	092-1001/31

Catalog n°			
Black colour	Red colour	Blue colour	White colour
092-1001/1	092-1001/2	092-1001/26	092-1001/22
092-1001/13	092-1001/14	092-1001/23	
092-1001/5	092-1001/6	092-1001/7	092-1001/8
092-1001/25	092-1001/4	092-1001/3	
092-1001/9	092-1001/10	092-1001/11	092-1001/12
	092-1001/16	092-1001/15	092-1001/24
092-1001/18	092-1001/19	092-1001/20	092-1001/21
092-1001/27			
092-1001/28			

Connection Accessories.



Nstat box (8 channels)

Bipot cable: dedicated to RRDE applications

Multi-Electrode Investigation Cables

Multi-electrode investigation cables				
	Channel	Nb of channels	Length/m	Catalog n°
Nstat box (for VSP, VMP3)	standard	4	1.5	092-16
		8	1.5	092-22/3
	Low current	8	1.5	092-22/4
Bipot cable (for VSP, BiStat, VMP3)	Standard	2	1.5	092-22/12
		2	3	092-22/12A
		2	4	092-22/12B
Option				
External power supply for the Nstat box (this option is needed if more than one Nstat box is connected to a VMP3 or if the user uses VSP)				092-16/1

Faraday Cages



Advanced Faraday cage

To avoid any external perturbations, especially for low current applications, we recommend using a Faraday cage. Please note that in order to ensure that the cage is fully functional, it must be earthed by connecting it to the ground (this is done via a green plug on the instrument's rear panel).



FC-45 Faraday cage

Advanced Faraday cage specifications

Temperature/°C	0 - 50
Power supply	100 VAC-240 VAC, 50/60 Hz
Fuse/A	1
Gas pressure	< 34 kPa max.
Interface	IEEE-488, RS-232C
Weight/kg	3.8



Faraday cages	Catalog n°
FC-45 Faraday cage, 450x450x450 mm	094-084/1
Stand for FC-45	094-084/2
Standard Faraday cage, 400x200x600 mm	NS-FAR600
Advanced Faraday cage, 286x230x320 mm	A-012033

Scanning Product Accessories.

Probes

A range of probes dedicated for use with our SECM, SVP, SKP and LEIS scanning probe applications are available for the M370 and M470 systems. SECM probes can also be used with SECM150

Probes	Catalog n°
Fused silica-based SECM 10 µm diameter Platinum disk	U-23/10
Fused silica-based SECM 15 µm diameter Platinum disk	U-23/15
Fused silica-based SECM 25 µm diameter Platinum disk	U-23/25
Capillary based SECM 1 µm diameter Platinum disk	U-P5/1
Capillary based SECM 2 µm diameter Platinum disk	U-P5/2
Capillary based SECM 5 µm diameter Platinum disk	U-P5/5
Capillary based SECM 10 µm diameter Platinum disk	U-P5/10
Capillary based SECM 15 µm diameter Platinum disk	U-P5/15
Capillary based SECM 25 µm diameter Platinum disk	U-P5/25
SKP 500 µm diameter	U-SKP370/1
HR SKP 150 µm diameter	U-SKP-150
LEIS	U-LEIS370/1
SVP	U-SVP370/1



VCAM3 Video Microscope System

The VCAM3 is a long working distance video microscope which allows users to view the positioning between probe tip and sample surface in many scanning probe electrochemistry techniques.



VCAM3 specifications

Working Distance	108 mm
Min illumination/lux	0.0003
Field of view/mm	1.4 (x4.5) to 8.6 (x0.7)
Operation temperature/°C	-30 to +70
Catalog n°	U-VCAM3

Scanning Product Accessories.

Cells

Three cells are available :

The TriCell™ is a large volume, wide scan range cell, dedicated to LEIS, SVP, SKP, SDS techniques.

The μ TriCell™ and its Shallow version are dedicated to SECM techniques (dc, ac and ic mode). The Shallow μ TriCell™ contains a slightly smaller volume of electrolyte than the μ TriCell™ is more accessible and ideal for ic-SECM.

The Foil Cell has been designed for use with flat, foil type samples, such as those used for battery electrodes. It has been designed to mount directly on the baseplate of the μ TriCell™ and Shallow μ TriCell™.

All TriCells accommodate samples mounted in a 32 mm diameter resin cylinder.


Cells	Volume (mL)	Catalog n°
TriCell™	700	U-TRICELL
μ TriCell™	7	U-uTRICELL
Shallow μ TriCell™	6	U-SuTRICELL
Foil Cell	1	U-uFoilCell



TriCell™



μ TriCell™



Shallow μ TriCell™



Foil Cell™

M470 Glovebox Cables

The M470 Scanning Electrochemical Workstation is supplied with a full set of standard cables for use with all techniques.

For applications requiring use of the M470 in a glove box, additional cell cables are available.

These are supplied as a set of internal, feedthrough and external cables to replace a single cable.



Hermetic scan stage cable for glove box

	Electrometer	Piezo Strain Gauge	Piezo Drive	Scan Stage	3300
Content :					
Feedthrough Type/pins	12	6	3	8	25
Inside glove box	Cable with connector to electrometer on one side and 12-pin Jaeger connector on the other side (length 1.5 m)	Cable with connector to piezo strain gauge on one side and 6-pin Jaeger connector on the other side (length 1.5 m)	Cable with connector to piezo drive on one side and 3-pin Jaeger connector on the other side (length 1.5 m)	Cable with connector to scan stage on one side and 8-pin Jaeger connector on the other side (length 1.1 m)	Cable with 4 mm connectors on one side and 25-pin Jaeger connector on the other side (length 1 m)
Outside glove box	Cable with connector to SCV470 on one side and 12-pin Jaeger connector on the other side (length 1 m)	Cable with connector to SCV470 on one side and 6-pin Jaeger connector on the other side (length 1 m)	Cable with connector to SCV470 on one side and 3-pin Jaeger connector on the other side (length 1 m)	Cable with connector to SCV470 on one side and 8-pin Jaeger connector on the other side (length 1.1 m)	Cable with connector to SCV470 on one side and 25-pin Jaeger connector on the other side (length 1 m)
Requirement:					
Hole to make in the glove box/mm	27	21	21	27	45
Max. Required	1	1	1	3	2
Catalog n°	U-HC470ELE	U-HC470PSG	U-HC470PD	U-HC470STG	U-HC3300CL

USB-PIO

The USB-PIO, designed for use with the M470, allows external devices to be switched on and read. It can be used to control up to four different channels individually or collectively using the M470 software.

The USB-PIO can interface directly to user supplied cables, or to the supplied breakout PCB using the DB25 pin female connector.

USB-PIO Specifications


Temperature/°C	0 - 50
Power supply	100 VAC-240 VAC, 50/60 Hz
Fuse/A	1
Gas pressure	< 34 kPa max.
Interface	IEEE-488, RS-232C
Weight/kg	3.8



Material Testing Accessories.

How to identify the correct solution

As the electrical properties of materials depend on temperature, a Temperature Control Unit and a Sample Holder are needed to control the temperature and hold the sample (solid, pasty and liquid sample) between the parallel plates.



Temperature Control Units

High Temperature Furnace (HTF-1100)

HTF-1100 is an horizontal laboratory tube furnace dedicated to the electrical characterization of materials and to heat treatment in the temperature range between the ambient and 1100 °C.



HTF-1100 & HTSH-1100



ITS

Intermediate Temperature System (ITS)

ITS is a compact temperature chamber dedicated to spectroscopy characterization of electrical properties of material by Impedance spectroscopy under controlled atmospheres in the temperature range between -35 °C and 150 °C.

Temperature Control Unit	Operating Temp.	Features	Catalog n°
HTF-1100	RT to 1100 °C	Heating rate adjustable K-type thermocouple	097-110
In-plane ITS	-35 to 150 °C	Temperature accuracy: 0.3 °C	097-140/11
Through-plane ITS	-35 to 150 °C	PT1000 probes	097-140/12





For more information visit
www.biologic.net



Application notes



White papers




Tutorials



Videos

BioLogic Science Instruments

serves its customers worldwide through its affiliate offices and an extensive distribution network



Headquarters

BioLogic SAS

4, rue de Vaucanson – 38 170 Seyssinet-Pariset
France

Phone: +33 476 98 68 31

Fax: +33 476 98 69 09

Affiliate offices

BioLogic USA, LLC

P.O.Box 30009 – 37 930 Knoxville
USA

Phone: +1 865 769 3800

Fax: +1 865 769 3801

BioLogic Pvt Ltd

Unit No.204, Odyssey IT Park, Road No. 9, MIDC
Wagle Estate, Thane West, Mumbai-400604 MH
India

Phone: +91 022 25829074

BioLogic Middle East

DMCC Office – No. 1206, Indigo Icon Tower, Cluster F,
Dubai JLT

United Arab Emirates

Phone: +971 455 151 82

Fax: +971 455 151 63

BioLogic Singapore

No.18, Boon Lay Way #02-105 (A), Tradehub 21,
Singapore

Phone: +65 92335838

BioLogic Science Instruments Ltd

Bradbury House, Market Street, Glossop
Derbyshire, SK13 8AR

United Kingdom

Phone: +44 (0) 333 012 4056

BioLogic Spain

Calle Diputació, 262, 1 1, Barcelona 08007
Spain

Phone: +34 656 876773

BioLogic Science Instruments GmbH

Rodeweg 20, D-37081 Göttingen
Germany

Phone: +49 (0)551 38266900