

Nitrogen/Protein Determination in Flours by the Thermo Scientific FLASH 2000 Analyzer using Argon as Carrier Gas

Dr. Liliana Krotz, OEA Product Specialist, Walter Galotta, OEA Product Support Engineer, and Dr. Guido Giazzi, OEA Product Manager, Thermo Fisher Scientific, Milan, Italy

Key Words

Argon, combustion, flour, nitrogen, protein

Goal

This application note presents data on nitrogen/protein determination in flour reference materials with different nitrogen concentrations to demonstrate the performance of the system using argon as an alternative carrier gas and show the reproducibility of the results obtained.

Introduction

Flour is a fine powder made by grinding cereal grains or other suitable edible plant matter which contains starch as polysaccharides and protein. It is most commonly made from wheat but also rye, barley, rice and corn (similar to “polenta”). Ground legumes and nuts, such as soy, peanuts, almonds, and other tree nuts are also called flours. It is the key ingredient of bread; a staple food in many countries, and its availability has often been a major economic and political issue.

Regulations in many countries require specific control of flour quality. One of the most important parameters is protein content, which can be monitored through the precise and accurate determination of the levels of nitrogen in the product to determine nutritional quality.



Helium supply has been beset with problems over the last few years, with world-wide shortages and subsequent price increases. It has therefore been necessary to test an alternative gas, argon, which is readily available.

The Thermo Scientific™ FLASH 2000 analyzer (Figure 1) copes effortlessly with a wide array of laboratory requirements such as accuracy, day by day reproducibility and high sample throughput. It operates according to the dynamic flash combustion of the sample and typically uses helium as the carrier gas.



Figure 1. FLASH 2000 Elemental Analyzer

Methods

Samples are weighed in tin capsules and introduced into the combustion reactor via the Thermo Scientific™ MAS™ 200R autosampler together with the proper amount of oxygen. After combustion, the resultant gases are carried by an argon flow to a second reactor filled with copper, then swept through CO₂ and H₂O traps, a GC column and are finally detected by a thermal conductivity detector (TCD). The analytical configuration as well as the TCD detector are the same as those used with helium as the carrier gas (see Figure 2).

A complete report is automatically generated by the Thermo Scientific™ Eager Xperience dedicated data handling software and displayed at the end of the analysis. The Eager Xperience software provides a new option AGO (Argon Gas Option) which modifies the argon carrier flow during the run to optimize the analysis.




Figure 2 – FLASH 2000 Nitrogen configuration

Analytical conditions

Combustion Furnace Temperature:	950 °C
Reduction Furnace Temperature:	840 °C
Oven Temperature:	50 °C (GC column inside the oven)
Argon Carrier Flow:	60 ml/min
Argon reference Flow:	60 ml/min
Oxygen Flow:	300 ml/min
Oxygen Injection End:	30 sec for standard, Category B (OxyTune) for flours and 30 sec for soy bean meal
Sample Delay:	10 sec
Run time:	10 mins

Official methods requirements for Protein determination by combustion

AOAC (Method 992.23) and AACC (Method 46-30, 1999) indicate that a suitable fineness of grind must be determined (for each different material analyzed) to achieve precision that gives RSD of $\leq 2\%$ for 10 successive determinations of nitrogen.

Results

To evaluate the performance of the system for flour analysis, samples with different nitrogen content were chosen. Instrument calibration was performed with approximately 50 - 70 mg of Aspartic acid standard (10.52 N %), using K factor as calibration method. To evaluate the calibration, Thermo Scientific™ Pasta Reference Material was analyzed. The protein factor 6.25 was used to calculate the protein content.

Table 1 shows the nitrogen/protein data obtained of the analysis of Pasta Reference Material. The sample weight was 150 – 160 mg. The certified N % is 2.227 and the uncertainty declared by the supplier is 0.097.

Table 1 – Nitrogen/Protein data of Pasta Reference Material

N %	RSD %	Protein %	RSD %
2.21		13.81	
2.23		13.94	
2.19		13.69	
2.20		13.75	
2.22		13.87	
2.21	0.542	13.81	0.535
2.20		13.75	
2.22		13.87	
2.21		13.81	
2.20		13.75	

Five Flour Reference Materials from 1.36 to 7.55 N % were chosen to correlate the experimental results to the expected values. Table 2 shows the certified N % and the uncertainty declared by the supplier.

Sample Description	Supplier Specification	
	N %	Uncertainty (\pm)
Wheat Flour Reference Material	1.36	0.25
Rice Flour Reference Material	1.38	0.05
Barley Flour Reference Material	1.9	0.04
Oatmeal Reference Material	1.9	0.10
Soy Bean Meal Reference Material	7.5	0.05

Table 3 shows the experimental N/Protein data obtained which were very satisfactory. The sample weight was 140 – 150 mg.

In all cases the RSD % was less than 2 % as requested in the Official Methods.

Table 3 – Experimental N % of Flours Reference Material

Sample Determination	Wheat Flour		Rice Flour		Barley Flour		Oatmeal	
	N %	Protein %	N %	Protein %	N %	Protein %	N %	Protein %
%	1.37	8.57	1.36	8.49	1.89	11.80	1.88	11.76
	1.36	8.51	1.37	8.59	1.91	11.94	1.88	11.73
	1.34	8.36	1.40	8.73	1.91	11.91	1.86	11.61
	1.35	8.41	1.40	8.74	1.91	11.95	1.90	11.90
	1.36	8.52	1.40	8.75	1.94	12.10	1.86	11.64
	1.36	8.49	1.38	8.65	1.91	11.91	1.85	11.56
	1.33	8.28	1.38	8.65	1.87	11.66	1.92	11.97
	1.36	8.51	1.38	8.62	1.89	11.80	1.87	11.67
	1.33	8.30	1.41	8.83	1.90	11.89	1.92	12.01
	1.36	8.52	1.39	8.69	1.92	11.97	1.89	11.79
Average %	1.35	8.45	1.39	8.67	1.90	11.89	1.88	11.76
RSD %	1.034	1.210	1.130	1.106	0.997	0.997	1.302	1.307

Table 5 – Experimental N % of Soy Bean Meal Reference Material (130 – 190 mg)

Weight (mg)	N %	Av. N %	RSD %	Protein %	Av. Protein %	RSD %
130.48	7.50	7.50	0.077	46.86	46.87	0.077
132.67	7.51			46.91		
129.90	7.50			46.84		
140.57	7.49	7.49	0.534	46.81	46.81	0.513
141.41	7.53			47.05		
140.28	7.45			46.57		
161.30	7.48	7.52	0.682	46.73	47.02	0.703
160.29	7.51			46.95		
162.44	7.58			47.38		
170.53	7.47	7.50	0.407	46.66	46.88	0.431
171.02	7.53			47.06		
171.22	7.51			46.91		
182.89	7.49	7.53	0.413	46.81	47.09	0.417
182.54	7.56			47.23		
180.50	7.55			47.22		
91.57	7.54			47.12		

Conclusions

Good repeatability, accuracy and precision was obtained with the FLASH 2000 Nitrogen Analyzer using argon as the carrier gas. The RSD % obtained was less than 2% according to the AOAC and AACC Official Methods Performance Requirements.

No memory effect was observed when changing the type of sample, indicating complete combustion and detection of the element.

www.thermoscientific.com

©2014 Thermo Fisher Scientific Inc. All rights reserved. AOAC is a trademark of The Association of Official Analytical Chemists. AACC is a trademark of The American Association of Cereal Chemists. ISO is a trademark of the International Standards Organization. All other trademarks are the property of Thermo Fisher Scientific and its subsidiaries. This information is presented as an example of the capabilities of Thermo Fisher Scientific products. It is not intended to encourage use of these products in any manners that might infringe the intellectual property rights of others. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

**COMPANY WITH QUALITY MANAGEMENT
SYSTEM CERTIFIED BY DNV
= ISO 9001:2008 =**

*Thermo Fisher Scientific S.p.A. Milan, Italy
is ISO 9001:2008 Certified.*

Africa +43 1 333 50 34 0	Denmark +45 70 23 62 60	Japan +81 45 453 9100	Singapore +65 6289 1190
Australia +61 3 9757 4300	Europe-Other +43 1 333 50 34 0	Latin America +1 561 688 8700	Spain +34 914 845 965
Austria +43 810 282 206	Finland +358 9 3291 0200	Middle East +43 1 333 50 34 0	Sweden +46 8 556 468 00
Belgium +32 53 73 42 41	France +33 1 60 92 48 00	Netherlands +31 76 579 55 55	Switzerland +41 61 716 77 00
Canada +1 800 530 8447	Germany +49 6103 408 1014	New Zealand +64 9 980 6700	UK +44 1442 233555
China 800 810 5118 (free call domestic) 400 650 5118	India +91 22 6742 9494	Norway +46 8 556 468 00	USA +1 800 532 4752
	Italy +39 02 950 591	Russia/CIS +43 1 333 50 34 0	

Thermo
SCIENTIFIC

Part of Thermo Fisher Scientific