Application proof note | 002340



**Environmental** 

## Fast determination of anions in municipal drinking water using a compact ion chromatography system with a chemically regenerated suppressor

#### Authors

Hua Yang, Jeffrey Rohrer Thermo Fisher Scientific, Sunnyvale, CA, USA

#### **Keywords**

Dionex IonPac AS22-Fast-4µm column, suppressed conductivity detection, Inuvion system, Safe Drinking Water Act (SDWA) and Clean Water Act (CWA), U.S. EPA Method 300.1, Dionex ACRS 500 Anion Chemically Regenerated Suppressor

#### Introduction

lon chromatography (IC) is a well-accepted technique for monitoring inorganic anions in water, including surface, ground, drinking, and wastewaters. In the U.S., water quality is regulated through the Safe Drinking Water Act (SDWA) and the Clean Water Act (CWA) and enforced through the United States Environmental Protection Agency (U.S. EPA)<sup>1</sup>. Ion chromatography methods have been approved by the U.S. EPA for compliance monitoring of inorganic anions in drinking water since the 1980s through U.S. EPA Method 300.0 that was updated in 1997 to U.S. EPA Method 300.1<sup>2</sup>. Various IC methods for water analysis have been demonstrated in Thermo Scientific application notes using standard or microbore flow rate columns with both carbonate/bicarbonate and hydroxide eluents.<sup>3</sup>

This work demonstrates the determination of inorganic anions in drinking water by IC using a Thermo Scientific<sup>™</sup> Dionex<sup>™</sup> IonPac<sup>™</sup> AS22-Fast-4µm column<sup>4</sup> set on a Thermo Scientific<sup>™</sup> Dionex<sup>™</sup> Inuvion ion chromatography system with a peristaltic pump used for regenerant flow to a chemical suppressor. The Dionex Inuvion IC system is an integrated, single-channel compact IC system. Coupled to the Thermo Scientific<sup>™</sup> Dionex<sup>™</sup> AS-DV autosampler and using a chemically regenerated suppressor, it provides a fast and economical choice for routine water analysis.

# thermo scientific

#### Method

#### Reagent, sample, and standard

- Degassed deionized (DI) water, 18 MΩ·cm resistance or better. •
- Drinking water •
- Thermo Scientific<sup>™</sup> Dionex<sup>™</sup> Combined Seven Anion Standard II (P/N 057590)

#### System preparation and setup


Figure 1 shows the flow diagram of a Dionex Inuvion IC system using the peristaltic pump to deliver suppressor regenerant.


#### Instrument method parameters


Instrument	Dionex Inuvion system (P/N 22185-60108)	
Autosampler	Thermo Scientific <sup>™</sup> Dionex <sup>™</sup> AS-DV autosampler (P/N 068907) with 5 mL Thermo Scientific <sup>™</sup> Dionex <sup>™</sup> PolyVials <sup>™</sup> and filter caps (P/N 038141)	
Columns	Thermo Scientific <sup>™</sup> Dionex <sup>™</sup> IonPac <sup>™</sup> AS22-Fast-4µm and AG22-Fast-4µm, 4 mm i.d. column set (P/N 088487, 088486)	
Eluent	4.5 mM Na <sub>2</sub> CO <sub>3</sub> / 1.4 mM NaHCO <sub>3</sub> (Prepared from Thermo Scientific <sup>™</sup> Dionex <sup>™</sup> AS22 Eluent Concentrate, P/N 063965)	
Eluent flow rate	2.0 mL/min	
lnj. volume	10 μL (full loop)	
Column temp.	30 °C	
Detection	Suppressed conductivity, Thermo Scientific <sup>™</sup> Dionex <sup>™</sup> ACRS 500 Anion Chemically Regenerated Suppressor, 4 mm (P/N 085090)	
Regenerant	50 mM sulfuric acid (prepared from Thermo Scientific™ Dionex™ Anion Regenerant Concentrate (P/N 057555)	
Regenerant flow rate	2.0 mL/min (delivered by the peristatic pump in the Dionex Inuvion system)	
Background conductance	~21 µS/cm	
System backpressure	~1,850 psi (100 psi = 0.6894 MPa)	
Run time	5 min	
Software	Thermo Scientific <sup>™</sup> Chromeleon <sup>™</sup> Chromatography Data System (CDS) software version 7.3.2	

#### **Results**

Figure 2 shows a separation of inorganic anions within 5 min using the Dionex IonPac AS22-Fast-4µm column. As this figure shows, seven inorganic anions are well resolved. The Dionex IonPac AS22-Fast-4µm column can be used for compliance monitoring of inorganic anions in water. This 5-min method allows additional samples to be run every hour compared to most IC methods used to determine anions in drinking water.







AS-DV autosampler

Figure 1. Illustration of the IC system flow diagram

Figure 3 shows the determination of inorganic anions in a drinking water sample. Among the common anions, fluoride, nitrate, and nitrite are regulated with the maximum contaminant level (MCL) for fluoride at 4 mg/L, nitrite at 1 mg/L, and nitrate at 10 mg/L. The results show that the drinking water sample contains chloride (9.5 mg/L), nitrate (1.2 mg/L), and sulfate (16.5 mg/L) with less than 1 mg/L of fluoride (0.6 mg/L), and nitrite (0.2 mg/L), which meet the regulatory criteria.

Sample:	Municipal dr	inking water
Peaks: 1. Fluoride 2. Chloride 3. Nitrite	Min mg/l 1.2 0.6 1.6 9.5 2.0 0.2	
4. Unknow 5. Nitrate 6. Phospha 7. Sulfate	2.6 1.2	
2		5
	-0.2	3.0
		7

### Figure 3. Determination of inorganic anions in a drinking water sample

34

2



9

uS/cm

0

0

### Learn more at thermofisher.com/inuvion

5

Minutes

3

6

**General Laboratory Equipment – Not For Diagnostic Procedures.** © 2023 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. This information is presented as an example of the capabilities of Thermo Fisher Scientific products. It is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details. **AP002340-EN 0623S** 

#### Conclusion

This work shows the determination of inorganic anions in a drinking water sample using a Dionex Inuvion IC system with a peristaltic pump. The integrated Dionex Inuvion IC system coupled with the Dionex AS-DV autosampler, Dionex IonPac AS22-Fast-4µm columns, and chemically regenerated suppressor provides a fast, simple, and Iow-cost instrument setup for routine determination of inorganic anions in drinking water.

#### References

7

4.0

5

6

- 1. National Primary Drinking Water Regulations (Accessed May 31, 2023.)
- U.S. EPA Method 300.1. The Determination of Inorganic Anions in Water by Ion Chromatography; rev 1.0; U.S. EPA, Office of Water: Cincinnati, OH, 1997. (Accessed May 31, 2023.)
- 3. Inorganic Anions Analysis by EPA 300.0 & 300.1. (Accessed May 31, 2023.)
- 4. The Dionex IonPac AS22-Fast-4µm IC columns. (Accessed May 31, 2023.)

### thermo scientific