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
MAX-OXT

Generate “Interference Spectra” from sample gas by selectively removing selected hydrocarbons

Does not destroy H_2O , CO_2 and CH_4 (primary atmospheric interferences)

“Interference Spectra” added to analysis algorithm

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
Novinka 2022: MAX ANALYTICAL

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<https://patents.justia.com/patent/20200292447>


Background & sample spectrum calculated from same Interferogram

- High resolution sample (2 cm^{-1})
- Low resolution background (8 cm^{-1})
- No independent background or high purity zero gas is required
- Eliminates baseline drift and bias drift in low level measurements

Novinka 2022: MAX ANALYTICAL

Kombinace: MCT + Starboost + Autoref


Ethylene Oxid - EtO

**± 1µabs
Residual**

**~170 ppt EtO
or ~600 pg**

**H₂O Intf.
~12,000x**

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- H₂ as target fuel for renewable/green energy applications
 - Advantages over batteries for many applications
 - Industrial/power plant energy supply via pipeline (CNG + H₂)
 - Fuel cell power for semi trucks, buses, cars, (train/marine/air?)
 - Refillable chemical batteries
- Currently most hydrogen generated from petrochemical reactions – steam reforming
 - “Blue” hydrogen uses carbon capture from petrochemical production
 - Goal – use solar/wind energy to split H₂O as “green H₂”
- Europe has goal to be carbon free by 2050
 - Banning diesel/gasoline vehicles in cities by 2030?
 - Significant financial commitments –
 - EU 3.8B Euro investment annually
 - Source, European commission communication
 - Note: picture Toyota Mirai 2

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According to 78 % of automotive executives polled, Fuel Cell Electric Vehicles (FCEV)s are “the real breakthrough for e-mobility” and could represent 32 % of fuel demand by 2050. **Current barriers to mass implementation of hydrogen in transport arise from European Directive 2014/94/EU and International organisation of legal metrology (OIML) recommendations that must be met by all European hydrogen refuelling stations (HRS).**

Hydrogen

>> Hydrogen purity measurements according to ISO 14687-2 and risk assessment for fuel cells

>> Analytical methods for performing hydrogen purity testing to enable the full implementation of the revised ISO 14687-2 standard

>> Development and validation of traceable methods for mass measurements of hydrogen absorbed in metal hydrides

End May 2019

HYDRAITE

Activities & Achievements

Establishing 3 EU laboratories for Hydrogen quality analysis

- ✓ Capable ISO 14687
- ✓ Gas sampling
- ✓ Comparability

End Dec 2020

METROLOGY for HYDROGEN VEHICLES 2

- Flow metering: accurate metering of dispensed hydrogen through calibration with secondary standards
- Hydrogen fuel quality control: provide reference materials for laboratories to be compliant with ISO 14687 and ISO 21087
- Provide scientific input to standardisation of sampling of hydrogen from HRS
- Validation of online sensors for fuel quality control at HRS
- Harmonized methodology for evaluation of impact of impurities on FC stacks

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MAX-IR Starboost enhancement of detection limits
No need of liquid nitrogen or purge

Performance Specification	Procedure / Reference Document	CO2	CO	CH4	THC (C1)
Instrument Detection Limit (IDL)	3 σ for n=9 analyses of blank sample	0.2ppb	1.7ppb	0.7ppb	6.3ppb
Method Detection Limit (MDL _s)	EPA 821-R-16-006	0.7ppb	0.6ppb	1.1ppb	7.6ppb
Accuracy (Average % of Expected)	4 replicates at 8 levels (0-5ppm), randomized order	1.25%	2.91%	2.15%	4.34%
Linearity (R ²)	4 replicates at 8 levels (0-5ppm), randomized order	0.9999	0.9999	0.9998	0.9998
Precision (RSD)	30 replicates at 5-10 x MDL	0.70%	2.43%	0.88%	4.67%

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MAX-IR Starboost enhancement of detection limits

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