


WOW!



**EASY-Spray C18,
75 cm x 75 µm column
(P/N ES805)**

- Higher peak capacity
- Highest resolution for protein identifications
- Faster loading and equilibration

The Thermo Scientific™ EASY-Spray™ C18, 75 cm x 75 µm column is for those who want a simple, easy-to-use format to obtain the highest resolution sequence coverage in proteomics.

Questions

- Why should I use EASY-Spray columns?
- I already use a 50 cm column – why do I need a 75 cm column?
- Do I need to run up to 1200 bar?

THEY ARE SIMPLE TO SET UP AND USE, WITHOUT TRICKY CONNECTIONS. THEY ALSO INCLUDE INTEGRATED TEMP CONTROL.

THE HIGHER PEAK CAPACITY GIVES BETTER SEPARATION-HELPING TO SEE LOWER ABUNDANCE PROTEINS USING THE SAME GRADIENT!

YOU CAN RUN UP TO 1200 BAR TO CUT DOWN LOADING AND EQUILIBRATION TIMES, BUT IT WILL ALSO WORK AT LOWER PRESSURES.

**THE EASY-SPRAY C18,
75 CM x 75 µm COLUMN –
RAISING THE BAR IN
PROTEOMICS!**

Find out more at thermofisher.com/EASYSpray

ThermoFisher
SCIENTIFIC