

Viscotester iQ

VISCOTESTER IQ

NOVÉ
Viscotester iQ

VISCOTESTER IQ

NOVÉ
Viscotester iQ

INDIVIDUÁLNÍ. INTUITIVNÍ. INTELIGENTNÍ.



NOVÉ

Viscotester iQ - INDIVIDUÁLNÍ

- robustní design
- výměnné temperační jednotky
- široký výběr geometrií
- RheoApp



NOVÉ

Viscotester iQ - INTUITIVNÍ

- 'smart lift'
- přesná aplikace vzorku
- jednoduchá výměna temperačních jednotek a geometrií



VISCOTESTER IQ

NOVÉ

Viscotester iQ - INTELIGENTNÍ

- touch screen
- průvodce měřením
- connect assist
- temperature assist



VISCOTESTER IQ

NOVÉ Viscotester iQ

- plug 'n play temperační jednotky: Peltier, kapalina, -20 až 180°C*
- standalone (**ČESKY**) vs. RheoApp vs. RheoWin
- koaxiální válce, deska-deska, double gap, lopatkové rotory
- univerzální držák, horní kryt, jednorázové desky



NOVÉ

Viscotester iQ - aplikační listy

- V267 Rheological Behaviour Of Shower Gel
- V268 Spreadability Of Cream Cheese
- V269 Flow Behaviour Of Chocolate Melts Working
- V270 Measuring Yield Stress To Correlate Slump
- V271 Investigation Of The effect Of Fat Content On The Yield Stress Of Mayonnaise
- V272 Yield Stress Of Jam
- V273 QC Tests On Cosmetic Emulsions With Viscotester iQ
- V274 Investigating The Thixotropic And Shear Recovery Behavior
- V275 How Much Squeezing Power Is Required To Get The Tooth Paste Out Of The Tube

www.pragolab.cz

<http://www.pragolab.cz/download>



[linkedin.com/company/pragolab-s-r-o-](https://www.linkedin.com/company/pragolab-s-r-o-)



twitter.com/pragolab