



**Thermo Scientific**  
Chromatography Columns  
and Consumables 2016-2017

# Connected chromatography solutions

**Thermo**  
SCIENTIFIC

## All-Plastic Disposable Syringes

Disposable syringes with polyethylene barrels and polypropylene plungers; use for all syringe filter applications

- Two-part, all-plastic construction eliminates the need for rubber or synthetic plunger gaskets
- No silicone or oil lubricant is required in the barrel
- Choose Luer-Slip or Luer Lock syringes, in capacities ranging from 1 to 50mL



### All-Plastic Disposable Syringes

Capacity (mL)	Cat. No.	Quantity
<b>Luer-Slip Syringes</b>		
1	<b>S7510-1</b>	100 Pack
3	<b>S7510-3</b>	100 Pack
5	<b>S7510-5</b>	100 Pack
10	<b>S7510-10</b>	100 Pack
20	<b>S7510-20</b>	100 Pack
30	<b>S7510-30</b>	50 Pack
50	<b>S7510-50</b>	30 Pack
<b>Luer Lock Syringes</b>		
3	<b>S7515-3</b>	100 Pack
5	<b>S7515-5</b>	100 Pack
10	<b>S7515-10</b>	100 Pack
20	<b>S7515-20</b>	100 Pack

## 750µL Micro-Centrifugal Filters, Non-Sterile

Filter volumes as low as 50µL with low hold-up volume

- Filter volumes as low as 50µL up to 750µL with low hold-up volume
- Use with any laboratory microcentrifuge
- Virgin polypropylene filter housing with tapered 2mL, capped receiver tube
- 10,000xG maximum centrifugal force



### 750µL Micro-Centrifugal Filters, Non-Sterile

Material	Pore Size (µm)	Cat. No.	Quantity
Cellulose Acetate	0.2	<b>F2517-1</b>	100 Pack
Cellulose Acetate	0.45	<b>F2517-2</b>	100 Pack
Nylon	0.2	<b>F2517-3</b>	100 Pack
Nylon	0.45	<b>F2517-4</b>	100 Pack
PVDF	0.2	<b>F2517-5</b>	100 Pack
PVDF	0.45	<b>F2517-6</b>	100 Pack
Regenerated Cellulose	0.2	<b>F2517-7</b>	100 Pack
Regenerated Cellulose	0.45	<b>F2517-8</b>	100 Pack
PTFE	0.2	<b>F2517-9</b>	100 Pack
PTFE	0.45	<b>F2517-10</b>	100 Pack

## 2mL Centrifugal Filters, Non-Sterile

- Filter sample volumes up to 2mL
- Virgin polypropylene filter housing with tapered 5mL, capped receiver tube
- Use with benchtop or floor model centrifuges
- 5,000xG maximum centrifugal force



### 2mL Centrifugal Filters, Non-Sterile

Material	Pore Size (µm)	Cat. No.	Quantity
Cellulose Acetate	0.2	<b>F2520-1</b>	25 Pack
Cellulose Acetate	0.45	<b>F2520-2</b>	25 Pack
Nylon	0.2	<b>F2520-3</b>	25 Pack
Nylon	0.45	<b>F2520-4</b>	25 Pack
PVDF	0.2	<b>F2520-5</b>	25 Pack
PVDF	0.45	<b>F2520-6</b>	25 Pack
PTFE	0.2	<b>F2520-7</b>	25 Pack
PTFE	0.45	<b>F2520-8</b>	25 Pack

## 25mL Centrifugal Filters, Non-Sterile

- Filter sample volumes up to 25mL
- Virgin polypropylene filter housing with conical receiver
- Use with benchtop or floor model centrifuges
- 2,500xG maximum centrifugal force



### 25mL Centrifugal Filters, Non-Sterile

Material	Pore Size (µm)	Cat. No.	Quantity
Cellulose Acetate	0.2	<b>F2519-1</b>	50 Pack
Cellulose Acetate	0.45	<b>F2519-2</b>	50 Pack
Nylon	0.2	<b>F2519-3</b>	50 Pack
Nylon	0.45	<b>F2519-4</b>	50 Pack
PVDF	0.2	<b>F2519-5</b>	50 Pack
PVDF	0.45	<b>F2519-6</b>	50 Pack

# Resources

## for Chromatographers



### Chromatography Resource Center

Our web-based resource center provides technical support, applications, technical tips and literature to help move your separations forward.

Visit [www.thermoscientific.com/chromatography](http://www.thermoscientific.com/chromatography)

### How to order

#### USA and Canada

800 332 3331  
865 354 4616 fax  
[chromatography@thermofisher.com](mailto:chromatography@thermofisher.com)

#### France

+33 (0)1 60 92 48 34  
+33 (0)1 60 92 49 00 fax  
[Consommableslscgcf@thermofisher.com](mailto:Consommableslscgcf@thermofisher.com)

#### Germany

+49 (0) 2423 9431 -20 or -21  
+49 (0) 6103 408 1640 fax  
[lscg.consumables.DE@thermofisher.com](mailto:lscg.consumables.DE@thermofisher.com)

#### United Kingdom

+44 1928 534110  
+44 1928 534001 fax  
[lscg.consumables@thermofisher.com](mailto:lscg.consumables@thermofisher.com)

#### China

+86 21 6865 4588

#### India

+91 22 6742 9494  
+91 27 1766 2352  
+91 22 6742 9495 fax  
[analyze.in@thermofisher.com](mailto:analyze.in@thermofisher.com)

#### Japan

+81 3 5826 1615  
+81 3 5826 1611 fax  
[analyze.jp@thermofisher.com](mailto:analyze.jp@thermofisher.com)

#### Thermo Fisher Scientific Australia Pty Ltd

1300 735 292 (free call domestic)  
1800 067 639 fax  
[AUinfo@thermofisher.com](mailto:AUinfo@thermofisher.com)

#### Thermo Fisher Scientific New Zealand Ltd

0800 933 966 (free call domestic)  
0800 329 246 fax  
[NZinfo@thermofisher.com](mailto:NZinfo@thermofisher.com)

#### All Other Enquiries

+44 (0) 1928 534 050  
+44 (0) 1928 534 049 fax  
[salesorders.columns.uk@thermofisher.com](mailto:salesorders.columns.uk@thermofisher.com)

### Technical Support

**For advice and support,  
please visit our website:**

[www.thermoscientific.com/chromexpert](http://www.thermoscientific.com/chromexpert)

For more information visit: [www.thermoscientific.com/chromatography](http://www.thermoscientific.com/chromatography)

© 2015 Thermo Fisher Scientific Inc. All rights reserved.  
Specifications, terms and pricing are subject to change. Not all products are available in all countries.  
Please consult your local sales representative for details.

**Thermo**  
SCIENTIFIC

A Thermo Fisher Scientific Brand