

Thermo Scientific Precision and Lab-Line Water Baths

optimal reliability
outstanding performance

Limited Time Offer!¹
Save 10%

on Thermo Scientific™ Precision
and Lab-Line Water Baths



► **Choose from 29 Thermo Scientific
Precision and Lab-Line Water Baths –**

now with CE certification and more to come in 2014



Thermo
SCIENTIFIC

Contact your Thermo Scientific representative to save

10% off your order TODAY!*



Thermo Scientific Precision Water Baths



Thermo Scientific Lab-Line General Purpose Digital Water Bath

Thermo Scientific Precision Water Bath Specifications

Precision models come standard with a stainless steel gable cover²

	Part No.	Model	Volume	Maximum Temperature
Precision Analog General Purpose Water Baths	2824	180	1.5L	95 °C
	2828	181	2.5L	
	2832	182	5.5L	99.9 °C
	2836	183	12L	
	2840	184	19.5L	
	2844	185	18L	
	2848	186	43L	
	2852	188 (Dual Chamber)	12L/12L	
Precision Internal Circulating Water Baths	2865	260	19L	99.9 °C
	2867	265	34.5L	
	2869	270	89L	
Precision Shaking Water Baths	2871	Reciprocal 25	14.5L	99.9 °C
	2873	Reciprocal 50	26.5L	
	2875	Shallow Form	14.5L	
	2877	Dubnoff Reciprocal	14.5L	

² Part number 2828 comes standard with polypropylene gable cover

Thermo Scientific Lab-Line Water Bath Specifications

Lab-Line models come standard with acrylic gable cover for a more cost effective product. Stainless steel gable covers are available, please contact your local sales rep.

	Part No.	Volume	Maximum Temperature
Lab-Line Analog General Purpose Water Baths	18050A-1CEQ	2L	100°C
	18020A-1CEQ	2L Shallow	
	18000A-1CEQ	5L	
	18800A-1CEQ	5 /10L (Dual Chamber)	
	18005A-1CEQ	10L	
	18100A-1CEQ	20L	
	18900A-1CEQ	28L	
Lab-Line Digital General Purpose Water Baths	18052A-1CEQ	2L	100°C
	18022A-1CEQ	2L Shallow	
	18002A-1CEQ	5L	
	18802A-1CEQ	5 /10L (Dual Chamber)	
	18007A-1CEQ	10L	
	18102A-1CEQ	20L	
	18902A-1CEQ	28L	

¹Offer valid through January 6, 2014 through March 31, 2014. Void where prohibited. Not valid with any other offer, contract or agreement. All offers subject to Thermo Fisher standard terms and conditions of sale. Not valid in connection with any other offer or discount. Void where prohibited by law or company policy. See your sales representative for additional information. Savings is based on published list price. Please note that the customer shall accomplish any legal requirements regarding to the disclosure of a discount or other reduction in price, which could apply following the applicable law.

www.thermoscientific.com/tc

© 2014 Thermo Fisher Scientific Inc. All rights reserved. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

PFTCP LLPROMOEU 0114

WE OFFER
a full complement of
Temperature Control
products and support
solutions.

▶ Contact a sales representative today to learn how your success can circulate at every degree:
+49 721 4094 444 or email
Sales.tc.eu@thermofisher.com

Thermo
SCIENTIFIC

Part of Thermo Fisher Scientific