

## Comfortable rotor exchange with Sigma g-lock®



### Order information of g-lock components

The new g-lock rotor system from Sigma is available for the Sigma 3-16, 3-18 and 3-30 collection. To benefit from the comfortable g-lock system, you can order a new device with g-lock adapter plus g-lock rotors. Furthermore, you have the possibility to upgrade an existing device with the g-lock adapter and order the novel g-lock rotors.

Article-No.	Max. capacity	Suitable for			Max. RCF *
		3-16	3-18	3-30	
91124	4 x 400 ml	●	●		4.248 x g
91373	4 x 200 ml	●	●		5.242 x g
91551	6 plates	●	●		3.285 x g
91355	12 x 2 ml			●	65.403 x g
91123	24 x 2 ml	●	●	●	61.973 x g
91371	36 x 2 ml			●	70.121 x g
91520	48 x 2 ml	●	●	●	25.115 x g
91527	10 x 10 ml	●	●	●	57.438 x g
91534	12 x 15 ml	●	●	●	24.903 x g
91365	6 x 30 ml	●	●	●	60.628 x g
91375	6 x 50 ml	●	●	●	41.415 x g
91377	6 x 50 ml	●	●	●	25.248 x g
91532	4 x 94 ml	●	●	●	40.695 x g
91119	6 x 94 ml	●	●	●	26.323 x g

\* dependent on chosen centrifuge

### Reordering

To order a g-lock centrifuge, please specify in addition to the article number of the centrifuge the **article number 17999**.

Please see article numbers of g-lock rotors in the table and check Sigma 3-16, 3-18 and 3-30 brochures for further, more detailed information.

### g-lock upgrade

To upgrade your Sigma centrifuge with the g-lock system, please order the adapter with **article number 17997**, which can be installed at your site by a qualified service technician. Consequently, the centrifuge runs exclusively with the new g-lock rotor.



g-lock video  
 on the Web!  
 Scan the QR code or go  
 directly to the web page:  
[www.sigma-zentrifugen.de/  
 en/g-lock](http://www.sigma-zentrifugen.de/en/g-lock)

